## Explanation of Symbols

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## Explication des symboles

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## Verklaring van de symbolen

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Phillips-head screwdriver.

- Place one bed sensor pad or chair sensor pad within reach.
- Connect the monitor to the bed sensor pad or the chair sensor pad. Plug the cord end of the sensor pad into the bottom of the monitor. Pinch the cord together and route it through the cord protection recess on the monitor (near the plug input).
- Put pressure on the sensor pad to activate the pad and monitor. The monitor will beep two times to indicate that the sensor pad and fall monitor are in use. The "Monitor Status" light will blink while pressure remains on the pad. When pressure is removed, the fall monitor will sound the alarm. Press the "Reset" button to silence the alarm
- If the "Pad Lost Light" is on, check to make sure the cord is plugged in and reset the system.
- Test the system before use.

#### MOUNTING THE MONITOR

- The strap with buckle can be used to attach the monitor to a desired location.
- To mount the monitor, mount the screw onto the wall or headboard. Hang the monitor on the screw using the mounting hole on back of the monitor.

#### TURNING OFF OR RESETTING THE ALARM

- When the alarm is triggered, the alarm can be turned off and reset with either the reset button or Caregiver Key™ function.
- Select either "Reset Button ON" or "Caregiver Key™ ON" with the switch on the side of the monitor. Use the end of a pen to move the switch.
- If set on "Reset Button ON" position, the caregiver can reset the monitor by pressing the reset button located on the front of the monitor.
- If set on the "Caregiver Key ON" position, the monitor can only be reset by using the Caregiver Key provided with the monitor. The Caregiver Key must be placed over the red dot to reset the monitor.

#### VOLUME, TONE AND DELAY OPTIONS

- Adjust the volume by turning the dial located inside the battery compartment.
- Select an alarm tone by adjusting the switch inside the battery compartment to the desired setting.
- Select alarm tone delay time by adjusting the switch inside the battery compartment.

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#### OPTIONAL ITEMS

- AC Adaptor can be purchased separately.
- Nurse call port can be purchased separately if the patient's room does not supply one that accommodates this device.

#### SIZING

Universal

#### CARE INSTRUCTIONS

Wipe clean with disinfectant wipes. After cleaning, test device before each use.

**WARRANTY**  
DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

#### ESPANOL

Improving Care. Improving Business.™

**Part #0-1887**  
Revised 10/2014  
©2014 DeRoyal Industries, Inc.  
All Rights Reserved.

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