

## Cleaning and Maintenance

**WARNING:** Always unplug the unit prior to cleaning! Clean with a damp soft cloth. A non abrasive cleaning agent may be used for stubborn stains. Make sure no liquid enters the base of the unit.

The unit has been designed for use by people with special needs, however this does NOT mean it is unbreakable. You must ensure the unit is not hit with anything and that no one pulls or pushes the unit.

## Lamp Change

Not required! It uses an LED which lasts for years!

Please note: the maximum ambient room temperature in which this unit will operate manually is 27 degrees C

## Troubleshooting

Q. The unit doesn't work.

A1. Is it switched on? If so,

A2. Has the unit scrolled through its colors? If so,

A3. Please contact us for further assistance.



/Manufactured proudly in the UK by:

Experia



Acorn Phase 3, High Street, Grimethorpe  
Barnsley, South Yorkshire, S72 7BD  
United Kingdom

© Copyright 2010-2013 Experia  
All rights reserved.  
Experia is a registered trademark.



experia **USA**  
multisensory for everyone

For optimal results, please read and follow instructions carefully  
Please keep for future reference



Please follow all warnings on the product and within this manual. If upon receipt of this product, the user suspects that the product is defective, and/or requires service or repair due to malfunction or damage to the product, please contact the manufacturer immediately. Please be aware that improper installation and/or use or misuse of this product may result in injury, death, or adverse reactions to its user. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended uses of the product.

**Congratulations** on the purchase of your new Aurora LED projector!

The Aurora LED Projector has been designed specifically for people with special needs and should, if used in accordance with these guidelines, offer years of fun! The Aurora LED Projector offers interaction, calming, relaxation, theming and storytelling opportunities depending on the wheels used.

It uses a Super Bright LED

It is incredibly easy to use. Simply place on a table or mount securely on a wall, connect to an electrical outlet and switch on.

**Please ensure you read and follow these Instructions carefully before installing or operating the unit. Please keep this manual for future reference.**

## Features

- Ergonomic design
- Low noise level
- Latest LED technology. No changing bulbs EVER!
- Environmentally friendly – 3 x as efficient as lamp based products
- Cool running
- Low voltage – transformer supplied
- Low maintenance
- Durable & Easy to use
- Study desktop stand
- Floor, wall or ceiling mountable
- Works with all standard accessories
- Designed for all ages and abilities
- 1 Year manufacturer's guarantee

## Pack contents

- 1 x Aurora LED Projector
- 1 x Low Voltage Transformer


## Getting Started

1. Place the Aurora LED Projector on a flat solid surface or mount securely on a wall.
2. Connect the low voltage transformer connector to the socket on the back of the unit (right)..
3. Connect the standard 110V AC plug to a nearby outlet.
4. Switch on the unit using the switch on the back of the unit.
5. The unit will now illuminate.



## Fitting a Wheel Rotator & Effects Wheel

To enable images to be projected from the unit you will need to fit a wheel rotator and effects wheel – this is easily done.

1. First, locate your wheel rotator and choose the effect wheel you wish to use.
2. Remove both items from packaging. You will notice that on the front of the wheel rotator there is a metal spindle. On back of the effects wheel there is a plastic receptacle in the center of the wheel. Push the receptacle onto the metal spindle.
3. Looking at the projector you will see a slot towards the front of the unit. In here you will see 2 gates: 
4. Slide the wheel rotator into the gate nearest the back of the projector and connect the plug to the socket underneath.
4. Switch the unit on. The wheel should now rotate.
5. If the image that is being projected is blurred then simply adjust the black focusing lens on the front until it becomes clear.



If you are still not sure then visit our website at [www.experiainnovations.co.uk](http://www.experiainnovations.co.uk) and look in the Projector section or visit our YOUTUBE channel – ExperiaInnovations.

**Please note: Only use the transformer supplied. Use of a different transformer may result in damage to the product and void the warranty. For warranty information please see separate warranty card supplied.**