# Adjustable Stylus - #9086 & Replacement Tips #9084

**USER'S GUIDE** 

### Choose your length!

The Adjustable Stylus with Velcro straps will fit any hand and any positioning. Just adjust the length of the stylus and use on touch screen surfaces including smart phones, iPads and Android tablets. Replaceable tips available separately. Size: 1/2" Diameter x 2" to 19"L.

#### **Operation:**

1. The Adjustable Stylus has two adjustable Velcro straps, the Adjustable Stylus can be positioned between any finger and the straps can move up or down the length of Velcro attached to the actual Stylus. Once you have the position set you can now secure the Adjustable Stylus with the Velcro straps across the hand and wrist a shown in figure 1 & 2.



Fig.1

2. The Adjustable Stylus is on a telescopic rod you can adjust its length to fit your needs.



Fig.2

## **Replacement Tips:**

When it comes time to replace the Stylus tip you will first need to remove band of silver tape wrapped around the base of the Stylus. Next you will unscrew the tip and replace it with one of the replacement tips (sold separately). The replacement tips will include a new tape strip to re-apply to secure the tip to the Stylus. **Please Note:** The tape must be used to prevent the tip from becoming a choking hazard for young children.

#### Care of Unit:

The Adjustable Stylus can be wiped clean with any mild household multi-purpose cleaner and disinfectant.

**Do not submerge** the unit, as it will damage the contents of the top and the electrical components.