



UNNA-Z™ AND UNNA-Z™ STRETCH

ZINC PASTE BANDAGE


UNNA-Z™ and UNNA-Z™ STRETCH provide conformable, semi-rigid support and mild compression.




APPLICATION:



1. Position the foot at a right angle, so the foot is perpendicular to the leg. Begin wrapping without pressure just above the toes and continue upwards using a spiral technique with a 50% overlap.



2. Mold wrap to the contours of the leg. Finish wrapping 1-2 in. below the knee. Smooth out wrinkles as needed.



3. Secure by covering with a secondary dressing such as a self-adherent wrap, elastic bandage, or gauze roll.

CONTRAINDICATIONS: Advanced arterial occlusive disease, arterial ulcers, congestive heart failure, ankle-brachial pressure index (ABPI) less than 0.5 or greater than 1.3, or advanced sensory disturbance of the lower extremities.
WARNING: Do not use if you have a hypersensitivity to zinc or zinc oxide.

See product packaging for more information.

© 2020 Medline Industries, Inc. Medline is a registered trademark and Unna-Z is a trademark of Medline Industries, Inc.
MKT19W628799 / e20024

COMPRESSION SYSTEM