

Dynaform[®] Foot & Ankle Positioner and Shoe Holder Assembly Instructions


The Dynaform[®] Foot & Ankle Positioners and Shoe Holders are designed to provide the important benefits of proper ankle and foot positioning for a wide variety of users. We are committed to providing the very best service and will where deemed necessary by you, provide additional training in the assembly of the positioners. These devices are designed for positioning and this should be undertaken by suitably qualified and competent persons including rehabilitation engineers and occupational therapists.

Tools required: 1/8" (3mm) and 5/16" (4mm) Allen Keys, 1/4" Drill

Step 1: Determine the desired location for the device on the chair footplate.

Step 2: Mounting the Dynaform[®] Foot & Ankle Positioner


Once the proper location has been determined, using a 1/4" drill, drill through existing holes in footplate 'A' and through the chair footplate. After the first hole is drilled install a 1/4" flathead screw with washer and nylock nut to secure the footplate. Repeat the procedure for the second hole. Tighten the screws as specified.*



Dynaform[®] Foot Positioner

Step 3: Mounting the Shoe Holder

Once the proper location has been determined, lift the foam pad 'D' and using a 1/4" drill, drill through footplate 'B' and through the chair footplate. After the first hole is drilled install a 1/4" flathead screw with washer and nylock nut to secure the footplate. Repeat the procedure for a second hole, making sure that the second hole is separated from the first hole by 2 or more inches. Tighten the screws as specified.*



Shoe Holder

Transportation:

This device is not intended to take the place of or considered to be a substitute for other restraints required for transportation.

Limitations on use:

The product is compatible with a wide range of wheelchair systems. We will be glad to advise on appropriate usage.

*All screws to be tightened to 20in-lb (2.3 N).



Authorized Representative:
Medimark Europe BP2332
38033 Grenoble Cedex 2
France

This product is manufactured
to comply with the 'Medical
Device Directive'