
INSTRUCTIONS FOR USE

# Bariatric Slings

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

# Contents

Foreword .....	2	Attachment Methods.....	8
Customer Support.....	2	Metod 1 - Cross through (Legs crossed with cross-	
Definitions in this IFU .....	2	ing straps) .....	8
Intended Use .....	3	Metod 2 - Abduction, Leg Separation (Legs opened	
Patient/Resident Assessment .....	3	with non-crossing straps) .....	8
Expected Service Life .....	3	Applying the Sling.....	9
Safety Instructions .....	4	In Bed (22 steps).....	9
Safety Practises .....	4	In Chair/ Wheelchair (19 steps).....	10
Serious incident.....	4	On Floor (35 steps) .....	10
Parts Designation .....	5	Removing the Sling .....	12
Preparations .....	5	In Bed (12 steps).....	12
Before First Use (6 steps) .....	5	In Chair/Wheelchair (13 steps).....	12
Before Every Use (5 steps) .....	5	Cleaning and Disinfection .....	13
Select Sling Size.....	6	Cleaning Instruction (7 steps) .....	13
Using the Arjo Measuring Tape		Cleaning Chemicals .....	13
(4 steps) (Accessory) .....	6	Disinfection.....	13
Sizing Without a Measuring Tape (2 steps).....	6	Care and Preventive Maintenance .....	13
Sling Selection.....	6	Before and After Every Use.....	13
Transfer/Transport.....	6	When Soiled or Stained and Between Patients ...	13
Sling Loop Selection.....	6	Storage.....	13
Allowed Combinations.....	7	Service and Maintenance.....	13
Safe Working Load (SWL).....	7	Troubleshooting .....	14
Loop Attachment & Detachment.....	7	Technical Specifications .....	15
Attach the Loops (5 Steps).....	7	Label on the Sling.....	16
Detach the Loops (2 steps).....	7		

## WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product. Mandatory to read the Instructions for Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies. © Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

## Foreword

**Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly!**

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

### Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

### Definitions in this IFU

#### WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

#### CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

#### NOTE

Means: This is important information for the correct use of this system or equipment.

# Intended Use

---

The *Bariatric Sling* is a product intended for assisted transfer of bariatric patients/residents with limited ability to move and require full body support during transfer. The *Bariatric Sling* should be used together with Arjo lift devices in accordance with the Allowed Combinations specified in the Instructions for Use (IFU).

The Hammock model MAA8030 is required to be used for reclined position only.

The *Bariatric Sling* shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, and in accordance with the instructions outlined in the *Instructions for Use*.

The *Bariatric Sling* is intended to be used in hospital environments, nursing homes, other healthcare facilities and home care.

The *Bariatric Sling* should only be used for the purpose specified in this Instructions for Use. Any other use is prohibited.

## Patient/Resident Assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident sits in wheelchair
- No capacity to support herself at all
- Cannot stand unsupported and is not able to bear weight, not even partially
- Dependent on caregiver in most situations
- Physically demanding for caregiver

• Stimulation of remaining abilities is very important  
Or, for a patient/resident who:

- Is passive
  - Might be almost completely bed ridden
  - Often stiff, contracted joints
  - Totally dependent on the caregiver
  - Physically demanding for caregiver
  - Stimulation and activation is not a primary goal
- Safe Working Load (SWL) for Bariatric Sling:
- MAA8000, MAA8000A, MAA8010, MAA8020, MAA8030: 454 kg (1000 lbs).

The right type and size of slings should be used after proper assessment of each patient/resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

## Expected Service Life

The expected service time of Bariatric Sling is the maximum period of useful life.

The expected life of the sling is dependent on the actual use conditions. Therefore, before use, always make sure that the sling does not show signs of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking). If any such damage is observed, do not use the sling. If you have any doubts about sling safety, as a precaution and to ensure safety, do not use the sling.

The expected service life of the Bariatric Slings:

- Service Time: 2 years
- Shelf Time: 5 years

# Safety Instructions

## WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

## WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

## WARNING

To avoid injury, always assess the patient prior to use.

## WARNING

To avoid falling, make sure that the user weight is lower than the safe working load for all products or accessories being used.

## WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the *IFU*, it can cause weakness and breakage in the material.

## WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.



## Safety Practises

- At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.
- Patient with spasm can be lifted, but great care should be taken to support the patient's legs.

## Serious incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

# Parts Designation

1. Sling Label  
(Located on the outside of the sling)
2. Stiffener pocket
3. Stiffener
4. Shoulder strap with loops
5. Leg strap with loop

6. Leg flap
7. Body section
8. Head section
9. Pulling handle

Example of slings. Not all models are represented.

EN

## Preparations

### Before First Use (6 steps)

1. Check all parts of the sling, see section “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this IFU thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the IFU should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

#### NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

### Before Every Use (5 steps)

1. Patient assessment must always determine the methods used in everyday care, see “Intended Use” on page 3.

#### WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the IFU.

2. Check the allowed combinations against the sling, spreader bar and lift, see “Allowed Combinations” on page 7.

#### WARNING

To avoid injury, always make sure to inspect the equipment prior to use.

3. Check all parts of the sling, see “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
  - fraying
  - loose stitching
  - tears
  - fabric holes
  - soiled fabric
  - damaged loops
  - unreadable or damaged label

#### WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

4. Check that the sling is clean. If the sling is not clean see “Cleaning and Disinfection” on page 13.
5. Before using the sling, read the lift IFU for transfer and transport.

# Select Sling Size

## Using the Arjo Measuring Tape (4 steps) (Accessory)

The Arjo Measuring Tape is a sizing accessory only intended as a guide and is an approximation. Patient width and body proportions should also be considered when selecting the correct size of Arjo sling.

EN

1. If possible, the patient should be in a seated position.
2. Place the point marked A on the measuring tape over the patient's coccyx/against seat (position A). (See Fig. 1) and (See Fig. 2)



Fig. 2

3. Measure from coccyx/seat (position A) to the top of the head (position B) (See Fig. 1).
4. The coloured area on the measuring tape that is levelled with the top of the patient's head indicates the required sling size. If the sling size falls in between two sizes, it is recommended to select the smaller size. (See Fig. 2)

## Sizing Without a Measuring Tape (2 steps)

1. Place the sling over the patient's back.
2. Make sure the sling covers the patient from the top of the head (position B) to the coccyx (position A). (See Fig. 1)

## Sling Selection

The patient's physical disabilities, weight distribution and general physique needs to be taken into consideration when selecting a sling.

### Transfer/Transport

Article number	Safe Working Load	Product Description	Size
MAA8000	454 kg (1000 lbs)	Standard	M, L, XL, XXL
MAA8000A	454 kg (1000 lbs)	Standard	M, L
MAA8010	454 kg (1000 lbs)	Basic	M, L, XL, XXL
MAA8020	454 kg (1000 lbs)	Divided leg	M, L, XL, XXL
MAA8030	454 kg (1000 lbs)	Hammock Basic (Reclined position only)	M, L, XL, XXL

Article numbers with "A-suffix" (e.g. MAA8000A) has a BAA (Buy American Act) compliant country of origin.

## Sling Loop Selection

Body position according to loop selections	Shoulder loop	Leg loop
	3	1
	2	1
	1	1
	1	2

# Allowed Combinations

## WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

## Safe Working Load (SWL)

Always follow the lowest SWL of the total system. E.g. *Tenor* lift/spreader bar has a SWL of 320 kg (705 lbs) and the MAA8000 has a SWL of 454 kg (1000 lbs). This means that the *Tenor* lift/spreader bar has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

	Lift	Tenor	Maxi Sky 1000	Maxi Sky 2 Plus
	SWL	320 kg (705 lbs)	454 kg (1000 lbs)	454 kg (1000 lbs)
	Spreader Bar	4 point loop bar	4 point bariatric spreader bar	spreader bar
Sling	SWL	Size	Size	Size
MAA8000	454 kg (1000 lbs)	M, L, XL, XXL	M, L, XL, XXL	M, L, XL, XXL
MAA8000A	454 kg (1000 lbs)	M, L	M, L	-
MAA8010	454 kg (1000 lbs)	M, L, XL, XXL	M, L, XL, XXL	M, L, XL, XXL
MAA8020	454 kg (1000 lbs)	M, L, XL, XXL	M, L, XL, XXL	M, L, XL, XXL
MAA8030	454 kg (1000 lbs)	M, L, XL, XXL	M, L, XL, XXL	-

EN

## Loop Attachment & Detachment

**4 point spreader bar:** First attach the shoulder loops in the front hooks. Then attach the leg loops in the back hooks

## WARNING

Do not cross the shoulder straps.

## Attach the Loops (5 Steps)

## WARNING

Use the same loop length (loop colour) for the shoulder straps and the same length (loop colour) for the leg straps.

1. Place the loop over the spring loaded latch (A). (See Fig. 3)
2. Pull the loop down to force the latch to open.
3. Make sure that the spring loaded latch closes completely with the loop inside. (See Fig. 3)
4. Make sure that the latch is moving freely.
5. Make sure loops and straps are not twisted.

## NOTE

The loops lengths on the straps are used to adjust the patient's position and comfort.

- Loops closer to the shoulders; more forward position
- Loops further away from the shoulders; more reclined position
- Loops closer to the legs; raised legs/bottom
- Loops further away from legs; lowered legs/bottom

## Detach the Loops (2 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the loop.

### Method 1

1. Press the spring loaded latch open.
2. Remove the loop.

### Method 2

1. Pull one side of the loop over the hook and spring loaded latch. (See Fig. 4)
2. Pull down the loop. (See Fig. 5)

Fig. 3


Fig. 4




Fig. 5


# Attachment Methods

## Metod 1 - Cross through (Legs crossed with crossing straps)

This method is recommended for most general transfers.


Sling models: MAA8000, MAA8010


Sling model: MAA8020

## Metod 2 - Abduction, Leg Separation (Legs opened with non-crossing straps)

In this method, the patient's legs are separated, which is convenient for toileting and hygiene care.

Sling models: MAA8000, MAA8010

Sling model: MAA8020

### WARNING

Method 2 might not be suitable for patients with limited upper body control as they can slide down and almost out of the sling or pitch forward when the most seated position is reached.

# Applying the Sling

## In Bed (22 steps)

1. Place the lift nearby.
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. If possible, adjust the bed to a seated position to ease breathing and for easier attachment of the sling.
4. Check that the stiffeners are completely inside the stiffener pockets, if any.
5. Fold the sling lengthwise along the centre line, outside facing outside. The sling label is located on the outside.

### WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

6. Log roll the patient. (Use local technique.) If log roll is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*
7. Place the folded sling over the patient's side. Make sure that the centre line is aligned with the patient's spine, starting at the coccyx. Place the leg flaps towards the legs.
8. Fold down the top of the sling and tuck it under the patient's body. (See Fig. 6)
9. Log roll the patient to the other side.
10. Pull out the remaining part of the sling from under the patients body.
11. Return the patient to his/her back.
12. Place the leg flaps underneath the patient's legs. To avoid skin injuries, use an Arjo sliding sheet/tube when placing the leg flaps underneath the legs. Make sure the leg flaps go all the way around and are not twisted. See respective *Sliding sheet/tube IFU*
13. Make sure that the: (See Fig. 7)

- sling is centred and flat without creases,
  - sling's head support covers the neck/head area,
  - sling pieces are not twisted underneath the patient.
14. Cross the leg straps. Pull one strap through the other. (See Fig. 8)

### WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

15. Make sure the patient's arms are inside the sling.
16. Place the lift next to the patient and apply the brakes.

### WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.


17. Attach the sling.
18. Slightly lift the patient to create tension in the sling. (See Fig. 9)

### WARNING

To avoid the patient from falling, make sure that the sling attachments are securely attached before and during the lifting process.

19. Make sure that:
  - All loops are securely attached
  - All straps are straight (not twisted)
  - The patient lays comfortably in the sling.
20. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
21. Release the brakes on the lift.
22. Transfer according to the lift *IFU*. Only transfer for short distances e.g a couple of meters (feet) otherwise use a wheelchair, stretcher or bed.

EN

## In Chair/ Wheelchair (19 steps)

### NOTE


Do NOT use a *Hammock sling* in chair/ wheelchair. Apply these in bed.

1. Place the lift nearby.
2. Apply the wheelchair brakes.
3. Check that the stiffeners are completely inside the stiffener pockets, if any.
4. Stand in front of the patient and lean the patient forward. The patient can use the arm rest of the chair to help out.
5. Place the sling over the patient's back and head with the inside of the sling against the patient. The sling label is located on the outside. If the sling doesn't fit, change to a bigger size.
6. Make sure that the centre line is aligned with the patient's spine and coccyx.
7. Tuck in the end of the sling at the coccyx or against the seat. If needed, use an Arjo sliding sheet/ tube for placing the sling under the patient. See respective *Sliding sheet/tube IFU*.
8. Lean the patient back.
9. Place the leg flaps under the patient's legs. Make sure the leg flaps are not twisted. To avoid skin injuries, use an Arjo sliding sheet/tube when placing the leg flaps. Make sure the leg flaps goes all the way around. See respective *Sliding sheet/tube IFU*. (See Fig. 10)
10. Cross the legs straps. Pull one strap through the other. (See Fig. 8)
11. Make sure that the:
  - sling is centred and flat without creases,
  - sling's head support covers the neck/head area,
  - inside of the sling is against patient's body and
  - sling pieces are not twisted underneath the patient.
12. Make sure the spreader bar is in a seated position and place the lift in front of the patient with the lift legs wide open.
13. Apply the brakes on the lift.

### WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

Fig. 10

14. Make sure the patient's arms are inside the sling.

### WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

15. Attach the sling. Adjust the spreader bar if needed.

### WARNING

Make sure straps are not caught by wheelchair or lift castors.

16. Pay attention to the spreader bar and slightly lift the patient to create tension in the sling.

### WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

17. Make sure that:
  - All loops are securely attached
  - All straps are straight (not twisted)
  - The patient lays comfortably in the sling.
18. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
19. Transfer and transport the patient according to the lift IFU.

## On Floor (35 steps)

1. Make a clinical assessment of the patient. Do not use the sling if the patient has a head, neck, spine or hip injury. (Follow local routine.)
2. Place a pillow under the patient's head. (See Fig. 11)
3. Make sure the patient is in an area that is accessible for the lift. If needed, use an Arjo sliding sheet/ tube. See respective *Sliding sheet/tube IFU*.
4. Check that the stiffeners are completely inside the stiffener pockets, if any.
5. **If the patient is able to sit up, continue to next step. If not able to sit up, proceed to step 9.**

*Continue with the steps on the next page.*

Fig. 11

6. **Patient able to sit up:** Stand behind the patient and help him/her to a seated position.
7. Place the sling over the patient's back and head with the inside of the sling against the patient. The sling label is located on the outside.
8. Make sure that the centre line is aligned with the patient's spine, starting at the coccyx or against the seat. **Continue to step 16.**
9. **Patient not able to sit up:** Fold the sling lengthwise along the centre line, outside facing outside. The sling label is located on the outside.
10. Log roll the patient. (Use local technique) **(See Fig. 12)** If log rolling is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*.
11. Place the folded sling over the patient's side. Make sure that the centre line is aligned with the patient's spine, starting at the coccyx. Place the leg flaps towards the legs. **(See Fig. 13)**
12. Fold down the top of the sling and tuck it under the patient's body. Use an Arjo sliding sheet/tube for placing the sling under the patient. See respective *Sliding sheet/tube IFU*.
13. Log roll the patient to the other side.
14. Pull out the remaining part of the sling from under the patient's body.
15. Return the patient to his/her back.
16. Make sure that the: **(See Fig. 13)**
  - sling is centred and flat without creases,
  - sling's head support covers the neck/head area,
  - sling pieces are not twisted underneath the patient.
17. Position the lift sideways with the lift legs wide open. One caregiver needs to be close to the patient and pay attention to the patient's head and hair. The other caregiver needs to lift up the patient's legs.
18. Position one lift leg next to the patient's head and the other lift leg under the patient's legs. Make sure the spreader bar is positioned over the patient's shoulders. **(See Fig. 14)**
19. Lower the spreader bar.

## WARNING

**To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.**

20. Place the leg flaps underneath the patient's legs. To

avoid skin injuries, use an Arjo sliding sheet/tube when placing the legs flaps. Make sure the leg flaps goes all the way around and are not twisted. See respective *Sliding sheet/tube IFU*.

## WARNING

**To avoid injury, make sure the patient's arms are placed inside of the sling.**

21. Make sure the patient's arms are inside the sling.
22. Attach the shoulder loops.
23. Lift the patient's head up slightly.
24. Move the lift/spreader bar closer to the patient's legs.
25. Apply the brakes on the lift.
26. Cross the leg straps. Pull one strap through the other. **(See Fig. 8)**
27. Attach the leg loops.
28. Slightly lift the patient to create tension in the sling.

## WARNING

**To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.**

29. Make sure that:
  - All loops are securely attached
  - All straps are straight (not twisted)
  - The patient lays comfortably in the sling.
30. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
31. Make sure that: **(See Fig. 15)**
  - the spreader bar is in a reclined position,
  - one caregiver pays attention the patient's head and
  - that the other caregiver holds up the patient's feet to prevent them from hitting the lift leg.
32. Remove the pillow.
33. Release the brakes on the lift.
34. Bring the lift legs together.
35. Transfer and transport the patient according to the *IFU*.

Fig. 12


Fig. 13

Fig. 14


Fig. 15

# Removing the Sling

## In Bed (12 steps)

1. Position the patient over the bed. (See Fig. 16)
2. Apply the bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. Make sure the spreader bar is in a reclined position.

### WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

4. Lower the patient onto the bed. Make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
5. Detach the sling.
6. Move the lift away from the patient.
7. Remove the leg flaps from under the patient's legs by folding them backwards and pull them out gently. Use an Arjo sliding sheet/tube to remove the leg flaps. See respective *Sliding sheet/tube IFU*.

### WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

8. Log roll the patient. (Use local technique) If log rolling is not possible, use an Arjo sliding sheet/tube for removing the sling. See respective *Sliding sheet/tube IFU*.
9. Tuck the sling under the patient's body.
10. Log roll the patient to the other side and remove the sling.
11. Return the patient to his/her back.
12. Make sure the patient is in a comfortable and secure position in the bed.

## In Chair/Wheelchair (13 steps)

1. Apply the brakes on the wheelchair.
2. Position the patient over the chair/wheelchair. (See Fig. 17)
3. Lower the patient into the chair/wheelchair.
4. Make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
5. Stand behind the chair and pull in the patient towards the chair by the handles on the side of the sling.
6. Make sure that the patient's lower back is positioned all the way in against chair/wheelchair to avoid the patient from sliding out.

### WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

7. Detach the sling.
8. Move the lift away from the patient
9. Pull out the leg flaps from under the patient's legs gently and place them backwards along the patient's side. Use an Arjo sliding sheet/tube to remove the leg flaps. See respective *sliding sheet/tube IFU*. (See Fig. 18).
10. Stand in front of the patient and lean the patient forward. The patient can use the arm rest of the chair to help out.
11. Pull out the sling.
12. If the patient can't help out, fold the sling backwards behind the patient's back and pull out the sling gently. Use an Arjo sliding sheet/tube to remove the sling. See respective *Sliding sheet/tube IFU*.
13. Lean the patient back up to a comfortable and secure position.

Fig. 16

Fig. 17



Fig. 18

# Cleaning and Disinfection

## WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

## WARNING

To avoid material damage and injury, clean and disinfect according to this IFU.

- No other chemicals are allowed.
- Never clean with chlorine.
- Chlorine will deteriorate the surface of the material.

## WARNING

To avoid injury, always remove the sling before disinfecting the lift.

All *Bariatric slings* should be cleaned when they are soiled and stained and between patients.

### Cleaning Instruction (7 steps)

1. Detach the sling from the lift.
2. Remove the stiffeners from the stiffeners pockets, if any.
3. Close all buckles and hook and loop attachment on the slings and belts prior to washing.
4. Check the symbols on the product label for cleaning recommendations.

5. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection, please wash according to local hygiene legislation. If allowed according to the product label, tumble dry in low temperature, max 60°C (140°F).
6. Do NOT
  - wash together with other items that have rough surfaces or sharp objects.
  - use any mechanical pressure, pressing or rolling
  - use bleach
  - use gas sterilization
  - use autoclave
  - dry clean
  - ironing
7. Place the stiffeners if any, back in to the sling before use.

EN

### Cleaning Chemicals

For all sling material use common commercial detergent without optical brightener. No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

### Disinfection

The only allowed disinfection is disinfection by washing.

# Care and Preventive Maintenance

## WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

### Before and After Every Use

#### Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching
- Tears
- Holes
- Discolouration or stains from bleaching
- Sling soiled or stained
- Label unreadable or damaged
- Steam

### When Soiled or Stained and Between Patients

#### Clean/Disinfect

The caregiver should make sure that the sling is cleaned according to “*Cleaning and Disinfection*” on page 13, when it is soiled or stained and between patients.

#### Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

#### Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

# Troubleshooting

Problem	Action
<p>The patient is not positioned correctly in the sling (e.g. one strap is shorter than the other).</p>	<ul style="list-style-type: none"> <li>• Make sure that both shoulder loops are attached to the spreader bar at the same loop length.</li> <li>• Make sure that both leg loops are attached to the spreader bar at the same loop length.</li> <li>• Make sure the straps are not twisted.</li> <li>• Make sure the patient is seated on the inside of the sling. The outside has a label for reference.</li> <li>• Make sure the patient is centred in the sling.</li> </ul>
<p>The patient feels discomfort in the leg area when seated in the sling.</p>	<p>Make sure that there are no creases on the leg flaps of the sling.</p>
<p>The sling is difficult to apply (e.g. during log roll).</p>	<p>Use an Arjo sliding sheet/tube according to instructions in respective <i>IFU</i>.</p>
<p>There are difficulties applying the shoulder loops.</p>	<ul style="list-style-type: none"> <li>• Place the patient in a more seated position in the bed or on the floor. Place a pillow behind the patient's neck/back, outside of the sling.</li> <li>• Raise the back rest of the bed, if possible.</li> <li>• Adjust the spreader bar to a more reclined position. Pay attention to the patient's head.</li> </ul>
<p>There are difficulties applying the leg loops.</p>	<ul style="list-style-type: none"> <li>• Bend the patient's legs or place a pillow under the patient's knees. If the sling is applied on the floor, place the patient's legs on the caregiver's knee.</li> <li>• Adjust the spreader bar to a more seated position.</li> </ul>
<p>It is difficult to use the sling in combination with the lift and spreader bar.</p>	<p>Use the sling only in accordance with the combinations described in section "<i>Allowed Combinations</i>" on page 7.</p>
<p>The position of the patient in the sling should be more reclined or seated.</p>	<p>Use the loop lengths to adjust the patient's position.</p>
<p>The patient wishes to have extra support for the neck.</p>	<p>Place a pillow behind the patient's neck, inside the sling.</p>

EN

# Technical Specifications

EN

<b>General</b>	
Safe working load (SWL) = Maximum total load	See "Allowed Combinations" on page 7.
Service life – Recommended period of use	2 years* * Refer to "Expected Service Life" on page 3.
Shelf life – Maximum period of storing new unpacked product	5 years
Model and Type	See "Sling Selection" on page 6.
<b>Operating, Transport and Storage Environment</b>	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport
<b>End of Life Disposal</b>	
Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.
<b>Parts and Accessories</b>	
Stiffeners	GX21270 and TBS110
Measuring tape	MTA1000

# Label on the Sling

EN

Care and Washing Symbols	
	Machine wash 70°C (158°F)
	No bleaching allowed
	No tumble drying allowed
	Tumble drying allowed
	No ironing allowed
	No dry cleaning allowed

Certificates/Markings	
	CE marking indicating conformity with European Community harmonised legislation
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Symbol for Loop	
	Use a loop spreader bar.

Article Number	
REF XXXXXX-X	Article number with -X refers to the sling size.
REF XXXXXXX	If the article number is stated without the size letter at the end, it means the sling is a One Size.

Fibre Content	
PES	Polyester
PU	Polyurethane

Misc. Symbols	
	Safe Working Load (SWL)
	Patient's name symbol
	Record symbol
	Read the <i>IFU</i> before use
	Article number
	Manufacturing date and year
	Manufacturer name and address

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**HONG KONG**

Arjo Hong Kong Limited  
Room 411-414, 4/F, Manhattan Centre,  
8 Kwai Cheong Road, Kwai Chung, N.T.,  
HONG KONG  
Tel: +852 2960 7600  
Fax: +852 2960 1711

**ITALIA**

Arjo Italia S.p.A.  
Via Giacomo Peroni 400-402  
IT-00131 ROMA  
Tel: +39 (0) 6 87426211  
Fax: +39 (0) 6 87426222  
E-mail: Italy.promo@arjo.com

**MIDDLE EAST**

Arjo Middle East FZ-LLC  
Office 908, 9th Floor,  
HQ Building, North Tower,  
Dubai Science Park,  
Al Barsha South  
P.O Box 11488, Dubai,  
United Arab Emirates  
Direct +971 487 48053  
Fax +971 487 48072  
Email: Info.ME@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges.

With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden

[www.arjo.com](http://www.arjo.com)

**arjo**

04.SB.00-INT1

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Arjo Inc.