# Domus 3

User's Manual

Care for a Healthy Life

# CONTENTS

<b>ALTERNATING PRESSURE REDISTRIBUTION SYSTEM</b>	
USER'S MANUAL	1
<b>SISTEMA DE REDISTRIBUCIÓN DE PRESIÓN ALTERNA</b>	
MANUAL DE INSTRUCCIONES	20
<b>SISTEMA DE REDISTRIBUIÇÃO DE PRESSÃO ALTERNADA</b>	
INSTRUÇÕES DE UTILIZAÇÃO	39
<b>SYSTÈME DE REDISTRIBUTION À PRESSION ALTERNÉE</b>	
INSTRUCTIONS D'UTILISATION	57
<b>ALTERNIERENDES DRUCKUMVERTEILUNGSSYSTEM</b>	
GEBRAUCHSANLEITUNG	75
<b>ALTERNEREND DRUKSYSTEEM</b>	
GEBRUIKERSHANDLEIDING	94
<b>SISTEMA DI RIDISTRIBUZIONE A PRESSIONE ALTERNATA</b>	
MANUALE DI ISTRUZIONI	112

MODEL NO.: 9P-074540 , 9P-074520

PLEASE READ THE MANUAL BEFORE USE

# IMPORTANT SAFEGUARDS

## READ ALL INSTRUCTIONS BEFORE USING

**DANGER** - To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store this product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

**WARNING** - To reduce the risk of burns, electrocution, fire, or injury to persons:

1. Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
2. The product may be used for patients with spinal injury, but suggested to consult with physician before use. However, it should not be used for patients with unstable spinal fractures.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Apex Medical Corp. for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
8. Never drop or insert any object into any opening or hose.
9. Do not modify this equipment without authorization of the manufacturer.
10. Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.

**CAUTION** -

1. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

**NOTE, CAUTION AND WARNING STATEMENTS:**

**NOTE** - Indicate some tips.

**CAUTION** - Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

**WARNING** - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

**SYMBOLS**

	Authorized representative in the European community.
	Manufacturer
	Complies with standards protecting against electric shock for type BF equipment.
	Attention, you should read the accompanying information carefully!
	Consult operating instructions for use
<b>IP21</b> 	Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops
	Class II
	Temperature limitation
	Dry clean, Any Solvent Except Trichloroethylene
	Do Not Iron
	Tumble Dry, Normal, Low Heat
	Do Not Tumble Dry
	Do Not Bleach
	Do Not Dry Clean
	Machine wash, regular / normal, 95 degrees C (203 degrees F)
	Machine wash, regular / normal, 60 degrees C (140 degrees F)
	Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.


# 1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

## 1.1 General Information

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure ulcers.

The system has been tested and successfully approved to the following standards:

### FOR US AND CANADA ONLY

E228589  
53DG

Medical Equipment-Air Pump  
WITH RESPECT TO ELECTRICAL SHOCK,  
FIRE AND MECHANICAL HAZARDS ONLY  
IN ACCORDANCE WITH ANSI/AAMI  
ES60601-1 (2005,3rd ed.) and CAN / CSA C22.2 No. 60601-1  
(2008, **3<sup>rd</sup> ed.**) and IEC 60601-1-11

### The below description for is for Canada only

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### EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

## 1.2 Intended Use

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It also provides following purposes:

- to help and reduce the incidence of pressure ulcers while optimizing patient comfort.
- for long term care of patients suffering from pressure ulcers.
- for pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure ulcer.


**NOTE:** Equipments are not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

## 2. Product Description

### 2.1 Pump & Mattress System

1. CPR
2. Pump unit
3. Mattress System
4. Quick Connector


### 2.2 Pump Unit

Front

1. Quick Connector Slot
2. Display Panel


Rear

1. Power Switch
2. Power Cord
3. Bed Frame Hangers
4. Air Filter

## 2.3 Front Panel

### 1. Pressure Adjustment Knob

The pressure adjustment knob controls the air pressure in the mattress. Turning the knob clockwise will increase the pressure; counter-clockwise decreases the pressure. Higher pressures will support heavier patients. The pressure should be adjusted according to individual comfort preferences.

When the mattress is finished pressurizing to the desired level of firmness, the Normal Pressure LED indicator will illuminate.

**NOTE:** You can also check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. The air cells will alternately inflate and deflate. You should feel slight contact with the buttocks when the air cells beneath the buttocks deflate.

**NOTE:** Every time the mattress is initialized (inflated), it is recommended that you set the pressure knob to "Max" to hasten inflation. You can then adjust the air mattress to the desired firmness later.

### 2. Normal Pressure Indicator

When the green LED illuminates, the pressure inside of air mattress has reached the desired pressure setting (preset on the Pressure Adjustment Knob).

### 3. Low Pressure Indicator

When the mattress is initializing, the Low Pressure LED will illuminate until the appropriate pressure is reached (according to the Pressure Adjustment Knob); this is normal. Otherwise, the Low Pressure LED is a warning, indicating that the pressure in mattress is unusually low. Check to ensure all connections are secured and correctly installed as per the installation instructions.

**NOTE:** If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

#### 4. Alternate / Static Switch

The ALTERNATE/STATIC switch selects between Alternate Pressure mode and Static Pressure mode. With Alternate Pressure mode, alternating air cells are partially deflated and inflated, avoiding prolonged pressure on any single point beneath the patient; this is to prevent pressure ulcers. With Static Pressure mode, all of the air cells are equally inflated.


#### 5. Alarm Mute

The pump unit is equipped with an audible alarm, the alarm will sound and the low pressure indicator will light up when air pressure is below normal. Caregiver may check for possible air leaks at this time. To mute the alarm sound, simply press the "Alarm Mute"  button on the panel. However, the low pressure indicator still remain on until the low-pressure faulty condition is solved.

**NOTE:** To re-activate the audible alarm function, please press "Alarm Mute"  button on the panel again.

### 3. Installation

Unpacking the box to inspect for any damage, which may have occurred during shipment. If there are any damages, please contact your dealer immediately.

1. Place the DOMUS 3 mattress on top of a foam or mattress bed. Secure the mattress firmly by fixing the straps to the bed frame if available.

**⚠ NOTE:** Domus 3 mattress is an overlay system, so there must be a foam or mattress underneath when using.

2. Hang the pump onto the footboard and adjust hangers so the pump is secured in an upright position; or place the pump on a flat surface.
3. Connect air hose connectors from air mattress to the pump unit. A "click" sound will be heard when the connection is completed and secured.

**⚠ NOTE:** Check and ensure the air hoses are not kinked or tucked under the mattress.

4. Plug the power cord into a electrical outlet with the appropriate voltage.

**⚠ NOTE:** The pump unit must be suitable for the local power voltage.

5. Then turn the main power switch to ON position.

**NOTE:** You can unplug the unit to disconnect the device. Do not position the equipment so that it is difficult to operate the disconnecting device.

**CAUTION:** The pump supplied with the mattress can only be used for mattresses recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)


**NOTE:** After installation, make sure any excess cord or tubing is removed from possible foot traffic to avoid accidental tripping. All EQUIPMENT should be positioned to always allow unhindered patient access by physicians and caregivers.

6. After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents.

The EQUIPMENT should be firmly placed at position where users/doctors can access easily.

## 4. Operation


### 4.1 General operation

1. Switch on the main power switch found on the side of the pump.
2. Once the pump begins delivering air into the mattress it will take approximately 30 minutes to fully inflate the mattress.  
Every time the mattress is set up for use, it is recommended that the pressure first to be set to the max. The user / career can then adjust air mattress weight levels to the desired softness after set up has been completed.
3. The low pressure indicator (yellow LED) will light up when the mattress is not fully inflated.
4. When the appropriate pressure is reached, the low pressure indicator (yellow LED) will go off and normal pressure indicator (green LED) will light up automatically,
5. According to the weight and height of the patient, adjust the knob setting to the most comfortable level without bottoming out, then the pressure in mattress will slowly increase to a certain value after the air mattress is ready to use

**⚠ NOTE:** Check to see if the suitable pressure is selected by sliding one hand between the deflated air cells and the patient to feel patient's buttocks. Users should be able to feel the minimum contact.


## 4.2 Emergency CPR Operations

When there is a need to perform emergency CPR on the patient, pull the quick release CPR tag located at the head-end of the mattress on the right hand side.

The quick connector found from the pump unit can also be disconnected to hasten deflation. Be sure to reconnect the quick connector to the pump unit once restore the power supply.

**⚠ NOTE:** When there is a power outage or a need of transporting the patient, one can cap on the quick connector as shown on the right to slow down the air leakage. Depending on the setting and mattress height, the mattress can remain inflated for a period of time. If there is a bottom-out issue, please remove the patient from deflating mattress to another proper surface.

## 5. Cleaning

It is important to follow the cleaning procedures before first-time use with any patient; this is to avoid the possibility of transmission of pathogens between patients and caregivers.

Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit.

 **CAUTION:** Do not immerse or soak the pump unit in liquids.

Wipe-down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent. Avoid dust and proximity to dusty areas. The cover may also be cleaned using sodium hypochlorite diluted in water. All components should be air dried thoroughly before use.

 **CAUTION:** Do not use phenolic based products for cleaning.

 **CAUTION:** Do not dry the mattress in direct sunlight.

The carrying bag (if available) should be turned inside-out and completely wiped-down using disinfectant solutions. Allow it to air-dry thoroughly and then turn it back and wipe down the outside of the bag with disinfectant solutions.

Cover Material: Nylon/PU	    
Cover Material: Silver+ Nylon/PU	
Cover Material: Stretch	    

## 6. Storage

1. Lay the mattress on a flat surface, upside-down.
2. Roll-up the mattress starting from the head until you reach the foot of the mattress.
3. Use the straps at the end of the mattress and stretch them around the mattress to prevent it from unraveling.


 **NOTE:** Do not fold, crease or stack the mattresses.

## 7. Maintenance

### 7.1 General

1. Check the main power cord and plug for abrasions or excess wear.
2. Check the mattress cover for signs of wear or damage.
3. Disconnect the air tube from the mattress. Then check the airflow coming from the two air outlets on the pump. They should be alternately delivering air when the pump is set at "alternate" mode.
4. Check the air hoses for kinks or breaks. For replacement, please contact your local distributor.

### 7.2 Air Filter Replacement

1. Open the air filter cap located at the back of pump.
2. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
3. Check and replace air filter regularly if environment is dirty.

## 8. Expected Service Life:

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

## 9. Trouble Shooting

### **Q1 Power is not ON**

- Check if the plug is connected to mains.

### **Q2 Alarm is on (Visual)**

- Check if the power is suddenly shut down.
- Check if the CPR is sealed.
- Check if the connection between air tube connector to pump unit is tightly secured.
- Check if all tubing connections along mattress are secured.

### **Q3 Patient is bottoming out**

- Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort.

### **Q5 Mattress form is loose**

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

### **Q6 No air produced from some air outlets of the air tube connector**

- This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

If the above information does not solve your problems, please contact your local agent directly. They might require a technician to take care the problem.

## 10. Technical Specification:

Item		Specification		
Power Supply (Note: See rating label on the product)		AC 100-120V, 60Hz, 0.16A (For 120V System) AC 220-240V, 50 Hz, 0.06 A (For 230V System)		
Fuse Rating		T1A, 250V (120V); T1AL, 250V (230V)		
Cycle time		Fixed 8 min @ 60Hz; 9.6 min @ 50Hz		
Dimension (L x W x H)		29 x 18.5 x 12.6 cm or 11.4" x 7.3" x 5"		
Weight		2.2 Kg or 5 lb		
Environment	Atmospheric Pressure	700 hPa to 1013.25 hPa		
	Temperature	Operation: 10° C to 40° C (50° F to 104° F) Storage: -15° C to 50° C (5° F to 122° F) Shipping: -15° C to 70° C (5° F to 158° F)		
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10% to 90% non-condensing		
Classification		Class II, Type BF, IPX0 (120V); IP21 (230V) Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or AP6 protection)		
Mattress		Specification		
Model		4" Overlay	5" Overlay	5" Overlay +3" Foam pocket
Dimension (L x W x H)		190 x 80 x 10.2 cm 74.8" x 31.5" x 4"	200 x 85/90 x 12.7 cm 78.7" x 35.4" x 5"	200 x 85/90 x 12.7 cm 78.7" x 34.5/35.4" x 5"
Weight		3.1 kg or 6.9 lb	5.2 kg or 11.4 lb	5.6 kg or 12.5 lb

### NOTE:

1. Consult the distributor or EU representative for further technical documents.
2. These specifications are also applicable for other regions operating with the same power supply.
3. Mattress dimensions and weight is measured without foam.
4. The manufacturer reserves the right to modify the specifications without notice.

## Appendix A: Emc Information

### Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

#### Warning:

1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		120 / 230 V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.

<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz  385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz  385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10V/m</p>	<p>Recommended separation distance  <math>d = \sqrt{P}</math> 150kHz to 80MHz  <math>d = 0.6\sqrt{P}</math> 80MHz to 800MHz  <math>d = 1.2\sqrt{P}</math> 800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
---	--	---	--------------	---

NOTE 1:  $U_T$  is the a.c. mains voltage prior to the application of the test level

NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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
**⚠REMARQUE:** When there is a power outage or a need of transporting the patient, one can capped on the quick connector as shown on the right to slow down the air leakage. Depending on the setting and mattress height, the mattress can remain inflated for a period of time. If there is a bottom-out issue, please remove the patient from deflating mattress to another proper surface.

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Print-2018/All rights reserved
**Apex Medical (Kunshan) Corp.**