

The **NYOrtho CushionGuard Wheelchair Footrest Sleeve** offers an essential combination of comfort, safety, and hygiene for wheelchair users. Its foam-filled, non-slip design ensures user comfort and stability, while the antimicrobial fabric promotes cleanliness and health. Easy to apply, clean and durable, this footrest cover also allows personal customization, enhancing the wheelchair experience.

- **Comfort:** Foam-filled NYOrtho footrest cover provides substantial user comfort, reducing fatigue from extended periods in a wheelchair.
- **Safety:** The non-slip vinyl material prevents the user's feet from sliding off, enhancing safety and overall wheelchair operation. Meets CA 117-2013 standards.
- **Hygiene and Health:** The footrest cover is constructed with antimicrobial fabric, inhibiting the growth of bacteria, fungi, and other microbes, reducing the risk of infections and odors, and promoting a healthier environment.
- **Easy to Clean:** The vinyl cover can be quickly wiped clean, offering a practical and efficient method for maintaining cleanliness.
- **Customization:** The addition of a footrest cover allows users to personalize their wheelchair to better suit their needs and preferences.
- **Durability:** Made from durable vinyl, the NYOrtho no-slip footrest cover is designed to withstand daily wear and tear, thereby prolonging the life of the footrests.



INSTRUCTIONS

- Ensure the wheelchair is locked in place before attempting to attach the footrest cover.
- Align the NYOrtho No-Slip Footrest Cover with the wheelchair footrest, ensuring the foam side faces upward.
- Secure the cover onto the footrest using the hook and loop strap. Make sure the cover is firmly in place to avoid movement during use.
- Check the footrest cover before each use. If any signs of wear and tear are noticed, replace the cover immediately to ensure continued safety and comfort.
- The NYOrtho No-Slip Footrest Cover is designed to keep the patient's feet in position. Ensure the user's feet are comfortably placed on the cover during use.

Dimensions: 6.75x6.75 in.
UOM: PAIR

CLEANING INSTRUCTIONS

- Wipe Clean with mild detergent. **DO NOT** use quaternary/isopropyl based disinfectants on vinyls.



MODEL	UOM	SIZE
9542	PAIR	6.75 x 6.75 in.