
- Use a screwdriver or pen to slide the 10° stop out of both hinges. Reapply the desired extension stop and gently secure with the remaining screws (not over-tighten) (Fig. 5).
- Follow the three-step procedure for the addition of a flexion stop on the other side of each hinge.

Device Removal
Undo the three-step procedure. Grip the device at the level of the calf and pull it off in a downward direction or unwrap the device.

Accessories and Replacement Parts
Please refer to the Ossur catalog for a list of available replacement parts or accessories.

Cleaning and care

- Hand wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

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Note: Do not remove the yellow pull tab until after the initial fitting.

1. Entfernen Sie die Gurte.
2. Abhängend von der Version der device chosen, with the leg in full extension, either pull the device over the knee using the pull-tabs (sleeve version) or wrap the device around the knee (lower portion first) and secure the contact closures on the front of the knee (wrap version).

Position the device, such that the kneecap opening and the middle of the hinge are centered on the knee (Fig. 3).

Fasten the straps at the calf and thigh in the order shown by inserting them through the corresponding D-ring (Fig. 2). Straps should wrap inside the pull-tabs at the top of the device (Fig. 3).

Note: When the final fit is confirmed, remove the yellow pull tab to secure the strap in position (Fig. 3).
Note: Universal buttress accessory can be added to the inside of the sleeve (sleeve version) and adjusted to multiple positions to optimize patella support (Fig. 6).

Device Adjustments
The device comes with 10° extension stops installed and includes the following additional stops:

- Extension stops are labeled 0°, 10°, 20°, 30° and 40°.
 - Flexion stops are labeled 0°, 45°, 60°, 75°, and 90°.
- Aucune contre-indication connue.

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DISPOSAL
The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY
Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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