

FLEX-SUPPORT No. 326 HARRIS HEMI SLING

APPLICATION INSTRUCTIONS (SHOWN FOR LEFT ARM)



Place the elbow cuff to the affected side, wrist cuff to the unaffected side. Lengthen or shorten front and back straps for proper fit and support.



1) Place the square opening over the head with elbow cuff at affected side.



2) Put unaffected arm through the cuff and adjust straps behind shoulders.



3) Lift affected arm and place elbow in cuff and adjust.



4) Insert unaffected hand through cuff, pick up affected hand and place in cuff. Apply Velcro strap to maintain cuff location between elbow and wrist.



Check for proper fit given by therapist.