

INSTALLATION

For use on any size bed frame.

1. Loop bed ladder strap around the bed frame at the foot of the bed.
2. Snap the buckle to secure the strap and adjust the length to ensure the first handle meets waist lying down. Adjust as necessary.



HOW TO USE

Using a hand-overhand movement, sit up in bed by pulling on successive rungs of the ladder until reaching a seated position.

CLEANING

Wipe bed ladder strap and handles with a damp cloth and dry immediately.

SPECIFICATIONS

- Strap Length: Adjustable 48" – 72"
- Length between handles: 10.5"
- Handle Grip: 5"
- Weight Capacity: 300 lbs

WARNING

- Periodically inspect the bed ladder for signs of wear. Do not use the bed ladder if it appears weakened.
- Be sure the bed ladder is securely attached to the bed, adjusting as necessary.