





Fitting Protocol

Multi Sling
Model 25105



- | | |
|--|--|
| Strap A (Front)
<input type="radio"/> Red Strap
<input type="radio"/> Yellow Strap
<input type="radio"/> Green Strap | Strap B (Back)
<input type="radio"/> Red Strap
<input type="radio"/> Yellow Strap
<input type="radio"/> Green Strap |
| Strap C (Back)
<input type="radio"/> Red Strap
<input type="radio"/> Yellow Strap
<input type="radio"/> Green Strap | Strap D (Front)
<input type="radio"/> Red Strap
<input type="radio"/> Yellow Strap
<input type="radio"/> Green Strap |

Components of Multi Sling:




Choose the loop combination that will give the best height for the intended use. Match colours on both sides of the patient (2 front straps and 2 back straps of chest belt).



When fitting the chest belt portion of the Multi Sling, ensure that it rests on the lumbar curve above the coccyx.

If applying a waist belt, insert through the inner loops of the chest belt.



The Multi Sling should be applied while the patient is in the seated position.

Maximum User Weight: 150kg/330.7lb

Model 25105 - Polyester

Transfer Situation:

- | | | | |
|-----------------------|-----------------------|----------------------------|-----------------------|
| <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| To/From a Bed | To/From a Chair | Gait Training/Mobilization | Dressing/Undressing |

Available sizes:

Indicate size required

- S
 M
 L
 XL

Leg Placement:

- | | |
|--|---|
| <input type="radio"/>
Not applicable for this sling
Divided relaxed position | <input type="radio"/>
Not applicable for this sling
Open position |
| <input type="radio"/>
Not applicable for this sling
Divided position | <input type="radio"/>
Not applicable for this sling
Closed position |

Suspension:

- | | |
|----------------------------------|---|
| <input type="radio"/>
2-Point | <input type="radio"/>
4-Point |
| <input type="radio"/>
3-Point | <input type="radio"/>
Not available for this sling |

Accessories required:

- 50866 Upholstered Cover for Hanger Bar. **Do not use in wet environments.**
 - 25150 Neck Cushion
 - 25145 Waist Belt - Velcro (Inserted through inside loops of chest belt)
 - 25146 Waist Belt - Safety Lock (Inserted through inside loops of chest belt)
 - 25107 Hip Belt
- Waist Belt (Select Size)

Available Sizes:

 - 60 (60cm/23.6")
 - 70 (70cm/27.6")
 - 80 (80cm/31.5")
 - 90 (90cm/34.4")
 - 100¹ (100cm/39.4")
 - 110 (110cm/43.3")
 - 120 (120cm/47.3")
 - 130 (130cm/51.2")

¹Model 25146 only

Patient Information:

Name		
Facility	Room	
Lift Used	Hanger Bar Used	

The Multi Sling is designed for gait training and rehabilitation exercises. It can be used in many situations and together with a wide range of accessories can act as support during standing, walking and balancing exercises as well as for visits to the toilet or dressing and undressing. The Multi Sling is suitable for patients with some body stability and partial use of their legs and is approved for 2-point, 3-point or 4-point suspension. When used with a hip belt, it provides excellent support for patients with serious functional impairments such as stroke. Adding a groin support provides additional support and safety.

Additional Instructions

To staff/caregiver:

To patient:

Fitting Protocol Completed by:

Name		
Date	Position/Title	



Important! This fitting protocol does not replace the user manual. Read all instructions before use!