GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

USAGE

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

The device comes with a pre-installed Buttress (A), which can be removed or adjusted if necessary.

1. Undo the straps (**Fig. 2**).
2. With the leg in full extension, grasp the Pull Tabs (B) and pull the device over the knee (**Fig. 3**). Pull evenly over the knee joint until the Kneecap Opening (C) encloses the kneecap and the kneecap is free from pressure.
3. Fasten the Lower Strap (D) and then the Upper Strap (E) (**Fig. 4**).
4. Stand up and verify the correct position and fit of the device. Adjust the device positioning and strap pressure, if needed.

Device Removal

Undo the straps. Grip the device at the level of the calf and pull it off in a downward direction.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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MD Medical Device

INTENDED USE

The device is intended for external support or stabilization of the knee.

Indications for use

For patients suffering from knee conditions that may benefit from additional knee and patella support, such as: patellar tendonitis, subluxation dislocation, post-operation rehabilitation, post arthroscopic surgery, chondromalacia, bursitis, arthritis and tendonitis. No known contraindications.

Warnings and Cautions:

- Inspect the skin under the device's contact area before and after each use.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.