ENGLISH

| | |
|-----------|----------------|
| MD | Medical Device |
|-----------|----------------|

INTENDED USE

The device is intended for external support or stabilization of the knee.

Indications for use

For conditions that may benefit from realignment of the patella, such as:

- Lateral patellar subluxations and dislocations
- Lateral patellofemoral malalignment
- Patellofemoral maltracking and pain
- Chondromalacia due to lateral maltracking

No known contraindications.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities. The device should be used directly on intact skin and never over an open wound. The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple uses.

USAGE

Size Selection

Measure the calf and thigh circumference 15 cm (6”) above and below mid-patella. The sizing chart is on the packaging, if between sizes always select the smaller size for optimal fit.

Device Application

Note: For best results when contouring hinge arms, remove hinges and bend hinge arms as desired (**Fig. 1**). Re-apply hinges before application, ensuring that the hinge arms are completely seated and secured within the hinge pockets.

- Undo the straps. With the leg in full extension, grasp the edge and pull the device over the knee. Pull evenly over the knee joint until the lateral j buttress encloses the knee cap and the knee cap is free from pressure (**Fig. 2**).
- Affix straps in the order shown on the fastening locations to the inside of the knee (**Fig. 3**).
- Stand up and verify the correct position and fit of the device (**Fig. 4**). Adjust the device positioning and strap pressure, if needed.

Device removal

Undo the straps. Grip the device at the level of the calf and pull it off in a downward direction.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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