

USER MANUAL

OAKWORKS® Treatment Tables



Powerline Series



Seychelle™ Series



OAKWORKS®

Notice

The information contained within this document is subject to change without notice and should not be construed as a commitment by OAKWORKS[®], Inc.

OAKWORKS[®], Inc. encourages requests for technical specifications and the like documentation to ensure accuracy. The appropriate documentation is available upon request.

OAKWORKS[®], Inc. shall not be liable for incidental or consequential damages in connection with or arising out of the furnishing, performance, or use of this document and the program material which it describes.

Printed in U.S.A.

All rights are reserved. No part of this document may be photocopied, reproduced or translated to another language without prior written consent of OAKWORKS[®], Inc.

OAKWORKS[®] is a registered trademark of OAKWORKS[®], Inc.

TABLE OF CONTENTS

- Product Use Description 1
- Important Safety Instructions 1
- Product Description & Photo
 - Powerline Series..... 2-3
 - Seychelle Series..... 4-5
- Installation
 - Table Assembly 6
 - Inset Shelf 6
- Directions for Use
 - Paper Roll Holder Use..... 7
 - Wave Backrest Use 7
 - Table Extender 7
 - QuickLock Face Rest 8
- Cleaning & Disinfection
 - Recommended Cleaners/Disinfectants 9
 - Cleaning Process..... 9
- Inspections 10
- Warranty Information 10
- Model Number & Serial Number 11
- Specifications
 - Product Specifications 12
 - Environmental Conditions 12
- Contact Information back cover

PRODUCT USE DESCRIPTION / IMPORTANT SAFETY INSTRUCTIONS

PRODUCT USE DESCRIPTION

The Oakworks Powerline & Seychelle™ Tables are used to support the patient during examinations and procedures. They are intended to be operated by a professional in a exam/treatment environment. No special training is required but a review of the following manual is important for the safety of the operator and patient. The professional should read and understand this entire manual before use with a patient. There are no known contraindications to the use of this equipment.

IMPORTANT SAFETY INSTRUCTIONS



CAUTION

READ AND SAVE THESE INSTRUCTIONS



This symbol, when used in this manual and on product labels, represents a caution warning. Be sure to read and comply with all precautions and warnings.

The Backrest section is not designed to support the entire weight of the patient. Do not sit on the Backrest section.

The table is designed to be a stand alone table. This table must not be modified or incorporated into any other equipment.

There are potential pinch points around and underneath the tabletop. It is the responsibility of the operator of this equipment to insure that bystanders are not in the area below or around this equipment during operation.

Proper operation of this equipment is very important for the safety of the operator, patient, and any other individuals in the area of this equipment. Directions for use of this equipment are described in this manual. The operator should read these sections carefully.

Do not use outdoors.

Use this furnishing only for its intended use as described in this manual. Do not use attachments not recommended by Oakworks.

Weight Limit: (patient and accessories) 500 lbs./227 kg. Do not exceed.

PRODUCT DESCRIPTION

Powerline Flat Top



| STANDARD SPECIFICATIONS | |
|---------------------------------------|---|
| Height Range | Fixed 30" (76 cm.) |
| Table Capacity | 500 lbs. (226 kg.) (patient and accessories) |
| Widths | 27" (69 cm.) or 30" (76 cm.) |
| Length | 72" (183 cm.) |
| Padding | 2.5" (6 cm.) Semi-Firm |
| Upholstery | Terra Touch™ |
| Paper Roll Holder | Accommodates up to 21" (53 cm.) rolls (roll not included) |
| Table Weight* | 90 lbs. (41 kg.) |
| Warranty | 3 year parts, 1 year labor |
| Safety Listings | CE marked, FDA |
| OPTIONS & ACCESORIES | |
| QuickLock Face Rest | |
| Table Extender | |
| Inset Shelf (adds 25 lbs. (11 kg.)) | |
| Padding - 1.75" (4 cm.) Firm Response | |

PRODUCT DESCRIPTION


Powerline Wave Backrest




QuickLock Face Rest (accessory)



Table Extender (accessory)



Wave Backrest



| STANDARD SPECIFICATIONS | |
|---------------------------------------|---|
| Table Height | Fixed 30" (76 cm.) |
| Wave Backrest | Manual 0-60°, 6 set positions |
| Table Capacity | 500 lbs. (226 kg.) (patient and accessories) |
| Widths | 27" (69 cm.) or 30" (76 cm.) |
| Length | 72" (183 cm.) |
| Padding | 2.5" (6 cm.) Semi-Firm |
| Upholstery | Terra Touch™ |
| Paper Roll Holder | Accommodates up to 21" (53 cm.) rolls (roll not included) |
| Table Weight* | 90 lbs. (41 kg.) |
| Warranty | 3 year parts, 1 year labor |
| Safety Listings | CE marked, FDA |
| OPTIONS & ACCESORIES | |
| QuickLock Face Rest | |
| Table Extender | |
| Inset Shelf (adds 25 lbs. (11 kg.)) | |
| Padding - 1.75" (4 cm.) Firm Response | |

PRODUCT DESCRIPTION

Seychelle Flat Top



| STANDARD SPECIFICATIONS | |
|-------------------------------------|--|
| Height Range | Fixed 30" (76 cm.) |
| Table Capacity | 500 lbs. (226 kg.) (patient and accessories) |
| Widths | 27" (69 cm.) or 30" (76 cm.) |
| Length | 72" (183 cm.) |
| Padding | 3" (8 cm.) Comfort |
| Upholstery | Terra Touch™ |
| Table Weight* | 90 lbs. (41 kg.) |
| Warranty | 3 year parts, 1 year labor |
| Safety Listings | CE marked |
| OPTIONS & ACCESORIES | |
| QuickLock Face Rest | |
| Table Extender | |
| Inset Shelf (adds 25 lbs. (11 kg.)) | |
| Padding - 4" (10 cm.) Comfort | |

PRODUCT DESCRIPTION


Seychelle Wave Backrest




QuickLock Face Rest (accessory)



Table Extender (accessory)



Wave Backrest



| STANDARD SPECIFICATIONS | |
|-------------------------------------|--|
| Table Height | Fixed 30" (76 cm.) |
| Wave Backrest | Manual 0-60°, 6 set positions |
| Table Capacity | 500 lbs. (226 kg.) (patient and accessories) |
| Widths | 27" (69 cm.) or 30" (76 cm.) |
| Length | 72" (183 cm.) |
| Padding | 3" (8 cm.) Comfort |
| Upholstery | Terra Touch™ |
| Table Weight* | 90 lbs. (41 kg.) |
| Warranty | 3 year parts, 1 year labor |
| Safety Listings | CE marked |
| OPTIONS & ACCESORIES | |
| QuickLock Face Rest | |
| Table Extender | |
| Inset Shelf (adds 25 lbs. (11 kg.)) | |
| Padding - 4" (10 cm.) Comfort | |

INSTALLATION

TABLE ASSEMBLY

Tools Required:

9/16" / 14mm Wrench/Socket

SUPPLIED PARTS:

TOOLS


- 1 - Large Allen Wrench
- 1 - Small Allen Wrench

LEG HARDWARE

- 8 - Reinforcement Plates
- 16 - Spring Washers
- 4 - Long Large Bolts
- 4 - Lock Nuts



LEG BRACE

- 1 - Wood Leg Brace
- 2 - Short Bolts
- 2 - Barrel Nuts




The tables are shipped complete in one box. Follow the step below to fully assemble the table:

1. Remove the table from the box and place it down on the floor with the top down. Remove all hardware and the Leg Brace.
2. Unfold the legs to a vertical position and remove any protective wrappings.
3. Install the long large bolts, reinforcement plates, spring washers and lock nuts as shown, two bolts on each leg. Using the large allen wrench and 9/16" / 14mm wrench socket, tighten every bolt/nut connection -- 8 total -- to 12 ft-lbs minimum or until the spring washers are flat.
4. Install the barrel nut into the leg brace being sure the screw slit is horizontal. Install the short bolt through the leg rail. Start the bolt into the barrel nut using your fingers first. After the bolt has started, tighten the bolt using the small allen wrench. Make sure the leg brace is flush with the adjoining rails while tightening bolts. Repeat on opposite end of the leg brace.
5. Turn your table right side up. Get up on the table and rock back and forth to check that all bolts are tight. The table should be very stiff. Dependent on the environmental conditions of your table, it may be necessary to periodically re-tighten the hardware.



SHELF ATTACHMENT

If you have ordered a shelf, drop the shelf into place with the finished side up. The shelf has grooves at both ends that enable it to be held firmly in place by the legs. If necessary, the shelf can be further stabilized by installing the enclosed Velcro™ strips. Place strips on each corner of the shelf's bottom side, and on each end of the rails. Make sure not to mismatch Velcro™ strips; place all "hook" or all "loop" pieces on either the rail or shelf.




DIRECTIONS FOR USE

PAPER ROLL HOLDER USE




After the table is set-up, locate the hole in the leg.



Using the supplied screw & allen wrench, secure the paper roll holder bracket into place. Repeat step for the opposite side.



Locate the wooden dowel & feed it thru your paper roll. Place the wooden dowel in the groove located in each of the paper roll holder brackets.



To secure the wooden dowel, feed it into the pocket of the groove on both sides of the table.

WAVE BACKREST USE




CAUTION

Do not adjust the Backrest with a patient on the table.



TO RAISE THE BACKREST

Lift the backrest to the desired height and release. The Backrest support arm will lock in the nearest slot in the side rail.



TO LOWER THE BACKREST

Lift the backrest slightly, lift and hold the support arm out of the slot and lower the backrest to the desired position.


TABLE EXTENDER




CAUTION

Do not put excessive weight or pressure on the table extender.

The Table Extender lengthens the table for your tallest clients and can improve access to your patient.



Locate the Accessory attachment outlets.



Unfold the dowels and insert into the outlets.

DIRECTIONS FOR USE

QUICKLOCK FACE REST



* U.S. Patent #5,177,823



CAUTION

Do not put excessive weight or pressure on the Quicklock Face Rest. Do not extend the dowels of the Quicklock Face Rest 3" (7.6 cm.) from the table top. The Maximum Weight Capacity for the QuickLock Face Rest is 25 lbs. (11 kg.)

Double articulating action provides nearly limitless positioning options with the QuickLock™ Face Rest. The adjustable pad takes pressure off sensitive facial areas while providing superb support. The QuickLock™ Face Rest can be used in both prone and supine positioning and folds out of the way when not in use. Adjustments are made easily by loosening the cam locks, adjusting for height and angle and then locking the cam locks into the desired position.



Important Safety Note:
Maximum Distance from edge
of Table Top= 3" (7.6 cm.)



Insert dowels into
Accessory attachment out-
lets.



Loosen Cam Locks &
adjust face rest into place.



Tighten the Cam Locks.



The Face Rest can be used
for either prone or supine
positioning.

CLEANING & DISINFECTION

RECOMMENDED CLEANERS/DISINFECTANTS



DANGER

To reduce the risk of electric shock:

- Always unplug this furnishing from the electrical outlet before cleaning.

The following cleaners were tested by OAKWORKS® and found to be acceptably compatible with all materials used in construction of our tables:

1. Formula 409®
2. Fantastik®
3. Green Windex®

The following disinfectants were tested by OAKWORKS® and found to be acceptably compatible with all materials used in construction of our tables:

1. Parker Labs: Protex™ (EPA reg# 6836-152-82613)
2. Healthlink® E-Z Kill® (EPA reg# 59894-10-67454)
3. MADA Medical: Madacide-FD (EPA reg# 1130-15-11703)

The following disinfectants were tested by OAKWORKS® and found to be acceptably compatible with most materials, with some cautions as indicated:

1. 10% Sodium Hypochlorite (bleach) solution in water (EPA reg# 5813-100) (possible mild corrosion of plated metals)
2. Metrex Cavicide One (EPA reg# 46781-12) (possible slight deterioration of upholstery and wood lacquer)

OAKWORKS® recommends a prepackaged wipe for cleaners/disinfectants to ensure best distribution of disinfectant for the required kill time, without leaving excess residue and/or overexposing components therefore minimizing the potential for damage to materials. Please read and follow disinfectants manufacturers' directions for cleaning and disinfection.

The following disinfectants were tested by OAKWORKS® and found to be incompatible with some materials and are therefore NOT recommended for use:

1. Clorox Broad Spectrum Quaternary (EPA reg# 67619-20)(deterioration of upholstery, wood lacquer, plated metals)
2. Clorox Hydrogen Peroxide Disinfectant (EPA reg# 67619-24)(corrosion of metals, deterioration of upholstery and wood lacquer, leaves a hard-to-remove residue)
3. Virox Accel TB (EPA reg# 74559-1)(corrosion of metals, deterioration of upholstery and wood lacquer, leaves a hard-to-remove residue)
4. Clorox Healthcare Bleach Germicidal Cleaner (EPA reg# 56392-7)(deteriorates upholstery, wood lacquer, mild corrosion of metals, dissolves foam in pads)

OAKWORKS® does NOT recommend the use of cleaners/disinfectants containing Hydrogen Peroxide, Acetic Acid, or Phenolics. These chemicals can cause damage to the appearance and/or material integrity of various components. Also, while the recommended cleaners/disinfectants list includes products containing Quaternary Ammonium compounds ("quats"), not all products containing quats are approved for use. Some contain additional detergents and/or surfactants which can be detrimental to some materials.

A note on Bleach: While a 10% sodium hypochlorite (household bleach) solution (EPA reg# 5813-100 or equivalent) can be an effective disinfectant and is dilute enough to be benign to most materials, it alone is not an effective cleaner and a separate product must be used for the initial cleaning steps of the procedure. Because of possible chemical incompatibilities between various cleaning products and bleach, utmost care must be taken by the user to avoid potential exposure to harmful or toxic by-products of the combination. Also, because bleach leaves a potentially corrosive residue as it evaporates, it must be rinsed with clean water after disinfection. It is therefore NOT recommended by OAKWORKS®. Please note, one step cleaner/disinfectants using bleach and a detergent/ surfactant were found to damage materials and are NOT recommended.

CLEANING PROCESS



DANGER

To reduce the risk of electric shock:

- Always unplug this furnishing from the electrical outlet before cleaning.

Follow the cleaners/disinfectant manufacturers' directions for use. Please note that cleaning and disinfecting an OAKWORKS® table is a two part process. First it must be cleaned of any visible soil, then it can be disinfected. The table should be cleaned as soon as possible after use. Please follow these procedures for best results:

1. Using a fresh wipe or wipes, clean any visible soil off of the table, working from the top to the bottom of the table.
2. Throw the used wipe(s) away.
3. Using a fresh wipe or wipes, thoroughly wipe all surfaces of the tabletop and any high-contact areas such as handles, hand-sets, etc., making sure they remain wet for the disinfectant manufacturer's recommended contact time.
4. Throw the used wipe(s) away and allow the table to dry.

INSPECTIONS

RECOMMENDED REGULAR INSPECTIONS (monthly or local standard)

- Check for damage to the power, hand control(s) and foot control(s) on all cables.
- Visually inspect components for obvious damage that could cause problems during operation.

RECOMMENDED PERIODIC INSPECTIONS (yearly or local standard)


- Check for damage to the power, hands control and foot control cables and all visible wiring.
- Visually inspect components for obvious damage that could cause problems during operation.
- Check all mechanical functions using the hand control. Repeat using the foot control. Check for abnormal noises.
- Check that all fasteners are present and fastened securely.
- Check table grounding.
- Clean unusual buildup of dirt on the table and/or parts of the table not normally cleaned on a regular basis.
- Check for tears or cracks in the upholstery.

MODEL NUMBER & SERIAL NUMBER

MODEL & SERIAL NUMBERS



The model number and serial number are located under the top of the table.




Model Number

Serial Number

MODEL NUMBER DESCRIPTION

Model Numbers and Serial Numbers always start with a letter.

PWMPBRXXXXXXXXFRTTXX XX



Options

Not Used

Fabric: Terra Touch

Padding

FR=1.75" (4 cm.) Firm Response

SF=2.5" (6 cm.) Semi-Firm

PL=3" (8 cm.) Comfort

AC=4" (10 cm.) Comfort

Table Top Length

Minimum Table Height

Table Width

Top

FT=Flat Top

BR=Wave Backrest

Base

Table

PW=Powerline Series

SE=Seychelle Series

SPECIFICATIONS

| PRODUCT SPECIFICATIONS | |
|------------------------|--------------------------|
| Weight | 90-115 lbs. (41-52) kg. |
| Shipping Weight | 105-130 lbs. (48-59 kg.) |
| Weight Limit | 500 lbs. (226 kg.) |

| ENVIRONMENTAL CONDITIONS | | | |
|--------------------------|-----------------------------|---------------|----------------------|
| Conditions | Temperature | Humidity | Atmospheric Pressure |
| Normal Use | 50° (10°C) to 104° (40°C) | 20% to 60% RH | 98 to 105 kPa |
| Storage & Transport | -20° (-29°C) to 135° (57°C) | 20% to 95% RH | 98 to 105 kPa |

USER MANUAL

OAKWORKS® Treatment Tables



EMERGO EUROPE



Manual Part Number MMMNST0005-EN
Revision level: 2
Revision Date: 07/30/2018

English, Printed in USA



OAKWORKS®