

STAND ASSIST TRANSFER DISC

LVA3041BLK

vive




OWNER'S MANUAL

STAND ASSIST TRANSFER DISC

WHAT'S INCLUDED

- (A) Turning Plate
- (B) Bottom Pole with Knee Rests
- (C) Top Pole with Handles
- (D) Pole Locking Knob
- (E) Turning Plate Screws x3



INSTRUCTIONS

1. Take the (B) Bottom Pole with Knee Rests and insert it into the (A) Turning Plate so that all three screw holes all line up. Note that the pole needs to rise up just slightly from the floor of the base for the holes to line up.



2. Take the Base Plate Screws (3) and screw them into the holes on the base. Better to start with the high screw hole on the rear to balance the device in the middle and then slightly screw in the two side screws so they catch the thread.

Once all three screws are slightly threaded then tighten all three screws all the way down.

3. Insert the (C) Top Pole with Handles into the Bottom Pole.

Use the (D) Pole Locking Knob to securely fasten it into place.



Note that if you decide to break down your unit that there is a locking button on the read of the pole.



4. You can adjust the height of the leg rests up and down using the locking knob to loosen and then re-tighten the clamp around the pole where needed.



USE INSTRUCTIONS FOR THE CAREGIVER:

Place the device directly in front of the patient within one foot of the patient so their feet fully resting on the Turning Plate about 4 inches apart and that their knees are comfortably pressed against the knee rests.

Place your foot on the Turning Plate in between the two wheels (You will notice a raised slope with a V on it, Place your foot there). Shift your weight onto this foot so the device is firmly held down. Place your hands on the top of the handles.

Have the patient pull themselves up using the bottom handles or even the bottom part of the top oval-ish shaped handle (whichever is more comfortable or gives your patient the easiest pulling motion for their height and particular strength level). Again make sure you are using your weight and grip on the top handle to stabilize the device while your patient is making this motion.

Once your patient is standing and comfortably resting their knees against the knee rests you can turn them any degree to the left or the right to place them in the best position to seat them in the new device you are moving them to whether that is a bed, wheelchair, sitting chair, commode, etc.

Follow the same procedure in reverse as the sit back down in this new device again making sure you as the caregiver have both hands on the top bar and your weight properly placed down on the Turning Plate.

WARNINGS

- Make sure all screws and knobs are securely tightened before each use.
- 300lbs weight capacity.
- This device is to only be used by people capable of pulling themselves up and standing on their own powers. Assess the users ability to hold themselves into place before trying the device.

- This device under no circumstances is to be used as a dolly or transportation device to move a patient from point a to point b. This device is ONLY to be used to turn a patient while standing in place to help assist in the transfer of one device to another device within a 1 foot range of the Stand Assist Transfer Disc.
- The wheels on the device are to only move this device itself on its own, the wheels are not to be used with anything or anyone on the device.
- Never leave a patient unattended on the device as there is a tipping hazard if no one is there to keep the patient and the device stable.
- Use proper judgment to assess your weight, strength and capability as the caregiver. For instance it is not recommended for someone who weighs less than 130lbs with less than moderate strength to use this device with someone who weighs over 250lbs.