## **BED ALARM**

## **QUICK START GUIDE**

PLEASE NOTE: Before starting, read through all directions completely. The pager will automatically exit programming mode after 10 seconds of inactivity. If you are not prepared, you may be forced to start again. Also make sure that batteries are installed correctly.

## PROGRAMMING THE TRANSMITTER TO A WIRELESS PAGER

 Press and hold both the UP and DOWN buttons on the pager for more than 5 seconds until you see "1" blinking in red. This means the pager is ready for programming.



There are ten channels for programming.
While the light is flashing, press the button
on the alarm transmitter for 5 seconds. The
light on the transmitter will flash orange.
Release the button.



3. The pager will beep to signal a successful connection. The pager will then re-enter programming mode so that you can connect other devices. If you do not need to connect other devices, press and hold both the UP and DOWN buttons for more than 5 seconds to exit programming mode.

- 4. Connect the sensor pad to the transmitter and test the connection by applying pressure to the pad and then releasing. The alarm should sound on the pager and the connected channel number should flash in red.
- 5. Place the pad in the desired location, making sure the transmitter is mounted in an area that is easily accessible.

## TO SILENCE A TRIGGERED ALARM

When the alarm is triggered, the device number will flash on the pager. Press the STOP button to silence the alarm.



