

# HEATING PAD

LVA1060

vive

QUICK START GUIDE

## CONTROLLER FEATURES

Press button to set heating pad to your desired heat.



### Automatic Shut-off Feature:

This feature is designed to reduce the risk of burns and conserve electricity.

## INSERTING PAD INTO THE PAD COVER

Your heating pad has straps on the ends. To secure the pad in place inside the cover, close the buttons around the loops on the pads.



## MOIST HEAT THERAPY

Dampen the included sponge pad insert for a relaxing moist heat therapy treatment.

QUICK START GUIDE

# vive

Distributed by

vive  
health

## WASHING INSTRUCTIONS

- Hand wash and hang to dry cover only.
- Do not dry clean.
- Do not bleach.
- Do not machine dry.
- Insert the pad into the cover completely, the pad will lie flat and fill the cover completely. The pad should not be folded.
- The cord and switch should remain outside of the cover.

**Note:** Do not over-dampen or submerge your heating pad with water. Be careful and protect clothing and furniture from moisture.

## WARNINGS ⚠

- Do not use while sleeping.
- Do not use on infants or animals.
- Not to be used by children.
- Never use pad without cover in place.
- Do not sit on pad. Only place pad on top of body part.
- Do not tamper with pad in any way.
- Do not use this pad with any ointments.
- Do not submerge pad in water.