PressureGuard® Turn Select®









Customized Lateral Rotation – 24 Hours a Day!

Slowly, gently, and quietly rotates the user laterally up to 30 degrees, with dwell times as short as 1 minute and as long as 2 hours. Hand-held pendant allows convenient programming of the system's four operating modes. Provides caregiver the flexibility to select a standard turning program, or to create a pattern fully customized to the patient's needs.

Ultra-compact, microprocessor controlled pump is integrated directly into the mattress, leaving no bulky pump or external air lines to contend with. Just hour after hour of whisper quiet, trouble-free operation. Maintains proper inflation and support during power loss until power is restored and rotation can be resumed. Includes CPR mode (provides an extra-firm surface for performing CPR with a crash board), program lock-out, and CPR alarm.

Therapy Provided

✓ Programmable Lateral Rotation

OUICK FACTS

Indications:

- Aggressive treatment Stage I IV Pressure ulcers
- Manual turning prohibited due to pain
- Patient can't be positioned off the ulcer¹
- Patient requires custom turning pattern
- Prevention for high risk patients

Key Advantages:

- Pendant controlled, no external pump
- Standard or customized turn patterns
- Won't deflate in power outage
- Ultra-quiet
- PressureGuard® design: stability, safety, anti-shearing

Item #	Description
TS7535-29	PressureGuard Turn Select®, 75"L x 35"W x 7"H
* TS8032-29	PressureGuard Turn Select®, 80″L x 32″W x 7″H
TS8035-29	PressureGuard Turn Select®, 80"L x 35"W x 7"H
TS8435-29	PressureGuard Turn Select®, 84"L x 35"W x 7"H

Wider sizes available: See In-Home Styles – page 19.

 * $\,$ Allow up to 10 working days lead time for this size.

Medicare Approved: E0277 (Group II) Weight limit: 350 lbs.

Warranty: 2 years, non-prorated

Mattress Weight: 28 lbs. UL listed: Yes

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.



¹ In situations where positioning on the pressure ulcer cannot be avoided (e.g. the individual has multiple ulcers on multiple surfaces), limit the amount of time the individual is positioned on the ulcer. Source: NPAUP/EPUAP Pressure Ulcer Prevention and Treatment, Clinical Practice Guideline, p 66-67