



Solutions to Help Prevent
FALLS in Your Facility



Falls continue to be a problem within healthcare facilities today. Each year up to 1 million people fall in U.S. hospitals. Over 30% of those patients sustain an injury that can add on average 6.3 days to the hospital stay costing about \$14,056. An estimated 11,000 of falls each year are fatal. The Centers for Medicare & Medicaid Services (CMS) identifies falls as an event that is preventable and should never occur. Falls and fall-related injuries are considered Hospital-Acquired Conditions (HAC) and reimbursement is limited.*

DeRoyal offers a variety of products and solutions designed to help prevent falls while also helping to protect the patient in the event of a fall. A complete fall prevention program ensures patients will receive top quality care and protection, while also helping the facility reduce expenditures to treat fall-related injuries.

Additional product offerings are available.

Please contact your DeRoyal sales representative or visit www.deroyal.com for more product options.

* Health Research & Educational Trust. (2016, October). Preventing patient falls: A systematic approach from the Joint Commission Center for Transforming Healthcare project. Chicago, IL: Health Research & Educational Trust. Accessed at www.hpoe.org



FALL MONITORS

INDICATIONS: To alert staff and patients in the event the patient attempts to exit the bed, wheelchair, recliner, or toilet.

SAFETY AUTO-RESET™ MONITOR:

- Alarm automatically resets once it has been triggered.
- Patented Caregiver Key™ option allows for only the clinician to reset the monitor so that patients and family members cannot tamper with the device.



M2100-SL

CORDLESS MONITOR:

- Patented cordless technology eliminates all cords and any tripping or entanglement hazards seen with corded models.
- Patented Caregiver Key™ option allows for only the clinician to reset the monitor so that patients and family members cannot tamper with the device.



***M2100-C**

FALL RISK

* All DeRoyal® Fall Monitors come with a protective case, batteries, wall mount option, and wipeable strap mounting.

INTENDED USE

The DeRoyal® Fall Alarm Monitors and Fall Sensor Pads are intended to alert staff and patients in the event the patient attempts to exit the bed, wheelchair and/or recliner.

FALL MONITOR SENSOR PADS

INDICATIONS: To alert staff and patients in the event the patient attempts to exit the bed, wheelchair, recliner, or toilet.

BED AND CHAIR SENSOR PAD:

Weight-sensitive sensor pads allows for early detection of changes in pressure if the patient attempts to exit the bed or chair.



M2200-BP* (Bed)



M2200-CP* (Chair)

CORDLESS BED AND CHAIR SENSOR PAD:

Patented cordless technology eliminates all cords and any tripping or entanglement hazards seen with corded models.



M2200-CBP* (Bed)



M2200-CCP* (Chair)

TOILET ALARM SENSOR PAD

Sensor pad can be placed on top of the feet of the underside of the toilet seat. When the patient or resident attempts to get up from the toilet and pressure is removed from the sensor pad, the monitor will alarm to alert the caregiver.



M2200-TB*

⚠️WARNING- This product can expose you to chemicals including bisphenol A, which is known to the State of California to cause birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

CHAIR BELT SENSORS

INDICATIONS: To alert staff and patients in the event the patient attempts to stand from the seated position.

Fall prevention Chair Belt Sensors help reduce falls by triggering a DeRoyal® fall monitor when the hook and loop or easy release buckle closures are released. These self-releasing options allow for patients to have freedom of movement.



M2200-WHL (Universal)
M2200-WHLX (Extra Long)



***M2200-WCB (Universal)**

***Only the Easy-Release Buckle Belt (M2200-WCB) can be hard mounted to a wheelchair.**

Medical History:

What type of reaction did you have?
 Obtained from patient
 Obtained from family

Do you have a history of?
High Blood Pressure
Allergy to Iodine
Other Allergies
Asthma / COPD / Emphysema
Heart Disease
Myasthenia Gravis
Multiple Myeloma
Kidney Disease
Diabetes

What type of reaction to the Contrast?
 Yes No Unknown
 Yes No Unknown
 Yes No Unknown
 Yes No Unknown
 Yes No Unknown
 Yes No Unknown

If yes, what type?
If yes, what type?
If yes, what type?
If yes, what type?

Ordering p...
ed? Yes No



Fall Monitor Products & Accessories

Product #	Description	Quantity
M2100-C	Fall Monitor - Cordless	EA
M2100-SL	Fall Monitor - Safety Auto-Reset™	EA
M2300-AC05	AC Adaptor - Universal	EA
M2300-CGK	CareGiver Key™ - Universal	20/CS
M2300-NC	Nurse Call Cord - 6 ft	EA
M2300-NCS	Nurse Call Splitter - Universal	EA
M2300-WMP	Wall Mount Kit - Universal	20/CS
M2300-WS	Wipeable Strap - Universal	EA
M2200-BP-CS ^{†*}	Sensor Pad (10"x30") Bulk, Bed	15/CS
M2200-CBP ^{†*}	Sensor Pad (10"x30"), Cordless Bed	EA
M2200-CCP ^{†*}	Sensor Pad (10"x15"), Cordless Chair	EA
M2200-CP-CS ^{†*}	Sensor Pad (10"x15") Bulk, Chair	25/CS
M2200-TB*	Toilet Pad Sensor - Universal	EA
M2200-WCB	Wheelchair Belt Sensor - Universal	EA
M2200-WHL	Chair Belt Sensor - Universal	EA
M2200-WHLX	Chair Belt Sensor - Extra Long	EA

† All Sensor Pads are single patient use with a 30-day warranty.

⚠️*WARNING- This product can expose you to chemicals including bisphenol A, which is known to the State of California to cause birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: Failure to comply with these warnings may result in injury or death. This device is not suitable for all individuals. This device is not a substitute for visual monitoring by caregiver. The manufacturer does not claim that this device will stop elopement and/or stop falls. This device is designed to augment caregivers' comprehensive patient mobility management program. Test this device before each use. Do not use without reading the instructions.

These fall prevention devices are to be installed by a licensed caregiver. As such, it is the entire responsibility of the caregiver to ensure that the system is not designed to replace good caregiving practices including but not limited to:

1. Direct patient supervision
2. Adequate training for staff personnel for fall prevention and/or elopement
3. Testing of the system before each use.



Improving Care. Improving Business.®

DeRoyal, DeRoyal logo, Improving Care. Improving Business., are registered trademarks of DeRoyal Industries, Inc.

Cordless, Caregiver Key, Safety Auto-Reset, and TamperProof are trademarks or registered trademarks of Smart Caregiver.

Copyright© 2018, DeRoyal. All rights reserved.