

SystemRoMedic[®] simple solutions for great results.

OF EASY TRANSFERS

1984-2010

AIRTY YEA

During the 1980s, load injuries increased dramatically within healthcare. At the same time, our company was founded and the very first SystemRoMedic™ transfer aid, the sliding mat EasySlide, was launched. EasySlide represented an entirely new approach and a new way of thinking when facilitating a supine transfer. The positive market response told us we were on the right track.

30 years in the service of simplicity

In 2014, Handicare celebrated its 30th year as a company. For more than a quarter of a century we have been dedicated to developing clever and easy-to-use assistive devices for safe and easy transfers and to making everyday life and work easier for both patients and personnel in healthcare.

And today, with SystemRoMedic[™], we offer the market's broadest and most complete range of easy-to-use, effective and safe assistive devices for all types of patient transfers.

The challenge

Every day, healthcare staff find themselves in situations where they are pressed for time, short of personnel and lacking knowledge and access to simple and functional assistive devices.

Inadequate budgets in combination with the ambition and desire to provide a high standard of care contribute to the occurrence of various load injuries. The result is increased staff absenteeism and escalating costs. Putting an end to this negative spiral can be difficult.

The solution is simple, naturally

The transfer-assistive devices in the SystemRoMedic[™] range are based on simplicity, logical thinking, ergonomic principles and design with the patient and caregiver in mind.

Simplicity is the key. Attention to the basic principles results in user-friendly solutions. And, in fact, it is the simple, most time-saving devices that actually get used. SystemRoMedic™ works with the laws of nature – not against them. Friction, for example, is a factor that must be considered. With smart solutions and combinations of special, functional materials, we can both reduce the friction under the patient to facilitate transfer or repositioning, and increase the friction to minimize the risk of the patient sliding out of bed.



Activating the patient

The more strength and functional capacity the patient has, the more he is able to participate in his own transfer or repositioning. And the more active the patient is, the easier the transfer will be for the caregiver. The SystemRoMedic[™] transfer-assistive devices help the caregiver to activate the patient while stimulating a natural movement pattern. In addition to strengthening the patient's muscles, this also contributes to the patient's sense of independence. In this way, we can improve the quality of life for the patient while at the same time making the caregiver's job easier and safer.

We love easy transfers!

The number of people requiring transfer assistance is steadily increasing. Load injuries and accidents cause huge occupational health and safety problems that account for enormous costs in healthcare.

SystemRoMedic[™] offers simple solutions for great results. It's a matter of simplifying life and work and of preventing load injuries and accidents that occur in connection with patient transfers. It also has to do with activating and rehabilitating people, facilitating transfers, even for painsensitive patients, and preventing pressure sores.

Plus, it is a matter of improving the patient's sense of independence, dignity and self-esteem. With SystemRoMedic™, healthcare providers can take a step in the right direction. Our products are simple, safe and renowned for their cost-effectiveness. Our training programs in transfer technique and product use give caretakers the means to

choose the right working method and assistive device for every transfer situation.

We have seen what our devices can do for people, and this inspires us to continue developing new products. We love easy transfers!

System RoMedic

SystemRoMedic[™] is a complete solution that provides for the majority of patient transfer requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for transferring patients between two locations
- Positioning, assistive devices for repositioning of patients within the same location
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training
- Lifting, assistive devices for manual and mechanical lifting of patients

Through a combination of training and a complete range of effective transfer-assistive devices, SystemRoMedic[™] offers improvement both in the work environment and in the quality of care and, at the same time, achieves significant cost savings.

Find the right transfer aid for the right situation.

SystemRoMedic[™] offers solution for all transfer To make it eas the right transf situation, pleas matrix to com in four transfer 14 most comm transfers.

1

.

EasyBelt/ FlexiBelt/ ReadyBelt

DisposableBelt/ WalkingBelt

31

32

gint ontaation	•															
rstemRoMedic™ offers solutions r all transfer situations and needs.	•	₹7	P	P	Þ	6 1 /2	Ť	ł	~~* \$	~~#	8⇔¥	ĥ	ಸ⇒ಶ	⇔⇔∌	æ	PREVENTS
a make it easier to find and choose e right transfer aid for the right uation, please use the below atrix to combine the main products four transfer categories with the a most common activities in patient ansfers. Products	Activities	Supine transfer	Turning in bed	Higher up in bed	Lying to sitting	Up from the floor	Support for gait training	Standing training	From bed to wheelchair	From bed to commode	From wheelchair to toilet	Rising - sitting down	Positioning in chair	Into and out of vehicles	Positioning in car	Preventing pressure sores
Transfer Sliding aids for manual f	transfers of pat	tients betw	veen two	locations												₽.
EasyGlide/ ReadyGlide	6	~							~	~	V			V		
MoveMaster	7	· ·	V	V					V	V	V			Ť		
EasySlide/ ReadySlide	7	V	V	V									~			~
RollerSlide	8	V														
	0															
Positioning Sliding aids, turning plat	tes and other p	ositioning	aids and	accesso	ries for n	nanual re	position	ing of pa	itients wi	thin the s	ame loca	ation.				
OneWaySlide	9												V			
EasyGlide, oval (mini)	9												V			
AutoSlide	9													~	V	
EasySlide/ EasySlide, disposable	10-11	V	V	V					~				V			V
WendyLett	12, 13		V													V
WendyLett2Way	14		V		V											V
WendyLett4Way	12, 14		V	~	V											V
WendyLean/ ReadyLett/ ReadySh	eet 15	V	V	V	V											V
PositioningCover	16		V	V												
TurnSheet	17-18		V	V												V
TurnTableSoft	19				V										V	
SitWalk	19												V			
ComfortCushion	19			V												
AntiSlip	20			V								V				
FootStool	20								V	V	V		V			
TopSheet	21		V						-	-	-		-			V
IncoSheet	21		V	V												
Support Different types of sit-to-	stand aids, bel	ts, turning	plates a	nd other s	support a	ids for n	nobility s	upport e	.g., durinę	g sit-to-s	tand or g	ait train	ing.			
MiniLift200/ MiniLift160/ MiniLift12	25 22							V	~	V	V	V	V			
ThoraxSling	23, 72							V	V	V	V	V	V			
Rena440/ Rena600	24							V	~	V	~	V	V			
StandUpSling	24, 73							V	~	V	V	V	V			
QuickMove	25							V	V	~	V	V	V			
FlexiMove/ ReadyMove	25, 26			V	V				V		V	V				
FlexiGrip/ ReadyGrip	26				V				V							
StandUp	27					V										
ReTurn	28-29							V	V	V	~	~	~			
ReTurnBelt/ SupportStraps	30						V	V				V				
		1							1				1			

V

1

1

1

V

1

System**RoMedic**™

		₹7	P	Þ	P		Ť	ħ	~7# ≵	~~#	^{⊜⇔} ≯	ħ	えきた	⇔ಾಶ	æ	PREVENTS
	Activities	Supine transfer	Turning in bed	Higher up in bed	Lying to sitting	Up from the floor	Support for gait training	Standing training	From bed to wheelchair	From bed to commode	From wheelchair to toilet	Rising - sitting down	Positioning in chair	Into and out of vehicles	Positioning in car	Preventing pressure sores
Products	Page						Su		Ë		Fro			_n		Prev
RoWalker400	32						V	V								
TurnTable	33								V	V	V					
TurnSafe	33							V	V	V	V	V	V			
Lifting Manual lifting aids	and mechanical lifti	na aids lik	e ceilina	lifts. rail	svstems	. mobile i	ifts. liftiı	na slinas	and acce	essories f	or manu	al or me	chanical	lifting of	patients	
LiftSeat	36					~			V	\checkmark	V			V		
FlexiMove	36	V	V	V		V				-				-		
MiniStretcher	37	V	V	V		V										
RiseBasic/ RiseAtlas, stationary ceiling	lift units 38-41	V	V	V	V	V	V	V	V	V	V	V	V			
RisePorto, portable ceiling lift units	35, 42	V	V	V	V	V	V	V	V	V	~	V	V			
Castor/ Pollux/ Helena, freestanding I		V	V	V	V	V	V	V	V	V	V	V	V			
MilkyWay, rail system	43, 46-47	V	V	V	V	V	V	V	V	V	V	V	V			
Eva, mobile lifts	48-49	V	V	V	V	V	V	V	V	V	V	V	V			
Rino450/ Rino600/ Rino1000/ C	arina350 50	V	V	V	V	V	V	V	V	V	V	V	V			
Vega505EE	51	V	V	V	V	V	V	V	V	V	V	V	V			
HighBackSling	60-61					V			V				V			
ShellHBSling	61					V			V				V			
FlexibleSling	62, 69					~			V	V	V		V			
AmpHBSling	62					-			V				V			
HygieneHBSling	63									V	V		V			
BasicSling	64								V				V			
AmpSling	65								V				V			
LowBackSling	66								V				V			
HygieneSling	67									V	V					
HygieneLBSling	68									V	V		V			
ClassicHBSling	69								V				V			
ClassicSling	70								V	V	V		V			
BariSling	70								V							
StretcherSling	71	V		V		V										
SafeHandlingSheet	71	V	V	~	V	~										
TurningSling	71		V													
StandingVest	74							~								
WalkingVest	74				V		V	V				V				
BariVest/ LimbSling	75						V	V				V				

More content

 Selecting sling model, material and size.....56-59 Size guide for lifting slings, overview.......76-79

Sliding aids



EasyGlide 5020



EasyGlide 5060

EasyGlide

For sitting transfers from wheelchairs to beds, toilets, car seats and many other seated positions. It is especially suitable for leveling out height differences or for bridging short distances. Available in two models and several sizes.

Art. No.	Product name	Description	Unit
5020	EasyGlide	60 x 33 cm/ 23.6" x 13", three-part	1 Pce
5050	EasyGlide	75 x 33 cm/ 29.5" x 13", three-part	1 Pce
5060	EasyGlide	60 x 25 cm/ 23.6" x 9.8"	1 Pce
5070	EasyGlide	45 x 20 cm/ 17.7" x 7.9"	1 Pce

Accessories

4014

Storage pocket For 5020, 5050





ReadyGlide

ReadyGlide is a sliding board that has been designed to simplify transfers in a seated position, for example, between bed and wheelchair or between wheelchair and toilet.

Art. No.	Product name	Description	Unit
15030	ReadyGlide	60 x 25 cm/ 23.6" x 9.8", polyethylene	1 Pce
		CESTA	€⇔≵



EasyGlide

For supine transfers from one bed to another.

Art. No.	Product name	Description	Unit
5000	EasyGlide, extra long	120 x 60 cm/ 47.2" x 23.6"	1 Pce
5010	EasyGlide, extra long	170 x 59 cm/ 66.9" x 23.2"	1 Pce
Glide 5010			(E #7



EasyGlide 5042

EasyGlide

For supine transfers from one bed to another but also for positioning higher up in bed and for placing a lifting sling.

Art. No.	Product name	Description	Unit
5040	EasyGlide, oval	85 x 39 cm/ 33.5" x 15.4"	1 Pce
5042	EasyGlide, oval	85 x 39 cm/ 33.5" x 15.4"	2 Pcs

EasyGlide, oval is also available in a smaller model for positioning further back in chair, see page 9.





MoveMaster

EasySlide

MoveMaster is a slightly different, open sliding mat with handles at one end. This makes it possible to help by pulling patients in various transfers, such as moving them further up in bed. MoveMaster is particularly well suited for use against bare skin when transferring to shower chairs or to the toilet.

Art. No.	Product name	Description	Unit
2090	MoveMaster	47 x 120 cm/ 18.5" x 47.2", 3 handles	1 Pce
		œCE ፇ₽₽ ๙๙ฃ₽	ઝ <i>~</i> ેટ્ટે

EasySlide (190 x 60) is primarily used for supine transfers, for example between beds and x-ray or operating tables, shower trolleys, stretchers and similar. EasySlide is also a very effective aid for turning patients in bed or on x-ray or operating tables. Disposable covers for EasySlide provide for improved hygiene and less laundry.



Protective Cover

WallRack

Art. No.	Product name	Description	Unit
1020	EasySlide	190 x 60 cm/ 74.8" x 23.6"	1 Pce
1030	EasySlide (1020 + 1051)	190 x 60 cm/ 74.8" x 23.6" complete	1 Set
Accessories	5		
1051	Protective cover	For 1020, 1030	1 Pce
1055	Disposable cover	For 1020, 1030	2 x 50 Pcs
1151	Protective cover	For 1120*, 1130*	1 Pce
1155	Disposable cover	For 1120*, 1130*	2 x 90 Pcs
1085	WallRack	For storing of sliding mats and	1 Pce
		disposable covers for EasySlide and Rolle	rSlide
* 1120, 1130, Ea	ssySlide 100 x 60, discontinued		



ReadySlide

ReadySilde is a thin, pliable and very easy-to-use sliding mat that facilitates various types of lying transfers of patients, positioning in chairs and application of lifting slings.

Art. No.	Product name	Description	Unit
11250	ReadySlide	70 x 70 cm/ 27.6" x 27.6", polyamid	1 Pce
11230	ReadySlide	110 x 70 cm/ 43.3" x 27.6", polyamid	1 Pce
11235	ReadySlide	110 x 90 cm/ 43.3" x 35.4", polyamid, extra wide	1 Pce
11210	ReadySlide	195 x 70 cm/ 76.8" x 27.6", polyamid	1 Pce
11215	ReadySlide	195 x 70 cm/ 76.8" x 27.6", polyamid, with handles	1 Pce



Sliding aids





RollerSlide

For supine transfers of inert patients between a bed and an x-ray or operating table, shower trolley, ambulance stretcher or similar. Bridges gaps and moves with the patients during the transfer.

RollerSlide consists of a padded sliding board with a cover made from ULF™ (Ultra Low Friction), a sophisticated sliding material that slides round the board during transfer. RollerSlide is available in three models and several sizes. RollerSlide, flexi is a lighter and bendable variant of RollerSlide.

RollerSlide, light is a light-weight variant of RollerSlide, flexi.

RollerSlide, flexi is also available in two wider variants adapted for larger and heavier patients weighing up to 300 kg / 660 lbs, depending on the distance between the underlying surfaces.

Art. No.	Product name	Description	Unit
9145	RollerSlide	175 x 50 cm/ 68.9" x 19.7", foldable	1 Pce
9155	RollerSlide	85 x 50 cm/ 33.5" x 19.7"	1 Pce
9148	RollerSlide, flexi	175 x 50 cm/ 68.9" x 19.7, foldable	1 Pce
9158	RollerSlide, flexi	85 x 50 cm/ 33.5" x 19.7"	1 Pce
9120	RollerSlide, light	175 x 50 cm/ 68.9" x 19.7", foldable	1 Pce
9110	RollerSlide, light	85 x 50 cm/ 33.5" x 19.7"	1 Pce
9170	RollerSlide, flexi	175 x 60 cm/ 68.9" x 23.6", foldable	1 Pce
9175	RollerSlide, flexi	85 x 60 cm/ 33.5" x 23.6"	1 Pce





Accessories 9136 9156

9172 9177 9142

6	Sliding cover	For 9120, 9145, 9148	1 Pce
6	Sliding cover	For 9110, 9155, 9158	1 Pce
2	Sliding cover	For 9170	1 Pce
7	Sliding cover	For 9175	1 Pce
2	Sliding cover, wipeable	For 9120, 9145, 9148	5 Pcs
)	Sliding cover, wipeable	For 9110, 9155, 9158	5 Pcs
3	Sliding cover, wipeable	For 9170	5 Pcs
3	Sliding cover, wipeable	For 9175	5 Pcs
3	Disposable cover	For 9120, 9145, 9148	2 x 50 Pcs
3	Disposable cover	For 9110, 9155, 9158	2 x 115 Pcs
5	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce

Patent: AT-1648368, CH-1648368, DE-502004013222.7





OneWaySlide 2069



OneWaySide has many functions. The main one improves sitting positions in chairs or wheelchairs. The sliding mat effectively retards forward movement yet slides backwards, making it easy for the caregiver to help the patient move into a more comfortable position in the chair. Also available in a tube-shaped design without handles for the more active patient.

Art. No.	Product name	Description	Unit
2074	OneWaySlide	110 x 30 cm /43.3" x 11.8", 6 handles	1 Pce
2075	OneWaySlide	130 x 40 cm/ 51.2" x 15.7", 6 handles	1 Pce
2076	OneWaySlide	130 x 45 cm/ 51.2" x 17.7", 6 handles	1 Pce
2077	OneWaySlide	130 x 55 cm/ 51.2" x 21.7", 6 handles	1 Pce
2069	OneWaySlide, tubular	50 x 40 cm/ 19.7" x 15.7"	1 Pce
2070	OneWaySlide, tubular	50 x 45 cm/ 19.7" x 17.7"	1 Pce
			€ ಸ⇒ಸ



EasyGlide

EasyGlide are used in pairs to facilitate positioning further back in chair. The smaller size makes it easier to place the sliding boards under the patient's thighs and seat.

Art. No.	Product name	Description	Unit
5045	EasyGlide, oval	45 x 20 cm/ 17.7" x 7.9"	2 Pcs
			(€ ಸ್ಸ್



AutoSlide

Facilitates getting in and out of cars, allowing patients to turn easier and move sideways. The seat cushions are firmly attached to the backs of car seats, so patients cannot slide forward.

Art. No.	Product name	Description	Unit
8100	AutoSlide	46 x 46 cm/ 18.1" x 18.1"	1 Pce



Sliding aids



EasySlide

EasySlide (100 x 70) is mainly designed for patients requiring only some care who need to turn and move in bed. It performs best on a rough surface. EasySlide reduces the friction on both shoulder and hip, which in turn reduces the risk of pressure sores. Thanks to these features, EasySlide is very gentle for patients in pain but also for the caregiver, thus avoiding the risks associated with lifting.

Art. No.	Product name	Description	Unit
1220	EasySlide	100 x 70 cm/ 39.4" x 27.6"	1 Pce
1230	EasySlide (1220 + 1251)	100 x 70 cm/ 39.4" x 27.6", complete	1 Set
Accessor	ies		
1251	Protective cover	For 1220, 1230	1 Pce
1255	Disposable cover	For 1220, 1230	2 x 90 Pcs
1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce



EasySlide 1430

WallRack

1599

1085

EasySlide

EasySilde (140 x 90) is designed for normal-weight patients with limited to severely reduced abilities to participate in the turning process. These patients are often permanently confined to the bed, and it is therefore important to prevent pressure sores. EasySlide (200 x 90) has extra padding in order to further reduce the risks of pressure sores. The padding distributes the pressure over the entire body and reduces pressure, friction and the risk of shearing on the especially vulnerable pressure points: the shoulders, hips, head, and feet.

Art. No. **Product name** Description Unit 140 x 90 cm/ 55.1" x 35.4", 2 fastening straps 1420 EasySlide 1 Pce EasySlide (1420 + 1451) 140 x 90 cm/ 55.1" x 35.4", complete 1430 1 Set 1470 EasySlide 200 x 90 cm/ 78.7" x 35.4", 4 fastening straps 1 Pce EasySlide (1470 + 1481) 200 x 90 cm/ 78.7" x 35.4", complete 1480 1 Set Accessories 1451 Prot 1481 Prot



WallRack

5		
Protective cover	For 1420, 1430	1 Pce
Protective cover	For 1470, 1480	1 Pce
LockingStrap	Fastening straps, spare part for 1420, 1430, 1470, 1480	1 Pair
WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce





EasySlide

EasySlide (50 x 60), (50 x 70) and (25 x 60) are excellent aids for turning in bed, for moving patients higher up in bed or for transferring them from the bed to the wheelchair. EasySlide is a small, practical aid that encourages activation. For the caregiver, EasySlide decreases the number of assisted positionings and eliminates unnecessary lifting, which spares both back and shoulders.

Art. No.	Product name	Description	Unit
2010	EasySlide	50 x 60 cm/ 19.7" x 23.6", polyester	1 Pce
2020	EasySlide	50 x 60 cm/ 19.7" x 23.6", microfiber	1 Pce
2030	EasySlide (2020 + 2051)	50 x 60 cm/ 19.7" x 23.6", complete	1 Set
2011	EasySlide	50 x 70 cm/ 19.7" x 27.6", polyester	1 Pce
2021	EasySlide	50 x 70 cm/ 19.7" x 27.6", microfiber	1 Pce
2015	EasySlide	25 x 60 cm/ 9.8" x 23.6", polyester	1 Pce

Accessories 2051

Protective cover	For 2010, 2020, 2030	1 Pce
Protective cover	For 2011, 2021	1 Pce
Disposable cover	For 2010, 2020, 2030	100 Pcs
WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce





WallRack

Disposable cover

EasySlide, disposable

EasySlide, disposable is a thin, pliable and very easy-to-use personal sliding mat for short-term use. EasySlide, disposable can be used to facilitate turning and repositioning higher up in bed, for supine and seated transfers, and for repositioning in the seated position. EasySlide, disposable can also be used to lift the patient's legs in and out of bed, to assist dressing and undressing, and to facilitate application of a lifting sling. EasySlide, disposable is a new, unbelievably versatile sliding mat that makes many common transfers, even heavy ones, easier for both the patient and the caregiver. The sliding mat is made of environmentally friendly polyethylene (PE) plastic and the practical, hygienic package contains a roll with up to 60 individual EasySlide mats for short-term use. EasySlide, disposable is also available in a practical and hygienic single-pack which is very convenient to always have close at hand.

M	EasySlide and Ball
2	Prosection and
14	
1	
1	AAA

EasySlide, disposable, 2040



including practical wall and tabletop dispenser in steelEasySlide, disposable80 cm x 70 cm/ 31.5" x 27.6", single-packed140 PcsEasySlide, disposable120 cm x 70 cm/ 47.2" x 27.6", single-packed100 Pcs	Art. No.	Product name	Description	Un
including practical wall and tabletop dispenser in steelEasySlide, disposable80 cm x 70 cm/ 31.5" x 27.6", single-packed140 PcsEasySlide, disposable120 cm x 70 cm/ 47.2" x 27.6", single-packed100 PcsEasySlide, disposable200 cm x 70 cm/ 78.7" x 27.6", single-packed70 Pcs	2040	EasySlide, disposable	50 cm x 70 cm/ 19.7" x 27.6"	2 x 60 Pc
EasySlide, disposable120 cm x 70 cm/ 47.2" x 27.6", single-packed100 PcsEasySlide, disposable200 cm x 70 cm/ 78.7" x 27.6", single-packed70 Pcs	2045	EasySlide, disposable		2 x 60 Pc
EasySlide, disposable 200 cm x 70 cm/ 78.7" x 27.6", single-packed 70 Pcs	2041	EasySlide, disposable	80 cm x 70 cm/ 31.5" x 27.6", single-packed	140 Pc
	2043	EasySlide, disposable	120 cm x 70 cm/ 47.2" x 27.6", single-packed	100 Pc
CE # # # # 10 10 10 10 10 10 10 10 10 10 10 10 10	2042	EasySlide, disposable	200 cm x 70 cm/ 78.7" x 27.6", single-packed	70 Pc
	2042	, ,	200 cm x 70 cm/ 78.7" x 27.6", single-packed	

Sit-to-stand aids



ReTurnBelt







HeelStrap



Padding for leg support, 7215



Bag

Art. No.	Product name	Description	Unit
7500	ReTurn7500	Transfer platform, complete with standard-height rising ladder	1 Pce
7500i	ReTurn7500	Transfer platform, complete with standard-height rising ladder, with an opening for hooking of ReTurnBelt	1 Pce
7400 7600	ReTurn7400 ReTurn7600	Transfer platform, complete with low rising ladder Transfer platform, complete with wide rising ladder	1 Pce 1 Pce
7501	Rising ladder, standard height	For ReTurn7500 with serial numbers from 373317697500500001 and for ReTurn7400 with serial numbers from 373317697400500001	1 Pce
7501i	Rising ladder, standard height	With an opening for hooking of ReTurnBelt, for ReTurn7500 with serial numbers from 373317697500500001 and for ReTurn7400 with serial numbers from 373317697400500001	1 Pce
7401	Rising ladder, low	For ReTurn7400 with serial numbers from 373317697400500001 and for ReTurn7500 with serial numbers from 373317697500500001	1 Pce

Δ.	cce	ee/	orie	20
- A9	Jue	ວວເ	טו וכ	:5

Accessorie	S		
60316038	ReTurnBelt***	Size XXS-XXXL, for ReTurn with article numbers 7100****, 7500, 7500i, 7400, 7600	1 Pce
61336136	ReTurnBelt	Wipeable, size S-XL	1 Pce
62336236	ReTurnBelt	Disposable, non-woven, size S-XL	5 Pcs
6061	SupportStraps	Accessory for belts	1 Pair
7250	HeelStrap	For ReTurn7500 with serial numbers from 373317697500500001 and for ReTurn7400 with serial numbers from 373317697400500001	1 Pce
7215	Padding for leg support	For ReTurn with article numbers 7100****, 7500, 7500i, 7400 (standard equipment for 7500i) and for TurnSafe	1 Pair
7216	Padding for leg support	For ReTurn with article numbers 7600	1 Pair
0954	Bag	For ReTurn with article numbers 7500, 7500i, 7400	1 Pce
		4	

) Patent: SE-0401517-8 *) ReTurn7100 is discontinued (દ 🖌 🖓 અન્ય 🖓 🚽 🖌



Belts



ReTurnBelt



ReTurnBelt

Our patented* ReTurnBelt is an excellent complement to ReTurn when there is a need for greater support and stability; for example, when the patient's weight-bearing ability is questionable or varies. ReTurnBelt can also be used separately, for support during sit-to-stand and to stabilize the patient in an upright position, for example, during standing training, gait training and transfer. ReTurnBelt has an anti-slip inside of PU-coated nylon and can therefore be used in environments where moisture and dampness are present or when there are special demands on hygiene. The fabric also ensures that the belt stays in place. ReTurnBelt is also available in a wipeable material, which can be cleaned using disinfectants, and in disposable material which is disposed after use instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacteria.

Art. No.	Product name	Size	Measurements	Handles	Unit
6031	ReTurnBelt	XXS	65 - 80 cm/ 25.6 - 31.5"	4 pcs	1 Pce
6032	ReTurnBelt	XS	80 - 90 cm/ 31.5 - 35.4"	4 pcs	1 Pce
6033	ReTurnBelt	S	90 - 105 cm/ 35.4 - 41.3"	4 pcs	1 Pce
6034	ReTurnBelt	Μ	105 - 125 cm/ 41.3 - 49.2"	4 pcs	1 Pce
6035	ReTurnBelt	L	125 - 145 cm/ 49.2 - 57.1"	4 pcs	1 Pce
6036	ReTurnBelt	XL	135 - 155 cm/ 53.2 - 61.0"	4 pcs	1 Pce
6037	ReTurnBelt	XXL	155 - 175 cm/ 61.0 - 68.9"	4 pcs	1 Pce
6038	ReTurnBelt	XXXL	175 - 195 cm/ 68.9 - 76.8"	4 pcs	1 Pce
6133	ReTurnBelt, wipeable	S	90 - 105 cm/ 35.4 - 41.3"	2 pcs	1 Pce
6134	ReTurnBelt, wipeable	Μ	105 - 125 cm/ 41.3 - 49.2"	2 pcs	1 Pce
6135	ReTurnBelt, wipeable	L	125 - 145 cm/ 49.2 - 57.1"	2 pcs	1 Pce
6136	ReTurnBelt, wipeable	XL	135 - 155 cm/ 53.2 - 61.0"	2 pcs	1 Pce
6233	ReTurnBelt, disposable, non-woven	S	90 - 105 cm/ 35.4 - 41.3"	2 pcs	5 Pcs
6234	ReTurnBelt, disposable, non-woven	Μ	105 - 125 cm/ 41.3 - 49.2"	2 pcs	5 Pcs
6235	ReTurnBelt, disposable, non-woven	L	125 - 145 cm/ 49.2 - 57.1"	2 pcs	5 Pcs
6236	ReTurnBelt, disposable, non-woven	XL	135 - 155 cm/ 53.2 - 61.0"	2 pcs	5 Pcs

Accessories

6061 SupportStraps

Accessory for belts

1 Pair

(€∦†

*) Patent: SE-0401517-8



SupportStraps

SupportStraps is a multi-functional accessory providing extra support and security for patients during standing, gait training and transfers with Handicare's SystemRoMedic[™] belts, sit-to-stand aids and transfer platforms. SupportStraps can be used as:

- leg harness for EasyBelt/FlexiBelt,
- torso/side support when using QuickMove together with EasyBelt/FlexiBelt/FlexiMove
- side support when using ReTurn together with ReTurnBelt

Art. No.	Product name	Description	Unit
6061	SupportStraps	Accessory for belts	1 Pair



FlexiBelt / EasyBelt

FlexiBelt and EasyBelt are used for sit-to-stand transfers or for moving from the wheelchair to the bed, toilet or car. FlexiBelt has a soft innerside of polyester velvet which makes it very comfortable to wear. EasyBelt has an anti-slip inside of PU-coated nylon and can therefore be used in environments where moisture and dampness are present or when there are special hygiene requirements. EasyBelt is also available in a wipeable material which can be cleaned using disinfectants.

FlexiBelt, hug and EasyBelt, hug are new versions of Handicare's classic SystemRoMedic™ belts. FlexiBelt, hug and EasyBelt, hug feature a new patented* design with horizontal handles that slide freely in loops to make the belt "hug" on to the patient instead of sliding up during transfer situations. With FlexiBelt, hug and EasyBelt, hug the caregiver can provide even more support without using more effort. The new design actually does all the work by itself.

	Art. No.	Product name	Description	Unit
	6023	FlexiBelt	S, 60 - 100 cm/ 23.6 - 39.4", 7 handles	1 Pce
	6024	FlexiBelt	M, 70 - 120 cm/ 27.6 - 47.2", 9 handles	1 Pce
	6025	FlexiBelt	L, 100 - 160 cm/ 39.4 - 63", 11 handles	1 Pce
(Restant	6026	FlexiBelt	XL, 115 - 175 cm/ 45.3 - 68.9", 13 handles	1 Pce
FlexiBelt, hug				
	6072	FlexiBelt, hug	XS, 45 - 70 cm/ 17.7 - 27.6", 3 + 2 movable handles	1 Pce
	6073	FlexiBelt, hug	S, 60 - 100 cm/ 23.6 - 39.4", 4 + 3 movable handles	1 Pce
and the second se	6074	FlexiBelt, hug	M, 70 - 120 cm/ 27.6 - 47.2", 5 + 4 movable handles	1 Pce
	6075	FlexiBelt, hug	L, 100 - 160 cm/ 39.4 - 63", 6 + 5 movable handles	1 Pce
	6076	FlexiBelt, hug	XL, 115 - 175 cm/ 45.3 - 68.9", 7 + 6 movable handles	1 Pce
	6012	EasyBelt	XS, 45 - 70 cm/ 17.7 - 27.6", 5 handles	1 Pce
	6013	EasyBelt	S, 60 - 100 cm/ 23.6 - 39.4", 7 handles	1 Pce
	6014	EasyBelt	M, 70 - 120 cm/ 27.6 - 47.2", 9 handles	1 Pce
	6015	EasyBelt	L, 100 - 160 cm/ 39.4 - 63", 11 handles	1 Pce
EasyBelt	6002	EasyBelt, hug	XS, 45 - 70 cm/ 17.7 - 27.6", 3 + 2 movable handles	1 Pce
EdSyden	6003	EasyBelt, hug	S, 60 - 100 cm/ 23.6 - 39.4" 4 + 3 movable handles	1 Pce
	6004	EasyBelt, hug	M, 70 - 120 cm/ 27.6 - 47.2", 5 + 4 movable handles	1 Pce
1	6005	EasyBelt, hug	L, 100 - 160 cm/ 39.4 - 63", 6 + 5 movable handles	1 Pce
C- Ch	6006	EasyBelt, hug	XL, 115 - 175 cm/ 45.3 - 68.9", 7 + 6 movable handles	1 Pce
	6113	EasyBelt, wipeable	S, 60 - 100 cm/ 23.6 - 39.4", 7 handles	1 Pce
	6114	EasyBelt, wipeable	M, 70 - 120 cm/ 27.6 - 47.2", 9 handles	1 Pce
EasyBelt, hug	6115	EasyBelt, wipeable	L, 100 - 160 cm/ 39.4 - 63", 11 handles	1 Pce
	6116	EasyBelt, wipeable	XL, 115 - 175 cm/ 45.3 - 68.9", 13 handles	1 Pce



cessories

6061

SupportStraps

*) Patent: SE-536273, EP-2611409, Patent pending: JP-2013-527040

Accessory for belts



ReadyBelt

ReadyBelt is a belt that is designed to provide the patient with support and assistance during mobility training and different transfer situations. For the caregiver, ReadyBelt offers an ideal means of supporting the patient during mobility training or visits to the toilet.



).	Product name	Description	Unit
	ReadyBelt	Small, 105 x 10 cm/ 41.3" x 3.9", 7 handles	1 Pce
	ReadyBelt	Medium, 120 x 10 cm/ 47.2" x 3.9", 9 handles	1 Pce
	ReadyBelt	Large, 140 x 10 cm/ 55.1" x 3.9", 11 handles	1 Pce

