Enabling Devices Capability Switches Comparison Chart

Best Selling Switches

*Activation Force - the amount of pressure needed on a switch to turn on a device. For your reference, a nickel weighs 5g.

DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Gumball Switches #782R, Y, B, G, P, #7801		Auditory click, Tactile	25 g	Pressure	2½"D	3¼" x ¾" ¼ lb	Red, Yellow, Blue, Green, Purple
Jumbo Switches #787R, Y, B, G, P. #789		Auditory click, Tactile	20 g	Pressure	5 1/8"	5 1/8" x 1" ¾ lb	Red, Yellow, Blue, Green, Purple
Mini-Gumball Switches #783R, Y, B, G, P, #7802		Auditory click, Tactile	35 g	Pressure	1½"D	2" x ½" 1 oz	Red, Yellow, Blue, Green, Purple

	Hand, Finger, Foot & Body Switches										
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR				
Adjustable Angle Switch #1032		Auditory click, Tactile	25 g	Pressure	5"	6½" x 8" x 8" 1 lb	Yellow/Black				
Button Click Switch #775		Auditory click, Tactile	40 g	Pressure	1 7/8"	2½" x 2½" ¼ lb	Yellow/Blue				
Chin Switch #6001		Auditory click	40 g	Pressure	2" x 1"	2½"W x 2½" L ½ lb	Black				
Compact Switches (Set of 3) #745	PROPES 1	Auditory click, Tactile	100 g	Pressure	³⁄4" x 3/8"	³¼" x 3/8" ¹¼ lb	Red, Blue, Yellow				
Cushion Grip Switch (#3141 Small & #3137 Medium)		Tactile, Auditory Click	40 g	Pressure	5" x 2" x 2"	5" x 2" x 2" ¼ lb	Blue				
Eye Blink Switch #9008		Visual	N/A	Eye Blink	N/A	4½" x 2½" x 1" ½ lb	Black/Silver				

DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Finger Extension Pressure Switch #724X	The state of the s	Tactile	Adjustable from 1 g to 40 g	Pressure, Adjustable		9" ½ lb	Black
Finger Isolation Button #716		Tactile	100 g	Pressure	1½"	3" x 1½" 3 oz	Yellow
Finger Switch #717		Tactile	15 g	Pressure		1½" x ½" 1 oz	Blue
Foot Switch #6000	NEW CONTRACTOR OF THE PARTY OF	Auditory click	5 lbs	Pressure	3½" x 2½"	3½" x 2½" ½ lb	Black
Joystick #720		Tactile	50 g	Pressure	3" x ¾"	7" x 6½" x 4½" ¾ lb	Red/Black
Pinch Switch #800		Tactile	100 g	Pressure, Pinch	3" x ¼"	3" x ¼" 2 oz	Red

DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
The Twitch Switch #1061		N/A	Adjustable from 1 g to 40 g	Pressure Adjustable	1½" x ½"	½ lb	White
Thumb Switch #802		Tactile	100 g	Pressure	1½"	4½" x 1" 3 oz	Yellow
Tilt Switch #990		None	N/A	Tilt Movement		3" x 2" x 1" 2 oz	Black
Wobble Switch #1070		Auditory click, Tactile	10 g	Pressure	4" Spring	24" x 2" ½ lb	Red/Yellow

	Lights, Vibration & Music Switches									
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR			
Jumbo Switch with Vibration #791		Auditory click, Visual, Tactile. Vibration		Pressure	5"	5" x ½" 1 lb	Yellow			
Jumbo Switch w/Lights, Music and Vibration #768		Auditory click, Visual, Tactile, Vibration Lights	20 g	Pressure, Visual	5"	5" x ½" ¾ lb	Green			
Silly Pin Switch w/Lights, Vibration & Music #209		Auditory click, Visual, Tactile Vibration	60 g	Pressure	N/A	8¼" x 9" x 7½" 1½ lbs	Blue/Green			
*Switch & Glow #894 *Not suitable for those who are allergic to Latex.		Tactile, Visual Light, Music Vibration, Lights	150 g	Pressure, Visual	N/A	7" x 6½" x 7" 1 lb	Green/Blue			
Textured Orbit Ball Switch #887		Auditory click, Visual, Tactile, Lights	90 g	Pressure	5½" x 4½"	5½" x 4½" 1 lb	Yellow/Blue			

	Switches for Visually Impaired										
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR				
Bright Switch for the Visually Impaired #2045		Visual, Tactile Light	50 g	Pressure	N/A	7" x 3¾" 1 lb	Black				
Bright Red Switch #262		Auditory click, Visual, Tactile, Lights & Vibration		Pressure	N/A	5" x 2½" ¾ lb	Blue				
Large Lens Switch #957		Visual, Tactile, Lights, Music, Vibration, Auditory Click	90 g	Pressure	N/A	8" x 2" 1¼ lbs	Blue				
Tail Light Switch, Say It Play It #464		Lights. Visual, Tactile, Voice Output, Vibration, Auditory Click	90 g	Pressure	N/A	7¾" x 2½" 1 lb	Silver, Yellow				
Tail Light Switch w/Vibration & Music #465		Music Visual, Tactile, Lights & Vibration, Auditory Click	90 g	Pressure	N/A	7¾" x 2½" 1 lb	Silver, Green				

	Pillow Switches									
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR			
Mother's Third Arm w/Pillow Switch #9090		Auditory click, Tactile	40 g	Pressure	3"	26" x 3¾" ¾ lb	Black & Green			
Petite Pillow Switches #831, S, R, F		Auditory click, Tactile	40 g	Pressure	3"	3" x 2" 3 oz	Green, Red, Safari, Floral			
		Р	late Swit	ches						
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR			
Plate Switch #810		Tactile, Auditory Click	50 g	Pressure	3" x 5"	3" x 5" ½ lb	Yellow/Blue			
Sensory Plate Switch #886		Auditory click, Visual, Tactile, Lights, Vibration	75 g	Pressure, Visual	5" x 8"	11¾" x 8½" x 3" 1¾ lbs	Yellow/Black			
Vibrating Plate Switch #1040		Tactile, Auditory Click, Vibration	50 g	Pressure	5" x 8"	5" x 8" ¾ lb	Yellow/Blue			

	Saucer Switches										
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR				
Jumbo Switch with Latch Timer #792		Tactile	20 g	Pressure	5"	5" x ½" 1 lb	Red				
Mini-Saucer Switch #718-S		Auditory click, Tactile	40 g	Pressure	3¼" x 3¼" x 1"	2¼" ¼ lb	Blue/Yellow				
Modified Pancake Switches #769R, B, Y, G, #766		Auditory click, Tactile	40 g	Pressure	5 3/8"	7" ½ lb	Red, Blue, Yellow, Green				
On/Off Jumbo Switch #7805		Tactile	20 g	Pressure	5"	5" x ½" ½ lb	Yellow				
Saucer Switch #718		Auditory click, Tactile	25 g	Pressure	7" x 1"	6" ½ lb	Blue/Yellow				

	Sip, Puff & Breath Switches									
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR			
*Grip and Puff Switch #700 *Not returnable due to sanitary reasons noted by OSHA standards. Straws should be sterilized before use.		None	20 g	Pneumatic, Pressure or Puff	N/A	½ lb	White			
*Sip and Puff Switch #970 *Not returnable due to sanitary reasons noted by OSHA standards. Straws should be sterilized before use.		None	20 g	Pneumatic, Pressure, Sip or Puff	Straw	19" Gooseneck 2¼ lbs	Black			

	Textured Switches										
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR				
Gooshy Switch #634		Tactile, Visual, Auditory Click	90 g	Pressure	N/A	10" x 4½" x 2"	Blue/Yellow				
Gooshy Switch w/Vibration, Lights & Music #636		Tactile, Visual, Lights, Music, Vibration, Auditory Click	90 g	Pressure	N/A	10" x 4½" x 2" ¼ lb	Blue/Yellow				
*Koosh Switch #883 *Not suitable for those who are allergic to Latex.		Tactile, Visual Light, Music Vibration, Auditory Click	30 g	Pressure	N/A	8½" x 5¼" x 2½" 1 lb	Yellow				
Multi-Sensory Center #725, #726		Auditory click, Tactile Light, Music Vibration	40 g	Pressure		9" x 6" 1 lb 11" x 7" 1¼ lbs	Blue/Yellow				
Senso Dot Switch #682 & Set of 4 #683		Tactile, Vibration, Auditory Click	40 g	Pressure	6" x 2¾"	6" x 2¾" 1 lb	Yellow/Green, Color may vary				
Small Oval Textured Switches (Blue) #711-S		Auditory click, Tactile	40 g	Auditory click, Tactile	4½"	5" x 4½" x 1" ¾ lb	Black/Yellow				

			ACTIVATION	ACTIVATION	ACTIVATION	SIZE/	
DEVICE	РНОТО	FEEDBACK	FORCE*	TYPE	SURFACE	WEIGHT	COLOR
Soft Shaggy Switch #3133		Tactile	40 g	Presure	7"	7" x 5" x 2" ¼ lb	Green/Color may vary
Squeeze Switch Kit #6430		Tactile	20 g 20 g 10 g	Pressure	Varies	6" x 2½" x 2½" 2¼" x 2¼" x 2¼" 3½" Round ¾ lb	Yellow, Green
Super Bumpy Switch #736		Auditory click, Tactile	30 g	Auditory click, Tactile	8" x 4½"	8" x 4½" x 1/8" ¼ lb	Yellow
NEW Tactile Switch Set #3111		Auditory click, Tactile	40 g	Pressure	3"	4½" x 2½" ½ lb	Green/Blue/ White - Color may vary
Textured Saucer Switches (Set of 3) #722		Auditory click, Tactile	25 g	Pressure	7"	6" 1¾ lbs	Orange/Blue, Blue/Yellow, Yellow/Red
Textured Roller Switch #639		Tactile, Visual, Lights, Vibration, Music	15 g	Pressure	N/A	11" x 4½" 1¾ lbs	Yellow/Green

	Theme & Animal Switches									
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR			
Credit Card Switch #5580		Auditory Click	60 g	Pressure	4¼" x 3"	4¼" x 3" x 1/8"	Black			
Fun Animal Switches #3157M,S,P,V & #3158		Auditory Click	60 g	Pressure, Pinch	Varies	Varies	Multi			
Lady Bug Switch w/Lights, Music & Vibration #729	The state of the s	Visual, Tactile, Lights, Vibration	50 g	Pressure	6" x 6¼" x 4	9" x 6¼" x 4 1 lb	Black/Red/ Yellow			
Multisensory Funny Face #743		Auditory click, Tactile, Visual, Lights		Pressure	7¼" x 4½ "	11" 1½ lbs	Yellow/Blue			

Wheelchair & Bedside Switches							
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Headswitch #705		Auditory click, Tactile	30 g	Pressure	3¼"	3¾" ½ lb	Black
Mother's Third Arm w/Switch #9098		Auditory click, Tactile	30 g	Pressure	2½"	24" x 2" 1 lb	Black/Yellow
Mother's Third Arm w/Pillow Switch #9090		Auditory click, Tactile	40 g	Pressure	3"	26" x 3¾" ¾ lb	Black & Green
Ultimate Switch #1671		Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2½ lbs	Black
Ultimate Switches #1551, 1552	3.8	Auditory click, Visual, Tactile	10 g	Pressure	4" Spring	6" 2 lbs 12" 2½ lbs	Black
Wobble Switch w/19" Gooseneck and Clamp #1063		Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2¼ lbs	Black/Yellow

	Wireless,	Sound	& Motion	Activate	d Switch	es	
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Movement Sensor Switch #655		Visual	N/A	Motion, Adjustable sensitivity	N/A	4½" x 2½" x 1" ½ lb	White
Sound Activated Switch #1051	S		N/A	Sound	N/A	4½" x 2½" x 1" ½ lb	White
Wireless Gumball Switch #786		Auditory click, Visual, Tactile	30 g	Pressure	2½"	4" x 5" x 2" ¼ lb	Red
Wireless Ultimate Switch #2700		Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2¾ lbs	Black

Switches for Computers and iPads								
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR	
Computer Mouse Interface #1165		Auditory click, Tactile, Visual	20.9	Switch Input Pressure	5"	5" ¾ lb	Yellow	
Pancake Switch with Wireless Computer Interface #5071		Auditory click, Tactile	40 g	Pressure	5 3/8"	7" ½ lb	Green/White	
iPad Wireless Switch #1164		Auditory click, Tactile, Visual	20.9	Switch Input Pressure	5"	5" ¾ lb	Yellow	
Switch Modifiers								
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR	
Switch Modifier #605		N/A	N/A	Switch Input Pressure		4½" x 2½" x 1" ¼ lb	White	

Switches Kits for Assessment								
DEVICE	РНОТО	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR	
Assessment Kit for More Involved Clients #7755		Varies	Varies	Pressure	Varies	Varies	Varies	
Basic Switch Assessment Kit #7750		Varies	Varies	Pressure	Varies	Varies	Varies	
Classroom Switch Kit #7803N		Auditory click, Tactile	Varies	Pressure	Varies	Varies	Varies	
For further information please see the classroom kit section of our website.								