

## ActivaPatch® IontoGo™ technology.

No external power source to connect, no leadwires,  
no dose controller needed, power shut down feature and  
**no saline required.**

# ACTIVApatch®

## IontoGo™ 4.0

Extended wear patch with an average drug delivery dose of 80 mA\*min in approximately 4 hours.



**Active area:** 11.4 cm<sup>2</sup>  
**Dosage:** 80 mA\*min  
**Skin Adhesive:** Hypoallergenic  
**Fill Volume:** 2.0cc  
**Packaging:** 6 IontoGo™ 4.0,  
6 Alcohol Wipes  
**Treatment Time:** 4 Hours

NC89189

**NO  
SALINE  
REQUIRED**

**FULLY  
SELF-  
CONTAINED**



# ACTIVApatch®

## IontoGo™ 12.0

Extended wear patch with an average drug delivery dose of 80 mA\*min in approximately 12 hours.



**Active area:** 11.4 cm<sup>2</sup>  
**Dosage:** 80 mA\*min  
**Skin Adhesive:** Hypoallergenic  
**Fill Volume:** 2.0cc  
**Packaging:** 6 IontoGo™ 12.0,  
6 Alcohol Wipes  
**Treatment Time:** 12 Hours

NC89191

**NO  
SALINE  
REQUIRED**

**FULLY  
SELF-  
CONTAINED**



### HOW SMART IS YOUR WIRELESS IONTO?

## Smarter Iontophoresis™