

Bariatric Gel-Foam Pad



Reorder #	Size	Unit
<input type="checkbox"/> 915110	20" x 18" x 3"	Each
<input type="checkbox"/> 915120	22" x 18" x 3"	Each
<input type="checkbox"/> 915130	24" x 18" x 3"	Each
<input type="checkbox"/> 915135	26" x 18" x 3"	Each
<input type="checkbox"/> 915133	26" x 20" x 3"	Each
<input type="checkbox"/> 915134	28" x 18" x 3"	Each
<input type="checkbox"/> 915137	28" x 20" x 3"	Each
<input type="checkbox"/> 915139	30" x 20" x 3"	Each

Purpose

Skil-Care's Bariatric Gel-Foam Pad features a surface layer of water-based gel that relieves interface pressure while cooling the skin to reduce perspiration. The foam base adds both comfort and stability. This cushion is specially constructed to prevent bottoming-out and maximize weight distribution.

Important Notice

The Bariatric Gel-Foam Pad's cloth cover has safety straps that must be secured to the wheelchair to prevent the cushion (and possibly the user) from slipping off the wheelchair seat.

Easy Care

The inner vinyl cover is heat sealed to prevent urine and other fluids from seeping into the cushion. The vinyl cover may be cleaned with a spray-type cleanser and disinfected with a non-chlorine disinfectant. The outer cloth cover may be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.