

PORTABLE CEILING LIFT SOLUTION

Maxi Sky[®] 440



Our complete range of Maxi Sky ceiling lifts is designed to integrate with your care environment, helping to ensure a comfortable and dignified experience for everyone involved. The Arjo ceiling lifts allow caregivers to transfer patients or residents without manual lifting, thereby helping to reduce the risk of injury for themselves and the person in their care.

- **Simple for the caregiver**

The Maxi Sky 440 is an easy-to-operate portable ceiling lift system, which allows a single caregiver to perform transfers of patients, avoiding the stress and strain associated with manual lifting.

- **Flexibility to meet your care setting needs**

Designed to be easy to move around and to quickly attach to a gantry or track, this portable ceiling lift can be set up wherever patient handling is required. Maxi Sky 440 and its integrated foldable two point spreader bar is also compatible with the large range of Arjo loop slings.

ARJO MOBILITY GALLERY



D



E

© Arjo 2019

D. Unable to bear weight through feet. Is able to sit if well supported.

E. Difficulty maintaining sitting, hence might remain lying/partially lying in bed, and may only sit in a special supportive chair.

Please contact Arjo for further information.

PRODUCT INFORMATION

| | |
|--|-------------------------------------|
| Max. safe working load (SWL) | 200 kg / 440 lb |
| Unit weight (batteries included) | 6.5 kg / 14.3 lb |
| Lifting speed | 4.3 cm/s / 1.7 in/s |
| Folding sling hooks | Reduced space for storage/transport |
| Electronic microprocessor | Soft-start and stop motor control |
| Maintenance required indicator | Visual LED |
| Lift - protection class - shock prevention | Internally powered equipment |
| Strap length | 2200 mm / 88 in |
| Optional scale | External scale available |
| Noise level for either raising or lowering, with or without load | 54 dBA max |

Batteries

| | |
|-----------------------|---|
| Battery | 2x2.3 Ah will average 2 to 4 hours to a complete charge |
| Battery capacity | Provides up to 30 transfers with a load of 80 kg (176 lb) |
| Low battery indicator | Audible |
| Visual LED indicators | Flashing yellow (charging). Solid yellow (charged). |

Hand control

| | |
|--|------|
| Low force activation buttons, polyester membrane | IPX7 |
| ABS fire retardant | V0 |

Charger unit

| | |
|------------------------------------|---------------------------|
| Power indicator on charging module | |
| 100-240 Vac/50-60 Hz/24 VA | |
| Battery charger safety protection | Class 2, double insulated |

Safety features

| | |
|---|-------------------------------|
| Electrical up and down emergency buttons | Located on the motor cab |
| Emergency stopping device (push button) accessible below unit | Stop functions |
| Emergency brake | In case of mechanical failure |
| Manual emergency lowering device | Accessible from the motor cab |
| ABS Fire retardant casing | V0 |

SPACE REQUIREMENT

- Blue area shows the minimum working area required for the staff to be able to use the mechanical aids in an ergonomic way from one side.
- Light blue area shows required extension of working area to facilitate activities from either side to provide adequate access.

Arjo slings and spreader bars are a vital part of our patient lifting systems. A variety of slings and spreader bars provides a range of different choices when lifting and transferring patients.

