


APB™ - All Purpose Boot

A Post-Op Solution for ALL Seasons

Standard Surgical Shoe



APB™



New EVA Sole!


The DARCO APB™ embodies nearly perfect 4 season postoperative protection with a closed-toe design that keeps dressings and toes warm and dry. Our boot is designed to provide unparalleled versatility by accommodating both bulky dressings and most fiberglass casts. Now with a softer, more durable EVA Sole.



DARCO



Larger toe box allows for more space, accommodating bulky bandaging and the PQ PegAssist™



NEW EVA Sole is:

- › Softer
- › More comfortable
- › Longer lasting
- › Provides less pressure & more secure footing

APB™ Indications:

- › Post Surgical healing
- › Unna Boots, casts and other bulky dressings
- › Foot/Ankle trauma
- › Wounds or ulcerations present under metatarsal heads and toes.

APB™ Features and Benefits:

- › **NEW EVA sole** is softer, providing less pressure on the foot and increased patient comfort. It's more durable, lasts longer, provides more secure footing, and is less weight than the previous TPR sole.
- › **Square Toe Design** acts as a bumper to protect the toes when k-wires are present and provides better Left/Right fit.
- › **Zoned Outsole** has an aggressive tread under the mid-foot where traction is needed most.
- › **Removable Insole** features twice the padding of standard insoles and can be modified as needed. Insole can be removed to accommodate DARCO's PQ Series PegAssist™ off-loading insole.
- › **Use with PQ PegAssist™** to provide excellent off-loading in a closed toe shoe.

APB™ Sizing Chart:

Estimated Shoe Sizing:

Size:	Women:	Men:	Product Codes:
XS	3 - 5.5	-	APQ0B
S	6 - 7.5	5 - 6.5	APQ1B
M	8 - 9	7 - 8	APQ2B
L	9.5 - 11	8.5 - 10	APQ3B
XL	11.5 - 13	10.5 - 12.5	APQ4B

