

LINNER



■ Your Hearing

■ Your Life

■ Your LINNER

Pioneer Of OTC Hearing Aids



LINNER, established in 2016, stands as a beacon of innovation to a life well heard. We're dedicated to transforming how you experience sound, via blending cutting-edge technology with sleek design.

LINNER

About LINNER – 3 Key Pillars

■ Hardware & Software & Audiology Professionalism


“Empowering Better Hearing through Innovative solutions”

Intelligent Hardware



"Linner Hearing Aid products integrate cutting-edge audio technology, tailored for our customers' needs."

Interactive Software



"Advanced digital functionality for tailored, effective customization to individual hearing needs."

Service



"Direct access to audiology experts for the answers customers need."

Market Size of Hearing Aids (OTC)


■ Marker Share of Different Style & Coverage of Different Age Group

In 2022, the total market value of OTC hearing aids was **1.06 billion US dollars**.


*Statistic from GrandViewResearch

The global OTC hearing aids market is projected to experience steady growth, with a **Compound Annual Growth Rate (CAGR) of 7%** during the forecast period from 2023 to 2031.

*Statistic from **RESEARCHANDMARKETS**
THE WORLD'S LARGEST MARKET RESEARCH STORE



- ordered online/direct
- purchased OTC
- prescribed



65 years old
33%




75 years old
45%



85 years old
55%

Market Potential Regarding Hearing Aids Coverage

■ Ranging from 30 - 90% is Unaided Population in the USA




LINNER Target Audience


■ Around 43 Million with Hearing Loss & Over 90% are Mild & Moderate Hearing Loss in the US

43 MILLION ADULTS WITH HEARING LOSS IN THE UNITED STATES

Hearing Loss Prevalence by Severity



Hearing Loss Prevalence by Severity and Age



LINNER Target Audience

■ Those Who Suffer From Hearing Loss (Mild, Moderate to Severe)

Target Audience:

PPL Suffer from Hearing Loss

Age:

40-65 Years Old (Main Focus)


Gender:

Male & Female & etc.


Degree of Hearing Loss:

From Mild, Moderate to Severe


What is your hearing loss level?



What is your age?



How often do you wear hearing aids on a daily basis?




Mars



LINNER Mars Hearing Aids

LINNER




Focus Mode
Pick Up Sound From A Specific Direction And "Tune Out" Noises From Others



Enhance Sound Clarity & Speech Understanding
18 Digital Channels Sound Processor



Universal Bluetooth Connectivity
TV & Mic And Phone



Reduced Feedback Noise For Improved Sound Quality
Adaptive Feedback Control



Cosmic Blue **Stellar White**

Aesthetic & Ergonomic Design




IPX4
Sweat & Water-Proof



Customized Fit
User Confidence




Great & Discreet Comfort
Compact & Sleek Design



Smart & Personalize Your Hearing Experience
LINNER HA APP



9 Hours of Use
Battery Life For **All-Day** Use



Noise Reduction
Enjoy Speech Precision Even In Noisy Environments
Filters Out Background Noise




Premium Ambient Sound
Dynamic Speaker



Mercury

LINNER Mercury Hearing Aids

LINNER



More Discreet
Ergonomic Design for Great Comfort of Long Time Use



3 modes & 9 selectable volume to meet your needs possible



Echo Shield to Help Stay Focus
Via Adaptive feedback control



Ocean Blue


Beige Gold

Aesthetic Ergonomic Design (Traditional)



20-hour

Battery life



ACTIVE NOISE CANCELLATION


Stay Focus in Noise
Noise Cancellation



For safety & Security
Sweat & Water-Proof



Recharge on-the-go
Fast & Wireless Charging



16-CHANNEL PROCESSOR

Enhance Sound Clarity
Via 16 digital channels processing



LINNER Nova/Lite Hearing Aids

LINNER




Enhance Sound Clarity & Speech Understanding
Via 16 Digital Channels Processing

Safety & Hygiene Assured

Via UV-C Antibacterial

UV-C



Echo Shield to Help Stay Focus

Via Adaptive Feedback Control




4 Sound Modes to meet your needs possible
Dialogue & Restaurant & Commute & Outdoor

nova



Pearl White

Aesthetic & Ergonomic Design



Universal Bluetooth Connectivity
TV & Mic And Phone



For safety & Security
IPX5 Sweat & Water-Proof



Stay Focus In Noise
Noise Cancellation

Smart & Personalize Your Hearing Experience
LINNER HA APP



Precise & Flexible Amplification Capabilities
Via Dynamic Speaker



Recharge On-the-go
Fast & Wireless Charging



Sound Therapy & Healing
Meditation & Deep Sleep & Relaxation


*Nova is with Connect, Lite is not



venus

LINNER Venus Hearing Aids

LINNER



12 Channels Channel Processor:
Precise Sound Processing.

Simplifying The User Experience



Adaptive Feedback Control
Minimizes Feedback Noise,
Enhancing Comfort and Sound Quality.

Mist White **Onyx Black**

← knowles

Knowles Speaker:
Enhanced Sound Quality and
Speech Intelligibility.

3

Hearing Modes

Eartips for Better Fit
Customizable Eartips Offer Improved
Comfort and a Secure Fit.

Mild to Moderate Degree of Hearing Loss
Experience Improved Hearing.




24 Hours

Battery Life:
Allows for
All-day Use.
































LINNER Saturn Hearing Aids

LINNER




LINNER Product Spec

	Solar Series Availability	Fitting Range	APP Compatibility	Bluetooth	Speaker	Channel	LINNER Connect	Auto Feedback Control	Noise Cancellation	Auto On	Waterproof	Charging Case	Battery Life	FSA & HSA	Style
	Mars 	Mild to Moderate	 LINNER HA		 Dynamic	 18	✓	✓	Dual Mic Beam-forming	IR Sensor	IPX5	 Wireless Charging Type-C Quick Charge	Up to 8 Hrs	Eligible	ITE
	Nova/Lite 	Mild to Moderate	 LINNER HA		 Dynamic	 18	✓	✓	Dual Mic	IR Sensor	IPX5	 Wireless Charging Type-C Quick Charge UVC-Sterilization	Up to 8 Hrs	Eligible	OTE
	Mercury 	Mild to Severe	—	—	 Armature	 16	✗	✓	Single Mic	Delayed Start	IPX2	 Type-C Quick Charge	Up to 20 Hrs	Eligible	RIC
	Venus 	Mild to Moderate	—	—	 Armature	 12	✗	✓	Single Mic	Delayed Start	IPX4	 Type-C Quick Charge	Up to 24 Hrs	Eligible	ITE
	Saturn 	Mild to Moderate	—	—	 Armature	 6	✗	✓	Single Mic	Delayed Start	IPX2	 Type-C Quick Charge	Up to 20 Hrs	Eligible	CIC


LINNER HA App

LINNER


Enhance, Customize, Simplify:



Enhance



Customize



Simplify




Sound Therapy



Hearing Test




Media Control



Hearing Aids Control



Online Customer Support



Advanced
Personalized Setting


Online Customer
Support

Hearing Aid Control

Media Control

Sound Therapy

Hearing Test



LINNER Connect

LINNER

Expand Convenience & Functionality to the Maximum
& MSRP: \$99.99



TV Mode



LINNER Connect



AUX 3.5mm Audio Cable

Wireless Mic Mode



-  Dual Modes For Versatility
-  Extended Sound Capture Range
-  Ultra-Low Latency Interaction
-  Convenient Rechargeability




Dual Modes for Versatility:



Extended Sound Capture Range:



Ultra-Low Latency Interaction:



Convenient Rechargeability:



*1 Feet = 0.3 Meters

LINER

Thanks