



Ergonomically designed  
to encourage self-mobility

*The Pedal Rocker is an exceptional  
mobility chair that provides a gentle rocking motion.*

#### STANDARD FEATURES

- BRODA Comfort Tension Seating®
- 16 Gauge Powder Coated Tubular Steel
- 16", 18" or 20" Seat Width
- 17" or 18.5" Seat Depth
- 32" Back Height
- 14" to 19" Adjustable Seat Height (1" increments)
- Height Adjustable Armrests
- 18° of Infinitely Adjustable Tilt (up to 33° of tilt when combined with rocking)
- Rocking Function with Caregiver Lock
- Adjustable Rocking Tension
- 5" Front Total-Lock Casters
- 5" Rear Fixed Casters with Brakes
- 250lb Weight Capacity

#### STANDARD PADDING INCLUDES:

- Neckrest Pad
- Lower Side Pads

## PEDAL ROCKER MODEL: 48R

*When Comfort is Essential*

- Residents can self propel and rock in the same chair! Excellent for those who are mobile that may be agitated and have a tendency to rock.
- Rocking motion helps to calm individuals and helps prevent wandering - Excellent for Alzheimers. Rocking also helps absorb rigid movements to aid in resident safety and comfort.
- 18° of infinitely adjustable tilt in space cradles the body helping to prevent slumping, sliding and falls without the use of restraints. It also provides exceptional comfort while helping to reduce pressure.
- Front pivot seat tilt allows for a proper foot on floor position for foot propelling while tilted - Encourages safe mobility and independence.
- Adjustable seat to floor height from 14" to 19" accommodates both short and tall residents.



Pedal Rocker with full padding package.



Resident can use rocking  
function independently of tilt  
(rocking function shown)

# PEDAL ROCKER

## MODEL: 48R

*When Comfort is Essential*

### OPTIONS

- Full padding package including neckrest, wings, shoulder bolsters, lower side pads, seat and back with removable terry or vinyl covers. Full vinyl covers also available (right) as additional option.
- 20" or 22" front or rear mounted mag wheels
- 22" seat width
- 15.5" or 20" seat depth
- 33.5", 35" back length

### COMMON ACCESSORIES

- Swing away or elevating swing away footrests
- Solid or clear tray
- Padded tray cover
- Half tray
- Thigh belt
- IV Pole
- O2 Holder
- Molded laterals

The ergonomic design maintains a proper foot on floor position when using the seat tilt, helping to prevent falls.



18° of seat tilt



Full seat tilt combined with rocking 33°



# BRODA®



Shown with optional mag wheels, clear tray and full padding with full upgraded vinyl seat and back.

*When not self propelling, residents can use the tilt function along with the rocking function to achieve a deep seat tilt position. Excellent for napping!*

### BRODA's Comfort Tension Seating®

*The Key to Comfort is in BRODA's Strapping!* Memory retention allows straps to return to original shape within seconds.

- Redistributes pressure helping to prevent skin breakdown
- Creates air flow for heat and moisture exchange
- Contours to the body for long term sitting comfort

