Protection through Containment"

CareBag®

## **TOILET BOWL LINER**

### **With Super Absorbent Pad**



The single-use CareBag® Toilet Bowl Liner with super absorbent pad allows regular toilets to be used as dry toilets, with no flushing required. The Toilet Bowl Liner is recommended for disposal of urine or stools contaminated with biocides, radioactive or chemotherapy substances. It is also ideal in situations where specimens need to be collected for analysis.

The CareBag® Toilet Bowl Liner with super absorbent pad is designed to safely collect, transport, and dispose of body fluids. The super absorbent pad turns liquid waste (urine, diarrhea, blood, serous fluid) into a gel within seconds.

Placed in a bucket, the CareBag® Toilet Bowl Liner is ideal for safe disposal of used specimens from laboratory analysis or as a liner for the kick bucket in the Operating Room or Delivery Room.

- Protects all types of toilet bowls and buckets
- Hygienic, single-use product
- Safe and convenient collection, transport, and disosal of body fluids
- Reduces risks and costs associated with HAI
- Quickly absorbs liquids and semi-solid materials
- Eliminates odors when completely closed
- Saves time in routine cleaning and decontamination of toilet bowls by housekeeping or healthcare personnel

#### Ideal for Use:

- When plumbing/sewage is not available
- With patients undergoing radioactive seeding or therapy
- In homecare for chemotherapy treatments when elimination should not be done in any water system
- Also useful as a kick bucket liner in the Operating Room and Delivery Room.



# CAREBAG® TOILET BOWL LINER

## With Super Absorbent Pad

#### **Product Description**

The CareBag® Toilet Bowl Liner is placed directly over the bowl and under the seat of the toilet, thus providing a clean, dry toilet. The CareBag® Toilet Bowl Liner can also be used to protect any kind of basin, especially in Operating or Delivery Rooms. It is also useful when there are sewage or other water management issues, or when water is simply not available. The super absorbent pad turns body fluids into a gel.

#### **Benefits**

- · Protects all types of toilet bowls and basins
- · Quickly absorbs body fluids and semi-solid materials
- Prevents the dispersion of dangerous or infectious substances into the environment
- Single-use, prevents cross-contamination
- Hygienic, eliminates odors when completely closed
- Minimizes risk of cross-contamination

#### Instructions for Use

- 1. Open box, unroll the first liner and tear along the perforated line.
- 2. Open the liner and position it over the rim of the toilet how!
- 3. Place GelMax® super absorbent pad at or near the center of the toilet bowl.
- After patient uses toilet bowl liner, body fluids turn to gel within seconds.
- 5. To close: Pull ties tightly and knot the ties.

#### Disposal: like diapers

- If patient is not infectious, dispose of closed liner with regular waste.
- Discard infectious waste according to facility infection control guidelines.
- Do NOT flush.

#### **Precautions**

Single-use product. Do NOT reuse (2)
Keep out of reach of children
DO NOT FLUSH the CareBag® Toilet Bowl Liner

#### **Product Composition**

Bag: Polyethylene liner with colored straps Pad: Non-woven composite of pulp and acrylate polymer

#### **Absorption Capacity**

Absorbs up to 16 oz (450 ml) of organic liquids

#### Identification & Expiration Date

Each box is identified by a reference/lot number and expiration date [5 years]

#### Safety Data

- This product is not considered dangerous according to the European Union regulation
- MSDS sheets available upon request

#### Conformity

- Compliant with the FDA regulation Class I, 510(k) exempt medical device
- Compliant with the 93/42/CEE regulation for medical device – Class I
- Eco-friendly: Compliant with anti-pollution legislation on incineration risks. During incineration, the product does not release any chloride element or byproduct that might create dioxin.
- Customs code: 5603149090
- Company certification: ISO 9001:2008 and ISO 13485:2003

Reorder Number: 9538534

Unit	Box of 20 liners	Liner: 25 <sup>9/16</sup> " x 15 <sup>3/4</sup> "; Pad: 5 <sup>1/8</sup> " x 9 <sup>7/16</sup> "
Case	18 boxes of 20 liners (total of 360 liners)	16 <sup>1/4</sup> " (L) x 14 <sup>1/8</sup> " (W) x 14 <sup>7/16</sup> " (H); 30.86 lbs
Pallet	30 cases (540 boxes, total of 10,800 liners)	31 <sup>1/2</sup> " (L) x 47 <sup>1/4</sup> " (W) x 77 <sup>15/16</sup> " (H); 970 lbs

