BRODA

Comfort

Rehab Wheelchair



Shown with optional head support, Matrx Elite back support, Matrx Vi seat cushion, and angle-adjustable lower leg supports

Patient, Facility, & Caregiver Benefits

- Pressure redistribution that aids in maintaining skin integrity
- Adaptable for wide range of patients and care needs
- Durable frame to withstand heavy workload healthcare environments



Front pivot seat tilt maintains user comfort with smooth, effortless tilt



Posterior tilt optimizes pressure redistribution



Solid seat pan offers support and is compatible with a variety of seat cushion options



Range of simple adjustments accommodates a variety of needs



Comfort

Rehab Wheelchair

SIMPLIFIED DURABILITY FOR REHABILITATION CARE

SPECIFICATIONS

Degrees of Tilt	<i>Up to 45°</i>
Weight Capacity	Up to 250 lbs. (114 kg)
Angle Adjustable Back Canes	5°, 10°, 15°, 20°
Seat Widths	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths	Adjustable: 16", 17", 18", 19", 20 (40cm, 43cm, 45cm, 48cm, 50cm)
Seat to Surface Heights	Adjustable: 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19" (35cm, 38cm, 39cm, 42cm, 43cm, 44cm, 45cm, 48cm)
Wheels	5", 6" or 8" Front Caster Wheels with fork (13cm, 15cm, 20cm) 20", 22" or 24" Rear Mag Wheels (50cm, 55cm, 60cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Government Coding (US & CA)	HCPCS: E1161 ADP (ON): WMA5B-025