Plastic Cylinder Wrenches (E-Keys)

25/pk Also known as "E Keys".	MED
Used on Post Valve Oxygen Cylinder	
to Open and Close.	

Oxygen Regulator Connector Swivel

1

EK-1

(XMAS Trees), 25/pk For Use on Concentrator or Regulator with DISS Fitting in Order to Attach Cannula.



E-Key

Cylinder Wrench

18501

For Use with Small Cylinders used in Most Portable Oxygen Set Ups, 10/pkg



- Has place to store Viton rubber/ brass sealing washer.
- Used to open post valves.
- Can be used to assist in install or removal of both large wing nut and small t-bar regulators.

Concentrator Water Trap 25/bx

CON 900

FEATURES:

Easy inline downstream installation with small bore tubing. Can be cleaned with detergent or in a vinegar bath.



18501

Tubing Extension Connector

Swivel, 50/bg
Straight, 50/bg
Ridged, 50/bg
Y Connector, 10/bg

FEATURES:

Rigid plastic tubing connector allows extension tubing to be connected to present cannula for greater length(s).



COMPARE TO:		
	Hudson®	Salter®
CON 400	N/A	1220
CON 500	N/A	2005
CON 550	1420	1215
XM-1	2555	1205