

Prevent Falls & Wandering

Corded Pads & Mats For Bed & Chair Exit Alarms



Works with: Corded Smart Fall Monitors

MFR#	SMART PRODUCT DESCRIPTION
PPB-WI	Antimicrobial Bed Pad 20" x 30" (1 year pad life)
PPB-90W, PPB-45W	Antimicrobial Bed Pad 20" x 30" (90 day, 45 day pad life)
PPB-RI	Antimicrobial Bed Pad 10" x 30" (1 year pad life)
PPB-90R, PPB-45R	Antimicrobial Bed Pad 10" x 30" (90 day, 45 day pad life)
PPC-WI	Antimicrobial Chair Pad 10" x 15" (1 year pad life)
PPC-90, PPC-45	Antimicrobial Chair Pad 10" x 15" (90 day, 45 day pad life)
FM-03	Floor Mat 24" x 36" Black (1 year warranty)
FM-05	Floor Mat 24" x 36" Gray (1 year warranty)
FM-06	Floor Mat 20" x 30" Gray (1 year warranty, for under rug use)
FM-07	Floor Mat 24" x 48" Gray (1 year warranty)
WC-02	Weight Sensing Fall Mattress 70.5" x 29.5" x 1.5" (90 Day warranty)
LM-01	Corded Landing Mat 24" x 71" x 1", Grey (1 year warranty)

Simple early warning system to help prevent falls.

Bed and chair sensor pads alert caregivers whenever pressure is removed from the sensor pad (i.e. when a resident gets up from a bed or chair and pressure is removed from the pad).

Floor Mats & Landing Mats alert caregivers when pressure is applied to the sensor mat (i.e. when a resident steps onto a floor mat and pressure is placed on the mat).

Floor mats are ideal for bedside and doorway mobility monitoring. Bed and chair pads are ideal early warning alerts of resident mobility as they exit beds and chairs.

