

SHOULDER P.A.D. II

**Pillow
Abduction
Device**



 DeRoyal

 DeRoyal

SHOULDER P.A.D. II

Pillow Abduction Device



The Shoulder P.A.D. II provides immobilization for surgical repairs of the rotator cuff, anterior glenohumeral capsular reconstruction, dislocation/subluxation of G-H joint, S/P glenohumeral acromioplasty, Bankart repairs, soft tissue repairs, and shoulder sprains and strains.

SMALL	11630005
MEDIUM	11630006
LARGE	11630007

*Measure elbow crease to wrist

Features & Benefits

- Allows the arm to be placed in an optimal, neutral position for effective recovery
- Arm may be positioned more anteriorly or posteriorly depending on the post-operative requirements
- Orthowick™ liner for greater comfort and breathability
- Sling pouch may be closed with top strap to secure arm
- Adjustable waist and shoulder straps
- Thumb loop to prevent wrist drop
- Stress ball provided to encourage exercise and proper circulation
- Universal left or right



200 DeBusk Lane
Powell, TN 37849 U.S.A.
888.938.7828
865.938.7828
www.deroyal.com