



Molift UnoSling

Patient specific
for infection control

 **etac**[®]

Molift UnoSling HighBack

Molift UnoSling HighBack is suitable for most users in common lifting situations. The sling is a patient specific sling, intended for an individual user and for short-term use. An effective solution to prevent transmission of infection.

High requirements for infection control

The Molift UnoSling HighBack is a personal disposable product designed to meet the demands and challenges of today's hospital environment. Slings for short-term use are specific to a single user, and an effective solution for preventing the transmission of infection.


The UnoSling has a high back and is suitable for most users in the most common hoisting situations, such as between bed and wheelchair or to the toilet. The sling remains with the user during the entire length of their stay in hospital.

Sling for short-term use

The UnoSling is manufactured in a strong non-woven material which is not washable. The sling is discarded if it becomes dirty or damaged, or if the user no longer needs it.

- A personal disposable sling for short-term use
- Hoisting to and from lying or sitting position
- Support for the entire back and the head
- For environments with a short term need for a sling. Can be used several times by the same user.



Variant	Pcs/pack	User weight	Item no.	Disposable
S	20	25-50	3060020	Do not wash. The UnoSling should be discarded when the user no longer needs it, or if it becomes damaged or dirty.  max 300 kg max 660 lbs
M	20	45-95	3060030	
L	20	90-165	3060040	
XL	20	160-240	3060050	
XXL	10	230-300	3060060	



Features



Straps with 2 heights
Straps with 2 different height settings.

2-point or 4-point
The Molift UnoSling is suitable for use with either 2-point or 4-point sling bar.

Comfortable support
Support for the entire back and the head.


Label with user's name
At the back: large, clear label with practical information. The size is shown clearly and it is easy to add the user's name, which ensures safe use.

Disposable product
All-round sling that fits most users.

Accessories



Molift RgoSling Extension Loops
For easier adjustment. Used to extend the hoisting straps with 3 different loop positions. Order 2 pcs for a pair.



Molift EvoSling Soft Padding
For pressure-sensitive users
EvoSling Soft Padding consists of pressure-relieving padded sleeves which are slipped over the leg parts of slings with split-leg support. Washable at 40°C.



Disposable product supplied in packaging which makes it easy to store and take out a new sling when needed.



Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.