

The ENHANCER DRY FLOATATION Cushion offers you a unique option in wheelchair seating. By offering a contoured form with a two manifold system, the client/user can benefit from positive midline channeling of the legs along with increased lateral stability. In addition, the special design of the ENHANCER cushion inhibits sacral sitting while enhancing tissue viability.

The ENHANCER's unique contoured shape is the result of flexible air cells of varying heights. The two manifold system allows you to: 1) adjust the interior bucket seat section for reduced pressure on the ischia and coccyx while 2) adjusting the two outer perimeter rows for enhanced stability and positioning. The two design features combined offer you a cushion that shapes to the user's body providing superior pressure distribution and minimal shear forces that are characteristic of every ROHO[®] DRY FLOATATION cushion, while offering you the added benefit of progressive position.





nexus®

AirLITE™

Features and Benefits

Variable Cell Heights: Provides for channeling of the legs and enhances lateral stability and body positioning.

Washable, Two-Way Stretch Cover: Allows for proper immersion into the cushion and "breathes" with every shift in position.

Adjustable: Shapes to the individual and to their specific seating needs.

Ventilation: Individual air cells allow for air circulation between the cells.



Lightweight: The ENHANCER is approximately four (4) pounds making it easy to transport and handle and minimizes fatigue of the individual. **Durable Neoprene Rubber:** Flame resistant and resistant to urine, oils and medicines. Easily cleaned with soap or water and can be steam autoclaved. Product can be repaired extending the overall life of the cushion.

Dual-Manifold Design: Permits specific adjustments to allow for reduced pressures on boney prominences while independently promoting better positioning and increases stability.

24 Month Warranty

Specifications

The ENHANCER comes complete with hand inflation pump, repair kit and operating instructions.

Construction–flame resistant neoprene rubber Weight–approximately four (4) pounds Cleaning–soap and water, or can be steam autoclaved

MODEL CHAIR SIZE fits 15"-16" widths & 15"-16" depths ENH88C fits 15"-16" widths & 16"-17" depths ENH89C fits 16"-17" widths & 15"-16" depths ENH98C fits 16"-17" widths & 16"-17" depths ENH99C ENH910C fits 16"-17" widths & 18"-19" depths fits 18"-19" widths & 16"-17" depths ENH109C ENH1010C fits 18"-19" widths & 18"-19" depths Can be custom made to fit a variety of seating surfaces and sizes.

Applications

Most effective for users:

- at moderate risk for skin breakdown needing some degree of abduction/adduction
 - with poor sitting balance
 who sacral sit

