

Made for

# Motion

Foot & Ankle Product Guide



DeRoyal®



**“The greater the difficulty, the more  
the glory in *surmounting it.*”**

**— Epicurus**



## Elastic Ankle Sleeve

- Pull on cotton elastic design
- Provides snug fit for maximum comfort and compression
- Heel cut-out for patient comfort
- Low profile design fits easily inside shoe
- Universal right/left design saves shelf space and inventory costs

**Indications:** Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4004-01
Medium	8" - 10"	4004-02
Large	10" - 12"	4004-03
X-Large	12" - 14"	4004-04

\*Measure circumference from the base of heel around the ankle

## Closed Heel Elastic Ankle Sleeve

- Closed heel, elastic construction
- Provides snug fit for maximum comfort and compression
- Pull on design for ease of application
- Low profile design fits easily in shoe
- Universal right/left design saves shelf space and inventory costs

**Indications:** Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4005-01
Medium	8" - 10"	4005-02
Large	10" - 12"	4005-03
X-Large	12" - 14"	4005-04

\*Measure circumference from the base of heel around the ankle



## Neoprene Ankle Support


- 1/8" neoprene sleeve provides warmth, compression and support
- Increases proprioception to the injured ankle
- Open heel construction for patient comfort
- Universal right/left design saves shelf space and inventory costs

**Indications:** Mild ankle sprains; Chronic ankle instability; Osteoarthritis; Ankle edema; Joint effusion

Size	Men's Shoe Size	Women's Shoe Size	Product #
Small	4 1/2 - 6	5 1/2 - 7	NE7732-72
Medium	6 1/2 - 8	7 1/2 - 9	NE7732-73
Large	8 1/2 - 10	9 1/2 - 11	NE7732-74
X-Large	10 1/2 - 12	11 1/2 - 13	NE7732-75
XX-Large	12 1/2 +	13 1/2 +	NE7732-76

HCPCS Code: L1901





“Sometimes being pushed to the wall  
gives you the *momentum*  
necessary to get over it.”

— Peter de Jager



HCPCS Code: L1901

## Figure 8 Wrap Ankle Support

- Constructed of a durable, breathable mesh sock material
- Elastic sleeve provides compression and edema control
- Figure 8 wrap with hook and loop closure provides additional compression and support
- Universal right/left design saves shelf space and inventory costs

**Indications:** Mild ankle sprains and strains; Ankle edema; Peroneal and anterior tibialis tendonitis; Osteoarthritis of the talocrural and subtalar joints

Size	Ankle Circumference*	Product #
Small	7" - 8"	4016-01
Medium	8" - 10"	4016-02
Large	10" - 12"	4016-03
X-Large	12" - 14"	4016-04

\*Measure circumference from the base of heel around the ankle

## Canvas Lace-Up Ankle Brace

- Made of durable canvas with flannel lining for added support and comfort
- Anatomically designed with contours to accommodate the ankle
- Removable medial and lateral stays provide stability
- Lace-up design for a secure fit
- Universal right/left design saves shelf space and inventory costs

**Indications:** Acute and chronic ankle sprains; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Small	7" - 8"	8761-01
Medium	8" - 10"	8761-02
Large	10" - 12"	8761-03
X-Large	12" - 14"	8761-04

\*Measure circumference from the base of heel around the ankle



HCPCS Code: L1902



HCPCS Code: L1902

## R.I.C.E. Ankle Splint

- Constructed of perforated vinyl with foam lining
- Removable stays with malleoli cut-out for comfortable support
- Lace-up design with figure 8 strap for a secure fit
- Universal right/left design saves shelf space and inventory costs

**Indications:** Acute and chronic ankle sprains; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Child	6" - 7"	4050-00
Small	7" - 8"	4050-01
Medium	8" - 10"	4050-02
Large	10" - 12"	4050-03
X-Large	12" - 14"	4050-04

\*Measure circumference from the base of heel around the ankle





*“Adversity causes some men to break;  
others to **break records.**”*

— William A. Ward



HCPCS Code: L1902

## Leatherette Ankle Brace

- Constructed of durable leatherette material with flannel lining
- Contoured design for best anatomical fit
- Removable medial and lateral metal stays provide stability
- Loop lock closure for compression and support
- Universal right/left design saves shelf space and inventory costs

**Indications:** Mild to moderate ankle sprains and strains; Ankle edema; Ankle rehabilitation

Size	Ankle Circumference*	Product #
Small	7" - 8"	4013-01
Medium	8" - 10"	4013-02
Large	10" - 12"	4013-03
X-Large	12" - 14"	4013-04

\*Measure circumference from the base of heel around the ankle

## DeRoyal® Sports Ankle Brace

- Constructed of durable nylon for long lasting wear
- Figure 8 strapping configuration, in conjunction with medial to lateral stirrup straps, places foot in slightly everted position for maximum ankle stability
- Removable medial and lateral plastic stays with malleoli cut-out provides comfort and support
- Lace-up design for a snug fit
- Low profile design allows use in athletic shoe

**Indications:** Chronic ankle instability; Prevention or treatment of acute ankle sprains; Prophylactic talocrural support; Subtalar joint sprains; Syndesmosis sprain

Size	Men's Shoe Size	Women's Shoe Size	Left	Right
X-Small	4 - 6	5 - 7	AB2801-10	AB2801-20
Small	6 1/2 - 8	7 1/2 - 9	AB2801-12	AB2801-22
Medium	8 1/2 - 10	9 1/2 - 11	AB2801-14	AB2801-24
Large	10 1/2 - 12	11 1/2 - 13	AB2801-16	AB2801-26
X-Large	12 1/2 +	13 1/2 +	AB2801-18	AB2801-28



HCPCS Code: L1902

## DeRoyal® Sports Ankle Brace 2

- Constructed of durable nylon for long lasting wear
- Universal stabilization strap forces the ankle into inversion or eversion
- Cuff closure for added support and compression
- Lace-up design for a snug fit and primary stabilization
- Low profile design allows use in athletic shoe
- Universal right/left design saves shelf space and inventory costs

**Indications:** Chronic ankle instability; Prevention or treatment of acute ankle sprains; Prophylactic talocrural support; Subtalar joint sprains; Syndesmosis sprain


Size	Ankle Circumference*	Product #
X-Small	10" - 11"	15600004
Small	11" - 12"	15600005
Medium	12" - 13"	15600006
Large	13" - 14"	15600007
X-Large	14" - 15"	15600008
XX-Large	15" - 16"	15600009

\*Measure circumference from the base of heel around the ankle



HCPCS Code: L1902



A close-up photograph of a basketball player's lower legs and feet. The player is wearing maroon shorts with white piping, white socks, and white sneakers with a visible air cushion in the sole. A basketball is suspended in the air to the right of the player's feet. The background is a green chain-link fence. The bottom of the image features a decorative wavy border in yellow, orange, and maroon.

“The great thing in the world is not so much where we stand, as in what direction we are *moving*.”

— Oliver Wendell Holmes





HCPCS Code: L1906

## DeRoyal® Functional Ankle Brace

- Anatomically designed uprights for an intimate fit
- Hinged stirrup to maintain support throughout range of motion
- Padded foam inlay for high energy absorption during activity
- Polyurethane foam padding for patient comfort and shock absorption
- Vel-stretch® straps for maximum support and edema control
- Universal right/left design saves shelf space and inventory costs

**Indications:** Acute and chronic ankle sprains; Chronic instability; Prophylactic talocrural support

Size	Men's Shoe Size	Women's Shoe Size	Product #
X-Small	—	5 - 7	AB2000-10
Small	5 - 8	7 1/2 - 10	AB2000-12
Medium	8 1/2 - 12	10 1/2 - 13	AB2000-14
Large	12 1/2 - 16	13 1/2 +	AB2000-16

## Element™ Ankle Brace

- Patented heel control strapping system holds the calcaneus under the talus controlling subtalar inversion and talar rotation
- Rigid medial and lateral uprights prevent abnormal eversion and inversion and provide protection during activity
- Functional hinge design allows for utilization of the brace prophylactically, after acute injuries and during rehabilitation
- Rear foot stabilization is enhanced with heat moldable footplate to control pronation during loading and mid-stance phase of the gait
- Pulley system strapping design provides optimal compression enhancing proprioception
- Wide base of support at impact reduces ankle inversion
- Unique design provides syndesmosis stabilization
- Concise version is 1 1/2" shorter for lower profile fit
- Available in standard and wide calf cuff

**Indications:** Acute and chronic ankle sprains; Chronic ankle instability; Syndesmosis sprain; Posterior tibial tendonitis; Peroneal tendonitis; Osteoarthritis

Size	Men's Shoe Size	Women's Shoe Size	Standard		Standard Concise		Wide	
			Right	Left	Right	Left	Right	Left
Small	< 8	< 9 1/2	AB2110-10	AB2110-16	AB2130-10	AB2130-16	AB2120-10	AB2120-16
Medium	8 - 12	9 1/2 - 13 1/2	AB2110-12	AB2110-18	AB2130-12	AB2130-18	AB2120-12	AB2120-18
Large	12 +	13 1/2 +	AB2110-14	AB2110-20	AB2130-14	AB2130-20	AB2120-14	AB2120-20

\*Measure calf circumference 8" from ground  
Standard = < 10"; Wide = > 10"



HCPCS Code: L1971

A person wearing a white t-shirt, a denim skirt, and a white and blue leg cast is playing mini-golf. They are holding a red and black mini-golf club, positioned over a blue ball on a green artificial turf surface. In the background, another person is visible, blurred, and the scene is set outdoors with trees and a building.

**“You can  
*promote your healing*  
by your thinking.”**

— James E. Sweeney





HCPCS Code: L4350

## Confor® Ankle Brace

- Confor® open cell polyurethane foam provides a superior fit and has proven to be one of the most conforming and shock absorbing materials on the market
- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Elasticized Vel-stretch® straps provide compression to help reduce edema
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily inside shoe
- Universal right/left design saves shelf space and inventory costs
- Available in each or 6/case quantity

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Each	Case
Standard	10"	AB2340-00	AB2340-06
Small	8 1/2"	AB2341-00	AB2341-06

## DeRoyal® Air Ankle Stirrup

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air bladder allows for simple adjustment of air with inflation bulb
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left design saves shelf space and inventory costs
- Kit features: standard ankle stirrup, ice pack, elastic bandage, exercise band and self-management guide

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2360-00
Small	8 1/2"	AB2361-00
Pediatric	6"	AB2362-00
Ankle Sprain Kit	10"	SC1000-00



HCPCS Code: L4350



HCPCS Code: L4350

## DeRoyal® Air/Gel Ankle Stirrup

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air/gel bladder can be placed in freezer for cold therapy
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left design saves space and inventory costs

**Indications:** Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2350-00
Small	8 1/2"	AB2351-00



*“There are many ways of  
going forward  
but only one way of standing still.”*

— Franklin D. Roosevelt





HCPCS Code: L4360

## ▶ Pacesetter™ Air Walker

- Inflatable air bladder design for compression and stabilization
- Unique non-welded air bladder reduces chance of over-inflation to the point of bursting or tearing
- Wide foot plate for patient comfort and stability
- Rocker sole allows for a rolling gait pattern
- Heel strap helps prevent foot migration

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Product #
Small	3 - 6 1/2	4 - 7 1/2	15780005
Medium/Large	7 - 10 1/2	8 - 12 1/2	15780067
X-Large	11 - 16	13 +	15780008

## ▶ Inline Pneumatic Walker

- Features inflatable air bladder wrap for compression and stabilization
- Heavy duty flared uprights are color-coded for easy size identification
- Low profile rocker sole promotes natural gait
- Wide foot plate for patient comfort and stability
- Treaded sole for patient safety
- Available in open and closed heel design
- Offered in both standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel	Low Top Closed Heel
Small	3 1/2 - 7	5 - 8 1/2	15470005	15490005	15480005	15500005
Medium/Large	7 1/2 - 10 1/2	9 - 12	15470067	15490067	15480067	15500067
X-Large	11 +	12 1/2 +	15470008	15490008	15480008	15500008
Pediatric	8 - 12 (Child)		15470002	15490002	15480002	15500002



HCPCS Code: L4360

## ▶ Pacesetter™ II Walker

- Rocker sole helps promote a rolling gait pattern
- Wide foot plate for patient comfort and stability
- Double laminated foam wrap is breathable
- Hook and loop straps for easy application and removal
- Available in open or closed heel design

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel
Small	3 - 6 1/2	4 - 7 1/2	CSG53510-A	CSG53550-A	CSG53520-A
Medium/Large	7 - 10 1/2	8 - 12 1/2	CSG53510-B	CSG53550-B	CSG53520-B
X-Large	11 - 16	13 +	CSG53510-C	CSG53550-C	CSG53520-C



HCPCS Code: L4386

A person is shown from the waist down, wearing khaki shorts and a grey sneaker on their left foot. Their right leg is in a black, multi-buckle medical cast. They are holding a fishing rod with a wooden handle and a silver reel. The background is a blurred body of water and trees. The bottom of the image has a decorative wavy border in shades of red and orange.

**“Effort only fully  
releases its rewards after a  
person *refuses to quit.*”**

— Napoleon Hill





HCPCS Code: L4386

## Inline Walker

- Heavy duty flared uprights are color-coded for easy size identification
- Low profile, rocker bottom promotes natural gait
- Flat, wide foot bed provides patient comfort and stability
- Treaded sole for patient safety
- Available in open or closed heel design
- Offered in both standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel	Low Top Closed Heel
Small	3 1/2 - 7	5 - 8 1/2	15520005	15530005	15540005	15550005
Medium/Large	7 1/2 - 10 1/2	9 - 12	15520067	15530067	15540067	15550067
X-Large	11 +	12 1/2 +	15520008	15530008	15540008	15550008
Pediatric	8 - 12 (Child)	—	15520002	15530002	15540002	15550002

## Tracker EX Walker

- Lightweight design for patient comfort
- Low profile rocker bottom sole helps promote natural gait
- Wide foot bed for patient comfort and stability
- Closed heel design helps prevent foot migration
- Extra long wrap provides full coverage and protection
- Available in standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard	Low Top
X-Small	up to 5 1/2	up to 4 1/2	15590002	15580002
Small	6 - 8 1/2	5 - 8 1/2	15590005	15580005
Medium	9 - 10 1/2	9 - 11 1/2	15590006	15580006
Large	11 - 13	12 - 13 1/2	15590067	15580067
X-Large	13 1/2 +	14 +	15590008	15580008



HCPCS Code: L4386



HCPCS Code: L4386

## Grid Iron XR® Walker

- Rigid, metal uprights for maximum stability
- Low profile rocker bottom sole helps promote natural gait
- Wide foot bed for patient comfort and stability
- Closed heel design helps prevent heel migration
- Available in standard and low top versions

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard	Low Top
X-Small	up to 3 1/2	up to 4 1/2	15570004	15560004
Small	4 - 7	5 - 8 1/2	15570005	15560005
Medium	7 1/2 - 10	9 - 11	15570006	15560006
Large	10 1/2 - 12	11 1/2 - 13	15570007	15560007
X-Large	12 1/2 +	13 1/2 +	15570008	—

A close-up photograph of a person's lower leg and foot wearing a blue, multi-strapped medical brace. The brace is made of a dark blue material with several horizontal straps and buckles. The person is wearing blue jeans. The leg is standing on a patch of green grass. To the right, a golden retriever is partially visible, looking towards the camera. The background is a blurred green field under bright sunlight. The overall scene conveys a sense of recovery and companionship.

*“Hope*  
is one of the principle springs that  
keeps mankind in motion.”

— Thomas Fuller





HCPCS Code: L4386

## Three-D® Lowest Profile Walker

- Extra low heel height helps prevent external hip rotation
- Tri-radiate sole encourages natural gait pattern
- Offset uprights at the malleoli accommodate swelling
- Padded foam inlay for increased energy absorption
- Flat foot bed keeps foot in neutral position
- Criss-cross straps help eliminate forward foot migration

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Product #
X-Small	up to 5	up to 6	SL5002-01
Small	5 - 7	6 - 8	SL5002-03
Medium	7 - 10	8 - 11	SL5002-05
Large	10 - 12	11 - 13	SL5002-07
X-Large	12 +	13 +	SL5002-09

## Three-D® Low Profile Walker

- Low heel height helps prevent external hip rotation
- Tri-radiate sole encourages a natural gait pattern
- Offset uprights at the malleoli accommodate swelling
- Two layer memory foam inlay conforms to foot
- Flat foot bed keeps foot in a neutral position
- Criss-cross straps help to eliminate forward foot migration

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Boot Wrap	Open Wrap
X-Small	up to 5	up to 6	SL2001-01	SL2001-08
Small	5 - 7	6 - 8	SL2001-03	SL2001-10
Medium	7 - 10	8 - 11	SL2001-05	SL2001-12
Large	10 - 12	11 - 13	SL2001-07	SL2001-14
X-Large	12 +	13 +	SL2001-09	SL2001-16



HCPCS Code: L4386



HCPCS Code: L4386

## One-Piece Samson Walker

- Low profile anatomical design promotes natural gait pattern
- Two layer memory foam inlay conforms to foot
- Adjustable straps for even compression
- Wrap material allows skin to breathe and helps to prevent skin breakdown
- Available in standard, mid and low top designs

**Indications:** Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Standard	Mid	Low Top
X-Small	up to 5	up to 6	SL8001-01	—	—
Small	5 - 7	6 - 8	SL8001-03	SL8001-43	SL8001-33
Medium	7 - 10	8 - 11	SL8001-05	SL8001-45	SL8001-35
Large	10 - 12	11 - 13	SL8001-07	SL8001-47	SL8001-37
X-Large	12 +	13 +	SL8001-09	—	—



*“Inspire yourself to dig deep  
and give the world your very best”*

— Laura Pickford Ramirez





HCPCS Code: L4386

## Grid Iron XR® R.O.M. Walker

- Range of motion is adjustable from fixed to 45°
- Rigid metal uprights for maximum stability
- Low profile rocker bottom sole helps promote natural gait
- Wide foot bed for patient comfort and stability
- Closed heel design helps prevent foot migration

**Indications:** Achilles tendon repair; Moderate to severe ankle sprains and strains; Stable foot and ankle fractures; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Product #
X-Small	up to 3 1/2	up to 4 1/2	15510004
Small	4 - 7	5 - 8 1/2	15510005
Medium	7 1/2 - 10	9 - 11	15510006
Large	10 1/2 - 12	11 1/2 - 13	15510007
X-Large	12 1/2 +	13 1/2 +	15510008

## Articulating Walker

- Adjustable hinge at ankle joint allows for positioning in 10° increments
- Low heel height helps prevent external hip rotation
- Tri-radiate sole promotes natural gait pattern
- Offset uprights at the malleoli accommodate swelling
- Two layer memory foam inlay conforms to patient's foot

**Indications:** Achilles tendon repair; Moderate to severe ankle sprains and strains; Stable foot and ankle fractures; Post ORIF

Size	Men's Shoe Size	Women's Shoe Size	Product #
Small	5 - 7	6 - 8	SL1201-10
Medium	7 - 10	8 - 11	SL1201-12
Large	10 - 12	11 - 13	SL1201-14



HCPCS Code: L4386



## Original Darby Shoe

- Semi-rigid metatarsal shank provides controlled forefoot motion
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Men's Shoe Size	Product #	Women's Shoe Size	Product #
Small	6 - 8	DS1014-10	4 - 6	DS1013-10
Medium	8 <sup>1</sup> / <sub>2</sub> - 10	DS1014-12	6 <sup>1</sup> / <sub>2</sub> - 8	DS1013-12
Large	10 <sup>1</sup> / <sub>2</sub> - 12	DS1014-14	8 <sup>1</sup> / <sub>2</sub> - 10	DS1013-14
X-Large	12 <sup>1</sup> / <sub>2</sub> - 14	DS1014-16	—	—
Pediatric	12 - 1 (Child)	DS1014-00	—	—

## Aegis Post-Op Shoe

- Semi-rigid rubber sole for support
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Men's Shoe Size	Product #	Women's Shoe Size	Product #
Small	6 - 8	2044-01	4 - 6	2044-07
Medium	8 <sup>1</sup> / <sub>2</sub> - 10	2044-02	6 <sup>1</sup> / <sub>2</sub> - 8	2044-08
Large	10 <sup>1</sup> / <sub>2</sub> - 12	2044-03	8 <sup>1</sup> / <sub>2</sub> - 10	2044-09
X-Large	12 <sup>1</sup> / <sub>2</sub> - 14	2044-04	—	—



HCPCS Code: L3260  
(Not Covered by Medicare)



## Deluxe Wooden Post-Op Shoe

- Rigid wood sole with rubber bottom for support
- Perforated suede upper for comfort and durability
- Available in hook and loop or lace up closure design
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Men's Shoe Size	Hook/Loop Closure	Lace Up Closure	Women's Shoe Size	Hook/Loop Closure	Lace Up Closure
Child	2 - 4 <sup>1</sup> / <sub>2</sub>	2007-01	2006-01	2 - 3 <sup>1</sup> / <sub>2</sub>	2007-06	2006-06
Small	5 - 6 <sup>1</sup> / <sub>2</sub>	2007-02	2006-02	4 - 6	2007-07	2006-07
Medium	7 - 8 <sup>1</sup> / <sub>2</sub>	2007-03	2006-03	6 <sup>1</sup> / <sub>2</sub> - 8	2007-08	2006-08
Large	9 - 11	2007-04	2006-04	8 <sup>1</sup> / <sub>2</sub> - 10	2007-09	2006-09
X-Large	11 <sup>1</sup> / <sub>2</sub> - 13	2007-05	2006-05	—	—	—

HCPCS Code: L3260  
(Not Covered by Medicare)





HCPCS Code: L3260  
(Not Covered by Medicare)

## Darco® Med-Surg™ Shoe

- MetaShank™ protection provides rigid control under metatarsal heads
- Strapless forefoot closure expands to accommodate bulky dressings
- Square toe design serves as a bumper for additional protection
- Rocker sole reduces plantar pressure on forefoot and heel
- Ankle strap helps eliminate heel migration
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Men's Shoe Size	Product #	Women's Shoe Size	Product #
Small	6 - 8	DS1016-10	4 - 6	DS1015-10
Medium	8 1/2 - 10	DS1016-12	6 1/2 - 8	DS1015-12
Large	10 1/2 - 12	DS1016-14	8 1/2 - 10	DS1015-14
X-Large	12 1/2 - 14	DS1016-16	—	—

## Darco® Orthowedge™ Shoe

- Reduces weight bearing pressure from the forefoot with 15° dorsi-angle
- Removable 10 mm insole allows customization for offloading specific areas
- Square toe design serves as a bumper for additional protection
- Hook and loop closure for easy application
- Universal right/left

**Indications:** Post-operative foot protection and support; Minor foot fractures

Size	Men's Shoe Size	Women's Shoe Size	Product #
X-Small	—	4 - 7	HS1005-08
Small	6 - 8	7 1/2 - 10	HS1005-10
Medium	8 1/2 - 10	10 1/2 - 13	HS1005-12
Large	10 1/2 - 12	13 1/2 +	HS1005-14
X-Large	12 1/2 - 14	—	HS1005-16



HCPCS Code: L3260  
(Not Covered by Medicare)



## Comfo-Eze Rocker Cast Shoe

- Constructed of durable navy canvas material
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application
- Available in open and closed toe design

**Indications:** Lower extremity cast protection

Size	Length	Width	Open Toe	Closed Toe
Pediatric	5 1/2"	2 1/2"	2004-00	—
X-Small	7"	3 1/2"	2004-01	2038-01
Small	8 1/2"	4 1/2"	2004-02	2038-02
Medium	10"	4 3/4"	2004-03	2038-03
Large	11 1/2"	5 1/4"	2004-04	2038-04
X-Large	13"	5 1/2"	2004-05	2038-05



## **Darco® SlimLine Cast Boot**

- Square toe design serves as a bumper for protection and provides room for larger casts
- Stretchable upper material helps boot conform to cast allowing precise fit
- Rocker sole promotes rolling gait pattern
- Padded insole to cushion patient's foot
- Strapless forefoot closure adjusts easily to cast size

**Indications:** Lower extremity cast protection

Size	Men's Shoe Size	Women's Shoe Size	Product #	Size	Child Shoe Size	Product #
X-Small	—	3 - 5	2039-01	Child/Small	8 - 10 <sup>1/2</sup>	2039-98
Small	—	5 <sup>1/2</sup> - 8	2039-02	Child/Medium	11 - 1 <sup>1/2</sup>	2039-99
Medium	6 - 8	8 <sup>1/2</sup> - 10	2039-03	Child/Large	2 - 4	2039-00
Large	8 <sup>1/2</sup> - 11	10 <sup>1/2</sup> - 13	2039-04			
X-Large	11 <sup>1/2</sup> - 14	13 +	2039-05			

## **DeRoyal® Night Splint**

- Dorsiflexion design applies stretch for effective pain relief
- Lightweight and low profile design for patient comfort
- Adjustable straps to vary amount of stretch
- Removable foam liner for laundering
- Optional foam wedges provided for additional stretch
- Kit features: Night splint, ice pack, exercise band, heel cups, foot control strap and self-management guide

**Indications:** Plantar fasciitis; Achilles tendonitis; Tendonosis

Size	Men's Shoe Size	Women's Shoe Size	Splint	Kit
Small	4 - 6	5 - 7	NS2000-10	SC2000-10
Medium	6 <sup>1/2</sup> - 9	7 <sup>1/2</sup> - 10	NS2000-12	SC2000-12
Large	9 <sup>1/2</sup> +	10 <sup>1/2</sup> +	NS2000-14	SC2000-14



HCPCS Code: L4396



## **Dorsal Night Splint**

- Thermoplastic design holds foot in 90° neutral position
- Neoprene straps with hook and loop closure
- Plastic free bottom and non-skid transfer surface
- Low profile for a cooler, more comfortable fit

**Indications:** Plantar fasciitis; Achilles tendonitis; Tendonosis

Size	Men's Shoe Size	Women's Shoe Size	Product #
Small/Medium	4 - 9	5 - 10	NS2000-20
Large/X-Large	9 <sup>1/2</sup> +	10 <sup>1/2</sup> +	NS2000-21

HCPCS Code: L4396





## ▶ Silicone Heel Cups

- Constructed of medical grade silicone for long-lasting support and comfort
- Soft blue spot helps relieve heel strike fatigue and discomfort through shock absorption
- Packaged 1/pair

**Indications:** Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

Size	Product #
Small	HC2000-10
Medium	HC2000-12
Large	HC2000-14

## ▶ DeRoyal® Heel Suspension Boot

- Smooth or convoluted foam construction
- Posterior panel helps reduce friction against bed surface
- Ventilation holes help promote air circulation and dissipate heat
- Sloped heel suspension wedge pre-attached or free for precise placement
- Extra pad included for additional elevation and rotation protection
- Hook/d-ring strap closures
- Packaged 6/case

**Indications:** Pressure reduction of heel; Pressure protection; Hip rotation; Foot drop; Hip fractures (pre and post operative); Heel decubitus treatment and prevention

Size	Product #
Convoluted	M60-100
Smooth	M60-200
Smooth w/Attached Wedge	M60-300



HCPSC Code: L4396

## ▶ Ankle Contracture Boot

- Static ankle foot orthosis used for positioning and pressure reduction
- Softgoods completely cover splint leaving no exposed metal against skin
- "No strap" closures make it easy for patients to apply and remove boot
- Anti-rotation bar, padded toe piece and mesh laundry bag included
- Universal right/left

**Indications:** Hip rotation; Foot drop; Not recommended for severe plantar flexion contractures, open wounds, severe spasticity or thrombosis

Description	Vel-Foam	Fleece
Ankle Contracture Boot	4302	4307
Ankle Contracture Boot w/Boot Sole	4301	4306
Ankle Contracture Boot w/Boot Sole & Extra Liner	4300	4305

Accessories	Product #
Extra Liner - Vel-Foam	4300K
Extra Liner - Fleece	4305K
Hip Abduction Bar	4300/05AB
Boot Sole	4300/05WK

Sizing Guide	Measure calf circumference 8" up from bottom of heel	Measure foot circumference around widest part of foot
Add letter to end of product # to indicate size	A* = 8" - 12"	A* = 6" - 8"
	B = 10" - 14"	B = 7" - 9"
*Size A available in 4301 style only	C = 14" - 16"	C = 9" - 11"
	D = 16" - 21"	D = 11" - 13"
	E = 16" - 21"	E = 13" +



DeRoyal

 DeRoyal®



A close-up photograph of a person's lower leg and foot wearing a black and white orthopedic brace. The brace has several horizontal straps and a logo that reads "Three-D Orthopedic, Inc.". The person is wearing blue jeans and a white sock. In the background, a golden retriever is sitting on a grassy field under a clear blue sky.

*“Hope*  
is one of the principle springs that  
keeps mankind in motion.”

— Thomas Fuller





A person is shown from the waist down, wearing khaki shorts and a grey sneaker on their left foot. Their right leg is in a black, multi-buckle medical cast. They are holding a fishing rod with a wooden handle and a silver reel. The background is a blurred body of water and trees. The image has a decorative wavy border at the bottom in shades of red and orange.

**“Effort only fully  
releases its rewards after a  
person *refuses to quit.*”**

— Napoleon Hill

