

## Frequently Asked Questions / Answers

- 1. Can the DeClogger be reused?**  
The DeClogger is designed for single-use only.
- 2. Can our facility get an in-house demonstration?**  
The DeClogger is simple and easy to use; just follow the step-by-step protocol for use of the DeClogger.
- 3. Can the DeClogger harm the patient in any way?**  
The DeClogger has a small disc at the handle that prevents over-insertion. The design of the device with its standardized length and flexible plastic molding are added safeguards.
- 4. Can the DeClogger work on a nasogastric or nasojejunal tube?**  
No, these tubes are longer in length and have a smaller French size.