Hammock Sling



Multi-Purpose Seated Transfer Sling

- Provides support around the hips and sacral area and under the thighs up to the shoulders with optional head support
- · Multiple loop and strap options can accommodate various seated and semi-reclined positions
- Divided leg supports that can be applied and removed in a seated or supine position
- Leg straps can be applied in various positions for different levels of support
- · Built-in caregiver handles on each side add an additional level of control when moving and transferring user
- Small commode opening

Materials









Applications



Up from the floor



Transfer to tub



Higher up in bed

Positioning in chair



From bed to wheelchair



From bed to commode



From wheelchair to toilet

Hammock Sling

Part Numbers

Mesh / Poly with Head Support

Part No.	Size	SWL
8D2730	Child	450 lb (204 kg)
8D2630	Junior	450 lb (204 kg)
8D2530	S	625 lb (283 kg)
8D2430	M	625 lb (283 kg)
8D2330	L	625 lb (283 kg)
8D2230	XL	1,000 lb (454 kg)

PolySlip with Head Support

Size	SWL
Child	450 lb (204 kg)
Junior	450 lb (204 kg)
S	625 lb (283 kg)
M	625 lb (283 kg)
L	625 lb (283 kg)
XL	1,000 lb (454 kg)
	Child Junior S M L

Spacer with Head Support

Part No.	Size	SWL
8D4730	Child	450 lb (204 kg)
8D4630	Junior	450 lb (204 kg)
8D4530	S	450 lb (204 kg)
8D4430	M	450 lb (204 kg)
8D4330	L	450 lb (204 kg)
8D4230	XL	450 lb (204 kg)

Quilted with Head Support

		• •
Part No.	Size	SWL
537105	Junior	800 lb (363 kg)
537110	S	800 lb (363 kg)
537120	М	800 lb (363 kg)
537130	L	800 lb (363 kg)
537140	XL	800 lb (363 kg)
537150	XXL	800 lb (363 kg)