



Re-Order #	Description	Size	Unit
401165	Headboard/Wall Protector	3" x 6" x 38"	Each

Purpose

The headboard wall protector is easily attached to the headboard with Velcro[®] straps and protects the wall from costly repairs resulting from the bed bumping the wall.

Easy to Use


Place the headboard protector between the bed and wall and attach with the Velcro[®] straps. The notched out, non-slip surface should face the bed while aligning the notches with the bed hardware.

Easy Care

The nylon cover can be wiped clean with a spray and wipe cleanser. If disinfected, avoid using chlorine (bleach) chemicals.


VINYL HEADBOARD PAD

Re-Order #	Size	Unit
401160	13" x 2" x 36"	Each



VINYL FOOTBOARD PAD


Re-Order #	Size	Unit
401170	9" x 2" x 36"	Each



Added protection for your restless or agitated patients. Amply padded with inch-thick polyfoam, these HealthGrade vinyl pads are non-allergenic and conform to the flammability standards set forth in California Bulletin 117B.

EASY TO USE AND
CARE FOR.

Both pads are held in place by elastic bands that go around the headboard or footboard. (See illustration.) These vinyl pads clean easily with a damp cloth and mild detergent.



Skil-Care Corporation . . . When Safety Matters and Comfort Counts