



## VIVE HEAVY DUTY WHEELCHAIR

### MOB1078

Experience a durable yet lightweight mobility solution designed for a variety of individual needs. With a strong steel frame, this wheelchair can support up to 400 pounds. It can also be folded for easy transportation and storage. The 20-inch wide seat and padded armrests ensure comfort, while the adjustable footrests allow for a more personalized fit.


#### WHAT'S INCLUDED

- 1 Heavy Duty Wheelchair

#### COLORS:



## HEAVY-DUTY AND LIGHTWEIGHT





**Product Dimensions:**

Seat width:	20"
Seat height:	19.50"
Seat depth:	16"
Back height:	18"
Total width:	27"
Total height:	36.3"
Total length(with footrest):	42"
Max loading:	400LBS



**HEAVY-DUTY AND LIGHTWEIGHT**

The Vive Foldable Wheelchair Transport offers robust support, accommodating up to 400 pounds. With a lightweight steel frame, this wheelchair is portable and makes for an excellent travel companion. It also folds, for hassle-free storage when not in use.

**ADJUSTABLE FOR GREATER COMFORT**

Experience superior comfort with this wheelchair's 20-inch wide seat, providing extra space for a more relaxed sitting position. You can customize your comfort further with the adjustable footrests, which conveniently swing away if needed. The padded armrests ensure comfort, while the easy-to-clean material of the chair ensures a cozy and hygienic experience.

**ENHANCED MANEUVERABILITY**

Navigate easily thanks to the effortless push design and the durable wheels on the Vive Foldable Wheelchair Transport. The hand-breaking system enhances control and safety, for greater peace of mind when maneuvering the chair. This wheelchair fits through most doorways with a total width of 27 inches, so you can move between indoor and outdoor spaces.

**VERSATILE DESIGN FOR A VARIETY OF NEEDS**

With four sturdy wheels and an adjustable footrest, the Vive Foldable Wheelchair Transport is great for individuals with a variety of needs. Provide seamless transportation between rooms and from hospitals to your home, or use this wheelchair to add stability during travel. With footrests that swing open, this wheelchair is easy to access and exit, making it great for those who are recovering from surgery or injury.

