

The new NYOrtho ULTRA Hip Protector combines our proven seamless fixation technology with energy-absorbing ultra-foam technology, to provide superior and comfortable protection for the trochanter.

### Ultra Hip Protector

RECOMMENDED USE: Patients at risk for falling; those at risk for a hip fracture

- High impact-absorbing viscoelastic pads
- 100% latex-free inner wear briefs are **ANTI-MICROBIAL**, soft and comfortably worn under clothing
- Low-profile pads are virtually undetectable
- Standard sizes S-XXL
- Unisex

#### STANDARD SIZES

Model	Size	Waist	Hips
9962-S	Small	26" - 30"	33" - 37"
9962-M	Medium	30" - 34"	37" - 41"
9962-L	Large	34" - 38"	41" - 45"
9962-XL	XL	38" - 42"	45" - 49"
9962-2XL	2XL	42" - 46"	49" - 53"

To insure proper fit, compare the waist measurement to hip measurement and select the size that is larger.



### Hip Fractures

More than 95% of all hip fractures are caused by falls. The majority of hip fractures occur in adults over 65 years old. Fall-related injuries such as hip fractures are associated with increased mortality and a major cause of loss of functional capacity and impaired quality of life. Hip fractures are one of the largest health concerns in nursing homes as well as an extremely serious threat to the quality of life and the future health of elderly patients. Our Ultra Hip Protector offers distinctive advantages over competing products, through the use of unique materials.

### Laundering Instructions

The ULTRA Hip Protector may be washed up to 120° F and dried at medium temperature. Check foam pads after each wash for evidence of cracked or broken pads. If the pad is damaged, discontinue use.

#### WARNING:

If a patient falls, replace the brief before patient use. NYOrtho does not guarantee that hip trauma will not occur.

**ANTI-MICROBIAL**

# NYOrtho

## ULTRA Hip Protector