

How to Use DeCloggers

There are boatloads of products that you look at and marvel that no one has invented them before.

Bionix Tube DeCloggers are one of those products.

Using the DeCloggers and providing relief and care to your patients is so simple that we had to think about how to make the *How to Use* them section longer. But we couldn't. Here's what you do:

1. Select the right DeClogger (they're color coded so even that should be simple).
2. Insert the DeClogger to reach the blockage and slowly rotate clockwise.
3. Take the DeClogger out.
4. Discard the DeClogger after use (intended for single patient use only).

