


# Multiglide/Spilerglide

**immedia**<sup>®</sup>  
a part of Etac

## Short Instruction



### Material

100% nylon si/si  
100% nylon pu/si

Label on the outside



Class 1 - Medical device directive

MDD 93/42/EEC

Shrink 5%.

This product can be cleaned with a solvent-free cleaning agent with a pH level between 5-9, or a 70% disinfection solution. Do not add fabric conditioner.

### Item no

### EAN no

### Size in mm

IM00100	5708012375127	W1000 x L700
IM00108	5708012001347	W1080 x L800
IM00120/1	5708012000258	W1200 x L900
IM00130	5708012000241	W1300 x L800
IM00139	5708012003174	W1300 x L900
IM00150/70	5708012511761	W1430 x L700
IM00170	5708012003181	W1400 x L1200
IM00200	5708012001743	W2000 x L700
IM0040	5708012000449	W400 x L400
IM0066	5708012499199	W600 x L600
IM0077	5708012375080	W720 x L700
IM0080	5708012375097	W720 x L1000
IM0082	5708012001729	W700 x L2000
IM0083	5708012375110	W720 x L1300
IM100	5708012140619	W1000 x L700
IM100/140	5708012226061	W1000 x L1400
IM102	5708012002252	W1000 x L2000
IM103/1	5708012000326	W1000 x L1350
IM108	5708012000319	W1080 x L800
IM120	5708012001668	W1200 x L700
IM130	5708012000302	W1300 x L800
IM168	5708012000296	W1600 x L800
IM200	5708012000289	W2000 x L700
IM40	5708012120871	W400 x L400
IM62	5708012002254	W600 x L250
IM66	5708012499205	W600 x L600
IM77	5708012140602	W720 x L700
IM80	5708012128761	W720 x L1000
IM82	5708012254477	W700 x L2000
IM83	5708012169764	W720 x L1300
IM88	5708012226078	W800 x L800
IM2245	5708012389544	W1500 x L10000
IM2249	5708012001637	W1500 x L10000

# Multiglide/Spilerglide



## Check always:



## Contact / Download

Etac A/S  
info@etac.dk  
www.etac.dk

Etac GmbH  
info@etac.de  
www.etac.de

Etac A/S  
hovedkontor.norge@etac.com  
www.etac.com

R82 UK Ltd.  
R82uk@r82.com  
www.r82.co.uk

Etac AB  
info@etac.se  
www.etac.se

Etac Export  
etac.immedia@etac.com  
www.etac.com/exp

Etac R82 BV  
info.holland@etac.com  
www.etac.com

Snug Seat Inc.  
information@snugseat.com  
www.snugseat.com

## Manufacture

Etac Immedia A/S  
Egeskovvej 12  
DK-8700 Horsens  
etac.immedia@etac.com

**etac**<sup>®</sup>  
Creating possibilities