

PediTurn

Short Instruction

immedia[®]
a part of Etac



Material

Toughtek, 100% nylon, PU
ABS/TPU, HDPE, brass, steel

Class 1 - Medical device directive
MDD 93/42/EEC


Cleaning

Do not add fabric conditioner.
This product can be cleaned with a
solvent-free cleaning agent with a pH
level between 5-9, or a 70% disinfection
solution.



Item no	EAN no	Size in mm
IM99836GT	5708012169702	Ø 400
IM99700	5708012001507	Ø 380

PediTurn



Check always:



Contact / Download

Etac A/S
info@etac.dk
www.etac.dk

Etac GmbH
info@etac.de
www.etac.de

Etac A/S
hovedkontor.norge@etac.com
www.etac.com

R82 UK Ltd.
enquiries@etac.uk.r82.com
www.etac.com/uk

Etac AB
info@etac.se
www.etac.se

Etac Export
etac.immedia@etac.com
www.etac.com/exp

Etac R82 BV
info.holland@etac.com
www.etac.com

R82 Inc.
information@r82.com
www.r82.com

Manufacture

Etac Immedia A/S
Egeskovvej 12
DK-8700 Horsens
etac.immedia@etac.com

 **etac**[®]
Creating possibilities