

JOBST® Opaque

NEW!
Midnight Navy



**Softness and ease
of donning for
exceptionally
comfortable leg
therapy.**

**Medical Compression Stockings
15-20mmHg*, 20-30mmHg*, 30-40mmHg***

#1 physician recommended



Comfort, Health and Style!



Experience comfort.

Features

Benefits



Enhanced with a silicone blend softener

Unmatched softness and superior moisture wicking for cool and comfortable legs.

Contoured panty and anatomical foot design with a relaxed toe cap and enlarged heel pocket

Opaque is easier to don than any other stocking you have worn before.

Refined opacity

Ensures a perfect fit from hip to toe all day long for exceptional wearing comfort.

Selected, strong yarns and a reinforced heel and toe


Covers blemishes and scars giving your leg an elegant appearance.

Fashionable but durable backed by Our Total Satisfaction Guarantee.

What is JOBST® Advanced Comfort

JOBST® Advanced Comfort combines exceptional yarns with a unique finish to create unbelievable wearing comfort. This technology encompasses:

- Unmatched Softness: finest yarns and unique finishes are combined to give a perfect soft touch
- Superior Moisture Wicking: Jobst Advanced Comfort boosts the spreading speed of any moisture on the fabric so that it evaporates faster. This keeps your legs comfortable all day long.



#1 physician recommended



Comfort, Health and Style!



Want softer compression therapy? Here you go ...


Why JOBST® Opaque?

- JOBST® Opaque is the softest Opaque ever and JOBST® Advanced Comfort Opaque gives unbelievable wearing comfort.
- The contoured panty provides a perfect fit around the hips to make our waist highs even more comfortable.
- JOBST®'s silicone blend softener makes Opaque very easy to don.
- New elegant collection of colors to mix and match with all wardrobes of any season.

Women voted. We delivered.

Women have voted! Trials have proven that women clearly prefer the new Opaque. 90% of women prefer the new Opaque vs. the current Opaque due to its perfect soft touch*.

* Blinded preference trial carried out by JOBST®. Data available upon request.



A variety of colors and styles to choose from.

Colors



Classic Black



Midnight Navy



Espresso



Honey



Natural

Actual garment color may vary from image shown.

Styles



**Knee High
Closed Toe**



**Knee High
Open Toe**



**Thigh High
Closed Toe**



**Thigh High
Open Toe**



Pantyhose

Knee High ordering information:

**Knee High
(Closed Toe)**



| Size | Catalog Number/Color | | | | | Petite (Closed Toe) | Size | Catalog Number/Color | | Compression mmHg* |
|------|----------------------|--------|---------|----------|---------------|------------------------|--------|----------------------|------------|-------------------|
| | Classic Black | Honey | Natural | Espresso | Midnight Navy | | | Classic Black | Natural | |
| S | 115200 | 115681 | 115212 | 115685 | 115752 | S | 115600 | 115604 | 15-20 mmHg | |
| M | 115201 | 115682 | 115213 | 115686 | 115753 | M | 115601 | 115605 | 15-20 mmHg | |
| L | 115202 | 115683 | 115214 | 115687 | 115754 | L | 115602 | 115606 | 15-20 mmHg | |
| XL | 115203 | 115684 | 115215 | 115688 | 115755 | XL | 115603 | 115607 | 15-20 mmHg | |
| S | 115132 | 115456 | 115270 | 115460 | 115740 | S | 115608 | 115614 | 20-30 mmHg | |
| M | 115133 | 115457 | 115271 | 115461 | 115741 | M | 115609 | 115615 | 20-30 mmHg | |
| L | 115134 | 115458 | 115272 | 115462 | 115742 | L | 115610 | 115616 | 20-30 mmHg | |
| XL | 115135 | 115459 | 115273 | 115463 | 115743 | XL | 115611 | 115617 | 20-30 mmHg | |
| LFC | 115364 | | 115378 | | | | | | 20-30 mmHg | |
| XLFC | 115365 | | 115379 | | | | | | 20-30 mmHg | |
| S | 115168 | 115701 | 115282 | 115705 | | S | 115626 | 115632 | 30-40 mmHg | |
| M | 115169 | 115702 | 115283 | 115706 | | M | 115627 | 115633 | 30-40 mmHg | |
| L | 115170 | 115703 | 115284 | 115707 | | L | 115628 | 115634 | 30-40 mmHg | |
| XL | 115171 | 115704 | 115285 | 115708 | | XL | 115629 | 115635 | 30-40 mmHg | |
| LFC | 115368 | | 115370 | | | | | | 30-40 mmHg | |
| XLFC | 115369 | | 115371 | | | | | | 30-40 mmHg | |

#1 physician recommended



Comfort, Health and Style!

This panel size is 7 3/4" wide to allow for hole punch at the fold in

Knee High ordering information:

Knee High
(Open Toe)



| Catalog Number/Color | | | Catalog Number/Color Compression mmHg* | | |
|----------------------|---------------|---------|--|---------|-------------------|
| Size | Classic Black | Natural | Petite (Open Toe) Size | Natural | Compression mmHg* |
| S | 115335 | 115331 | S | 115069 | 15-20 mmHg |
| M | 115336 | 115332 | M | 115070 | 15-20 mmHg |
| L | 115337 | 115333 | L | 115071 | 15-20 mmHg |
| XL | 115338 | 115334 | XL | 115072 | 15-20 mmHg |
| S | 115384 | 115480 | S | 115620 | 20-30 mmHg |
| M | 115385 | 115481 | M | 115621 | 20-30 mmHg |
| L | 115386 | 115482 | L | 115622 | 20-30 mmHg |
| XL | 115387 | 115483 | XL | 115623 | 20-30 mmHg |
| LFC | 115376 | 115366 | | | 20-30 mmHg |
| XLFC | 115377 | 115367 | | | 20-30 mmHg |
| S | 115388 | 115496 | S | 115638 | 30-40 mmHg |
| M | 115389 | 115497 | M | 115639 | 30-40 mmHg |
| L | 115390 | 115498 | L | 115640 | 30-40 mmHg |
| XL | 115391 | 115499 | XL | 115641 | 30-40 mmHg |
| LFC | 115380 | 115382 | | | 30-40 mmHg |
| XLFC | 115381 | 115383 | | | 30-40 mmHg |

Fiber contents:

15-20 mmHg*: 78% Nylon and 22% Spandex/Elasthan; 20-30 mmHg*: 78% Nylon and 22% Spandex/Elasthan; 30-40 mmHg*: 70% Nylon and 30% Spandex/Elasthan

Thigh High ordering information:

Thigh High
(Closed Toe)

(Silicone Dot Band)



| Catalog Number/Color | | | | | Catalog Number/Color Compression mmHg* | | | |
|----------------------|---------------|--------|---------|----------|--|--|---------|-------------------|
| Size | Classic Black | Honey | Natural | Espresso | Midnight Navy | Petite (Closed Toe) (Silicone Dot Band) Size | Natural | Compression mmHg* |
| S | 115504 | 115689 | 115508 | 115693 | 115760 | S | 115652 | 15-20 mmHg |
| M | 115505 | 115690 | 115509 | 115694 | 115761 | M | 115653 | 15-20 mmHg |
| L | 115506 | 115691 | 115510 | 115695 | 115762 | L | 115654 | 15-20 mmHg |
| XL | 115507 | 115692 | 115511 | 115696 | 115763 | XL | 115655 | 15-20 mmHg |
| S | 115144 | 115464 | 115274 | 115468 | 115748 | S | 115660 | 20-30 mmHg |
| M | 115145 | 115465 | 115275 | 115469 | 115749 | M | 115661 | 20-30 mmHg |
| L | 115146 | 115466 | 115276 | 115470 | 115750 | L | 115662 | 20-30 mmHg |
| XL | 115147 | 115467 | 115277 | 115471 | 115751 | XL | 115663 | 20-30 mmHg |
| S | 115180 | 115709 | 115286 | 115713 | | S | 115672 | 30-40 mmHg |
| M | 115181 | 115710 | 115287 | 115714 | | M | 115673 | 30-40 mmHg |
| L | 115182 | 115711 | 115288 | 115715 | | L | 115674 | 30-40 mmHg |
| XL | 115183 | 115712 | 115289 | 115716 | | XL | 115675 | 30-40 mmHg |

Thigh High ordering information:

Thigh High
(Open Toe)

(Silicone Dot Band)



| Catalog Number/Color | | | Catalog Number/Color Compression mmHg* | | | |
|----------------------|---------------|--------|--|--|---------|-------------------|
| Size | Classic Black | Honey | Natural | Petite (Open Toe) (Silicone Dot Band) Size | Natural | Compression mmHg* |
| S | 115556 | 115697 | 115552 | S | 115644 | 15-20 mmHg |
| M | 115557 | 115698 | 115553 | M | 115645 | 15-20 mmHg |
| L | 115558 | 115699 | 115554 | L | 115646 | 15-20 mmHg |
| XL | 115559 | 115700 | 115555 | XL | 115647 | 15-20 mmHg |
| S | 115560 | 115472 | 115544 | S | 115664 | 20-30 mmHg |
| M | 115561 | 115473 | 115545 | M | 115665 | 20-30 mmHg |
| L | 115562 | 115474 | 115546 | L | 115666 | 20-30 mmHg |
| XL | 115563 | 115475 | 115547 | XL | 115667 | 20-30 mmHg |
| S | 115564 | 115717 | 115548 | S | 115676 | 30-40 mmHg |
| M | 115565 | 115718 | 115549 | M | 115677 | 30-40 mmHg |
| L | 115566 | 115719 | 115550 | L | 115678 | 30-40 mmHg |
| XL | 115567 | 115720 | 115551 | XL | 115679 | 30-40 mmHg |


Fiber contents:

15-20 mmHg*: 72% Nylon, 18% Spandex, 10% Silicone; 20-30 mmHg*: 70% Nylon, 16% Spandex, 14% Silicone; 30-40 mmHg*: 62% Nylon, 28% Spandex, 10% Silicone

JOBST® Opaque

15-20 mmHg*, 20-30 mmHg*, 30-40 mmHg*







Pantyhose ordering information:

| Pantyhose | Catalog Number/Color | | | Compression mmHg* | |
|--|----------------------|---------------|---------|-------------------|-----------------------|
| | Size | Classic Black | Natural | | NEW! Midnight Navy |
|  | S | 115220 | 115232 | 115756 | 15-20 mmHg |
| | M | 115221 | 115233 | 115757 | 15-20 mmHg |
| | L | 115222 | 115234 | 115758 | 15-20 mmHg |
| | XL | 115223 | 115235 | 115759 | 15-20 mmHg |
| | S | 115156 | 115278 | 115744 | 20-30 mmHg |
| | M | 115157 | 115279 | 115745 | 20-30 mmHg |
| | L | 115158 | 115280 | 115746 | 20-30 mmHg |
| | XL | 115159 | 115281 | 115747 | 20-30 mmHg |
| | S | 115192 | 115290 | | 30-40mmHg |
| | M | 115193 | 115291 | | 30-40mmHg |
| | L | 115194 | 115292 | | 30-40mmHg |
| | XL | 115195 | 115293 | | 30-40mmHg |

Fiber contents:

15-20 mmHg*: 76% Nylon, 24% Spandex; 20-30 mmHg*: 76% Nylon, 24% Spandex; 30-40 mmHg*: 66% Nylon, 34% Spandex

Sizing chart:

| SIZE | ANKLE  | CALF  | THIGH  | HIP  |
|---|---|--|---|--|
| S | 7"-8 ¹ / ₄ " (18-21 cm) | 11"-15" (28-38 cm) | 15 ³ / ₄ "-24 ³ / ₈ " (40-62 cm) | 28"-46" (71-117 cm) |
| M | 8 ³ / ₈ "-9 ⁷ / ₈ " (21-25 cm) | 11 ⁷ / ₈ "-16 ¹ / ₂ " (30-42 cm) | 18 ¹ / ₈ "-27 ¹ / ₂ " (46-70 cm) | 30"-50" (76-127 cm) |
| L | 10"-11 ³ / ₈ " (25-29 cm) | 12 ¹ / ₂ "-18 ¹ / ₈ " (32-46 cm) | 21 ¹ / ₄ "-30 ³ / ₄ " (54-78 cm) | 32"-54" (81-137 cm) |
| XL | 11 ¹ / ₂ "-13" (29-33 cm) | 13 ³ / ₈ "-19 ⁵ / ₈ " (34-50 cm) | 23 ⁵ / ₈ "-32" (60-81 cm) | 40"-65" (102-166 cm) |
| LFC | 10"-11 ⁷ / ₈ " (25-30 cm) | 18"-24" (46-61 cm) | | |
| XLFC | 12"-14" (30-36 cm) | 18"-24" (46-61 cm) | | |
|  | Knee High Regular: Greater than 15" (38 cm) Knee High Petite: 15" (38 cm) or less | |  | Thigh High Regular: Greater than 27" (68.5cm) Thigh High Petite: 27" (68.5cm) or less |

For Knee High, Measure Ankle & Calf. For Thigh High, Measure Ankle, Calf & Thigh. For Waist High, Measure Ankle, Thigh & Hip.

Care instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.



15-20 mmHg*
Compression



20-30 mmHg*
Compression



30-40 mmHg*
Compression

* The mean compression of an average ankle size.
Please refer to the product label and /or package insert for full instructions on the safe use of these products.

