JOBST® ActiveWear

Features

Benefits

| dri-release® MICROBLEND PERFORMANCE PARRICE By Optimize My Optimize | Feels like cotton and offers superior moisture management to help keep feet dry, cool and comfortable | |
|--|---|--|
| Extra cushioned foot with a reinforced heel | Enhances wearing comfort and resists abrasion | |
| Seamless toe | Provides all-day wearing comfort | |
| All-day comfort knee band | Keeps socks up without binding or pinching | |
| Designed for both men and women | Accommodates a wider group of gradient compression stocking users | |

Unisex style

to suit both men and women



Colors

Cool White

for an active look

Cool Black

Actual garment color may vary from printed swatch sample.

Care instructions

Machine washable in lukewarm water.
Tumble dry in low heat.
Do not bleach, dry clean, or iron.

Sizing chart:

| SIZE | ANKLE — | CALF |
|------|---|---------------------------------------|
| S | 7"-8 ¹ / ₄ " (18-21 cm) | 11"-15" (28-38 cm) |
| М | 8³/s"-9 ⁷ /s" (21-25 cm) | 11 ⁷ /8"-16¹/2" (30-42 cm) |
| L | 10"-113/8" (25-29 cm) | 12¹/2"-18¹/8" (32-46 cm) |
| XL | 11½"-13" (29-33 cm) | 13³/s"-195/s" (34-50 cm) |
| LFC | 10"-11 ⁷ /8" (25-30 cm) | 18"-24" (46-61 cm) |
| XLFC | 12"-14" (30-36 cm) | 18"-24" (46-61 cm) |

Knee High ordering information:

Knee High (Closed Toe)



| Size | Cool White | Cool Black | Compression mmHg* |
|------|-------------------|-------------------|-------------------|
| | Catalog No. | Catalog No. | |
| S | 110479 | 110483 | 15-20 mmHg |
| M | 110480 | 110484 | 15-20 mmHg |
| L | 110481 | 110485 | 15-20 mmHg |
| XL | 110482 | 110486 | 15-20 mmHg |
| LFC | 110525 | 110531 | 15-20 mmHg |
| XLFC | 110526 | 110532 | 15-20 mmHg |
| S | 110489 | 110493 | 20-30 mmHg |
| M | 110490 | 110494 | 20-30 mmHg |
| L | 110491 | 110495 | 20-30 mmHg |
| XL | 110492 | 110496 | 20-30 mmHg |
| LFC | 110527 | 110533 | 20-30 mmHg |
| XLFC | 110528 | 110534 | 20-30 mmHg |
| S | 110051 | 110055 | 30-40 mmHg |
| M | 110052 | 110056 | 30-40 mmHg |
| L | 110053 | 110057 | 30-40 mmHg |
| XL | 110054 | 110058 | 30-40 mmHg |
| LFC | 110529 | 110535 | 30-40 mmHg |
| XLFC | 110530 | 110536 | 30-40 mmHg |

Fiber contents:

15-20 mmHg* and 20-30 mmHg*: 60% Polyester, 20% Nylon, 11% Cotton and 9% Spandex/Elasthan 30-40 mmHg*: 58% Polyester, 20% Nylon, 10% Cotton and 12% Spandex/Elasthan

Please refer to the product label and/or package insert for full instructions on the safe use of these products.









New Sizes Full calf Sizes JOBST ActiveWear



Effective leg therapy in an energizing athletic sock

Medical Compression Socks 15-20mmHg*, 20-30mmHg*, and 30-40mmHg*



© 2010 BSN medical Inc.

^{*} The mean compression for an average ankle size.

^{**} Tested against Sigvaris® Cushioned Cotton, Venosan® Sportsline and Mediven® Motion Sport.

Effective leg therapy in an energizing athletic sock

JOBST[®] ActiveWear


360° cushioned foot


JOBST® ActiveWear offers gradient compression socks that fit active lifestyles. They look and feel just like performance athletic socks! The ideal combination of therapeutic efficacy and Dri-release® yarn for superior moisture management.

Precise gradient compression in three compression classes provides effective relief and management of vascular disorders.

Superior moisture management

- Tests prove that JOBST® ActiveWear is superior in temperature control and moisture management.
- Dri-release®'s unique blend of natural and synthetic fibers helps keep feet dry, cool and comfortable.







More durable

- Tests prove the unique design of state-of-the-art polyester, nylon, cotton and spandex provides comfort and more durability than leading competitors.
- Extra cushioned foot with a reinforced heel and toe provide greater durability.

Easier-to-don

- Tests prove that JOBST® ActiveWear is easier to don at heel and foot than leading competitors.
- Easier donning ensures maximum comfort and compliance.





Knee high – satisfying in Unisex Design All-day comfort knee band every detail... keeps socks up without binding or pinching dri-release for superior moisture management IOBST® ActiveWear Seamless toe for increased wearing Reinforced heel comfort resists abrasion for greater durability

JOBST ActiveWear is superior to its competition:

| Attribute | Sigvaris® Cushioned Cotton | Venosan® Sportsline | Mediven [™] Motion <i>Sport</i> | JOBST® ActiveWear |
|---|-------------------------------|------------------------|---|--------------------------------|
| Toe Seam (picture of the inside of the foot, at the same angle and distance) | | | | seamless Toe! |
| Padding (picture of the inside of the foot, at the same angle and distance) | | | | 360° padding! |
| Knee Band (picture of the inside of the foot, at the same angle and distance) | | | | wide comfortable knee band! |
| Toe Temperature °F (after 30 min. workout) | 84° | 81° | 81° | 77° |
| Moisture Wicking (time to absorb water droplet in the inside of the sock) | 3 seconds | 66 seconds | 55 seconds | 1 second |
| Durability (abrasion: cycles to hole in the heel) | 800 cycles | 1,000 cycles | 188 cycles | 1,125 cycles |
| 3D-Donning (ankle/heel/foot force in kg) | 3.1/1.3/2.9 | 1.7/1.4/1.9 | 4.0/0.8/2.0 | 2.55/0.5/1.8 |

Dri-release® is an Optimer Inc. trademark for its microblend performance fabric. US PATENT #5,888,914.