










Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.


| Lift Types | | Spreader Bar | Sling SWL |
|--------------------------------------|-------------------------------|--|-----------|
| Advance Elara HPL450 HPL500 | HPL700 Presence Stature | Loop  | 500 lbs. |




| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | | | | | | | | | | |
|---|---|--------|--------------------------------|---------|---------|-----|-----|-----|-----|---------|-----|-----|--|--|--|--|--|--|--|--|--|--|--|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 | | | | | | | | | | | |
| Access  Material: Padded | XS | Brown | NA1600 | | | | | | | | | | | | | | | | | | | | |
| | S | Red | NA1601 | | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA1603 | | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA1606 | | | | | | | | | | | | | | | | | | |
| | XL | Blue | | | | | | | | NA1608 | | | | | | | | | | | | | |
| Access, Head Support  Material: Padded | S | Red | NA1626 | | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA1627 | | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA1628 | | | | | | | | | | | | | | | | | | |
| | XL | Blue | | | | | | | | NA1629 | | | | | | | | | | | | | |
| | Quickfit  Material: Padded | XS | Brown | NA1002P | | | | | | | | | | | | | | | | | | | |
| S | | Red | NA1000P | | | | | | | | | | | | | | | | | | | | |
| M | | Yellow | | NA1003P | | | | | | | | | | | | | | | | | | | |
| L | | Green | | | NA1006P | | | | | | | | | | | | | | | | | | |
| XL | | Blue | | | | | | | | NA1008P | | | | | | | | | | | | | |
| Quickfit  Material: Mesh | S | Red | NA1001 | | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA1004 | | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA1007 | | | | | | | | | | | | | | | | | | |

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | | | | | | | | | |
|---|------|--------|--------------------------------|---------|---------|-----|-----|-----|-----|---------|-----|-----|--|--|--|--|--|--|--|--|--|--|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 | | | | | | | | | | |
| Full Back  Material: Padded | XS | Brown | NA1067 | | | | | | | | | | | | | | | | | | | |
| | S | Red | NA1068 | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA1069 | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA1070 | | | | | | | | | | | | | | | | | |
| | XL | Blue | | | | | | | | NA1071 | | | | | | | | | | | | |
| Full Back  Material: Mesh | XS | Brown | NA1080 | | | | | | | | | | | | | | | | | | | |
| | S | Red | NA1081 | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA1082 | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA1083 | | | | | | | | | | | | | | | | | |
| | XL | Blue | | | | | | | | NA1084 | | | | | | | | | | | | |
| Full Back Disposable  Material: Disposable | XS | Brown | NA25620 | | | | | | | | | | | | | | | | | | | |
| | S | Red | NA25621 | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA25622 | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA25623 | | | | | | | | | | | | | | | | | |
| | XL | Blue | | | | | | | | NA25624 | | | | | | | | | | | | |
| Full Body Hammock  Material: Polyester | S | Red | NA3171 | | | | | | | | | | | | | | | | | | | |
| | M | Yellow | | NA3172 | | | | | | | | | | | | | | | | | | |
| | L | Green | | | NA3173 | | | | | | | | | | | | | | | | | |





Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | | Spreader Bar | Sling SWL |
|--------------------------------------|-------------------------------|--|-----------|
| Advance Elara HPL450 HPL500 | HPL700 Presence Stature | Loop  | 500 lbs. |

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|------|--------|--------------------------------|----|---------|---------|-----|-----|-----|--------|-----|-----|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| Quickfit Deluxe  Material: Polyester | XS | Brown | NA1048 | | | | | | | | | |
| | S | Red | NA1050 | | | | | | | | | |
| | M | Yellow | | | | NA1053 | | | | | | |
| | L | Green | | | NA1056 | | | | | | | |
| | XL | Blue | | | | | | | | NA1058 | | |
| Quickfit Deluxe  Material: Mesh | S | Red | NA1051 | | | | | | | | | |
| | M | Yellow | | | | NA1054 | | | | | | |
| | L | Green | | | NA1057 | | | | | | | |
| | XL | Blue | | | | | | | | NA1059 | | |
| Quickfit Deluxe  Material: Padded | S | Red | NA1050P | | | | | | | | | |
| | M | Yellow | | | | NA1053P | | | | | | |
| | L | Green | | | NA1056P | | | | | | | |


*All slings on this page may be used for amputees

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|------|--------|--------------------------------|----|--------|--------|-----|-----|-----|--------|-----|-----|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| Long Seat  Material: Polyester | S | Red | NA1100 | | | | | | | | | |
| | M | Yellow | | | | NA1103 | | | | | | |
| | L | Green | | | NA1106 | | | | | | | |
| | XL | Blue | | | | | | | | NA1109 | | |
| Long Seat  Material: Mesh | S | Red | NA1101 | | | | | | | | | |
| | M | Yellow | | | | NA1104 | | | | | | |
| | L | Green | | | NA1107 | | | | | | | |
| Long Seat Commode  Material: Polyester | S | Red | NA1125 | | | | | | | | | |
| | M | Yellow | | | | NA1128 | | | | | | |
| | L | Green | | | NA1131 | | | | | | | |
| Long Seat Commode  Material: Mesh | S | Red | NA1126 | | | | | | | | | |
| | M | Yellow | | | | NA1129 | | | | | | |
| | L | Green | | | NA1132 | | | | | | | |





NA2650 Extension Straps available add length to the existing sling straps

Hoyer® Clip Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | Cradle | Sling SWL |
|------------------|--|-----------|
| Presence Stature | Clip  | 500 lbs. |


| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|------|--------|--------------------------------|----|---------|---------|-----|-----|-----|---------|---------|-----|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| Comfort Standard  Material: Polyester | XS | Brown | NA25501 | | | | | | | | | |
| | S | Red | NA25502 | | | | | | | | | |
| | M | Yellow | | | | NA25503 | | | | | | |
| | L | Green | | | NA25504 | | | | | | | |
| | XL | Blue | | | | | | | | | NA25505 | |
| Comfort  Material: Mesh | S | Red | NA25512 | | | | | | | | | |
| | M | Yellow | | | | NA25513 | | | | | | |
| | L | Green | | | NA25514 | | | | | | | |
| | XL | Blue | | | | | | | | | NA25515 | |
| Comfort  Material: Padded | XS | Brown | NA25506 | | | | | | | | | |
| | S | Red | NA25507 | | | | | | | | | |
| | M | Yellow | | | | NA25508 | | | | | | |
| | L | Green | | | NA25509 | | | | | | | |
| Comfort Disposable  Material: Disposable | XS | Brown | NA25610 | | | | | | | | | |
| | S | Red | NA25611 | | | | | | | | | |
| | M | Yellow | | | | NA25612 | | | | | | |
| | L | Green | | | NA25613 | | | | | | | |
| XL | Blue | | | | | | | | | NA25614 | | |



| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|----------|----------|--------------------------------|------|-------------|---------|---------|-----|-----|-----|---------|---------|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| Comfort Access  Material: Padded | S | Red | NA25571 | | | | | | | | | |
| | M | Yellow | | | | NA25572 | | | | | | |
| | L | Green | | | NA25573 | | | | | | | |
| | XL | Blue | | | | | | | | | NA25574 | |
| Comfort Access Head Support  Material: Padded | M | Yellow | | | | | NA25577 | | | | | |
| | L | Green | | | | NA25578 | | | | | | |
| | XL | Blue | | | | | | | | | | NA25579 |
| Comfort Amputee*  Material: Polyester | S | Red | NA1201 | | | | | | | | | |
| | M | Yellow | | | | NA1202 | | | | | | |
| | L | Green | | | NA1203 | | | | | | | |
| Comfort Sling Accessory | | | Style | Size | Part Number | | | | | | | |
| Comfort Pad  | Padded | One Size | NA25517 | | | | | | | | | |
| | Standard | One Size | NA25516 | | | | | | | | | |


*Must be used with HOY_AMPUTEE_SB

Hoyer® Loop and Clip Style Sling Color/Size Chart - Standard Size


Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | Cradle | Sling SWL |
|--------------------------------|--|-----------|
| Ascend** Elevate Journey | Loop  | 440 lbs. |


| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | | |
|---|------|--------|--------------------------------|----------|----------|-----|----------|-----|-----|-----|-----|-----|--|--|--|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 440 | | | |
| Deluxe Stand-Aid  Material: Polyester | S | Red | | NA13499N | | | | | | | | | | | |
| | M | Yellow | | | NA13500N | | | | | | | | | | |
| | L | Green | | | | | NA13498N | | | | | | | | |
| Deluxe Transport  Material: Polyester | S | Red | | NA13502N | | | | | | | | | | | |
| | M | Yellow | | | NA13501N | | | | | | | | | | |
| | L | Green | | | | | NA13503N | | | | | | | | |

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | | |
|---|------|--------|--------------------------------|----------|----------|-----|----------|-----|-----|-----|-----|-----|--|--|--|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 440 | | | |
| Standing Disposable  Material: Disposable | S | Red | | NA13525N | | | | | | | | | | | |
| | M | Yellow | | | NA13526N | | | | | | | | | | |
| | L | Green | | | | | NA13527N | | | | | | | | |

**Ascend is discontinued - for reference only


| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | |
|--|------|--------|--------------------------------|----|---------|-----|---------|-----|-----|-----|-----|-----|--|--|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 440 | | |
| Stand-Aid Clip*  Material: Polyester | M | Yellow | | | NA13510 | | | | | | | | | |
| | L | Green | | | | | NA13511 | | | | | | | |




















*Elevate only. Must have 0Y0488 type spreader bar for clip style sling.

| Lift Types | Cradle | Sling SWL |
|------------|--|-----------|
| Elevate | Clip  | 440 lbs. |

Hoyer® Loop Style Sling Color/Size Chart - Bariatric Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | Spreader Bar | Sling SWL |
|----------------|--|-----------|
| Calibre HPL700 | Bariatric Loop  | 850 lbs. |



| Sling | Size | Color | Color Example | Part Number |
|---|------|---------------------------|--|-------------|
| Quickfit  Material: Polyester  | M | Yellow with White Stripes |  | NA1009P |
| | L | Green with White Stripes |  | NA1010P |
| | XL | Blue with White Stripes |  | NA1011P |
| | XXL | Black with White Stripes |  | NA1012P |
| Full Back  Material: Polyester  | M | Yellow with White Stripes |  | NA1072 |
| | L | Green with White Stripes |  | NA1073 |
| | XL | Blue with White Stripes |  | NA1074 |
| | XXL | Black with White Stripes |  | NA1075 |
| Hammock*  Material: Spacer  | XL | Red with White Stripes |  | NA1030 |
| | XXL | Yellow with White Stripes |  | NA1031 |
| | 3XL | Green with White Stripes |  | NA1032 |
| | 4XL | Blue with White Stripes |  | NA1033 |
| | 5XL | Black with White Stripes |  | NA1034 |




*Only smaller sizes should be used with the HPL700 due to clearance over beds.



** Due to the varying nature of bariatric patient shapes and sizes, it becomes increasingly difficult to provide a guide on sling size based on weight alone. Patients should be individually assessed to ascertain the correct size sling.



Hoyer® Loop Style Sling Color/Size Chart

Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | | Spreader Bar |
|----------------------|-------------------------|--|
| Advance Elara HPL450 | HPL500 Presence Stature | Loop  |
| Calibre | HPL700 | Bariatric Loop  |


| Sling | Item Number | Sling SWL | Sling | Item Number | Sling SWL | Sling | Item Number | Sling SWL |
|---|-------------|-----------|---|-------------|-----------|--|-------------|-----------|
| Repositioning, Standard  | NA25710 | 1000 lbs. | Twin Turner, Regular  | NA25705 | 625 lbs. | Limb, Multi Loop  | NA25700 | 400 lbs. |




| Sling | Sling SWL | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | | | | | | |
|--|-----------|------|--------|--------------------------------|----|-----|-----|-----|-----|-----|-----|-----|-----|--|--|--|--|--|
| | | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 | | | | | |
| Standing Harness  Material: Polyester  | 500 lbs | XS | Brown | NA26600 | | | | | | | | | | | | | | |
| | | S | Red | NA26601 | | | | | | | | | | | | | | |
| | | M | Yellow | NA26602 | | | | | | | | | | | | | | |
| | | L | Green | NA26603 | | | | | | | | | | | | | | |

| Sling | Size | Part Number |
|---|------|-------------|
| Slide Sheet  Material: Polyester  | S | NA25525 |
| | M | NA25526 |
| | L | NA25527 |
| | XL | NA25528 |
| | XXL | NA25529 |

Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

| Lift Types | | Spreader Bar | Sling SWL |
|------------------|-------------------|--|-----------|
| HML400 HPL402 | CHLA-2 CHLA-2T | Loop  | 600 lbs. |

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|------|--------|--------------------------------|----|-----|-------|-----|-----|-------|-----|-----|-----|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| U Sling Head Support  Material: Padded | S | Red | 70013 | | | | | | | | | |
| | M | Yellow | | | | 70012 | | | | | | |
| | L | Green | | | | | | | 70011 | | | |
| | XL | Blue | | | | | | | 70010 | | | |
| U Sling Head Support  Material: Mesh | S | Red | 50027 | | | | | | | | | |
| | M | Yellow | | | | 50026 | | | | | | |
| | L | Green | | | | | | | 50025 | | | |
| U Sling  Material: Padded | S | Red | 70003 | | | | | | | | | |
| | M | Yellow | | | | 70002 | | | | | | |
| | L | Green | | | | | | | 70001 | | | |
| | XL | Blue | | | | | | | 70000 | | | |

| Sling | Size | Color | Recommended Weight Range (lbs) | | | | | | | | | |
|---|------|--------|--------------------------------|----|-----|-------|-----|-----|-------|-----|-----|-----|
| | | | 35 | 75 | 100 | 125 | 150 | 175 | 200 | 275 | 350 | 500 |
| U Sling  Material: Mesh | S | Red | 50023 | | | | | | | | | |
| | M | Yellow | | | | 50022 | | | | | | |
| | L | Green | | | | | | | 50021 | | | |
| | XL | Blue | | | | | | | 50020 | | | |
| One Piece  Material: Polyester | S | Red | 70053 | | | | | | | | | |
| | M | Yellow | | | | 70052 | | | | | | |
| | L | Green | | | | | | | 70051 | | | |
| One Piece Commode  Material: Polyester | S | Red | 70058 | | | | | | | | | |
| | M | Yellow | | | | 70057 | | | | | | |
| | L | Green | | | | | | | 70056 | | | |
| | XL | Blue | | | | | | | 70055 | | | |

*One Piece slings may be used for amputees