

Latitude

Pedal Wheelchair

PEDAL
BRODA[®]



Shown with optional rear mag wheels and swing-away lower leg supports



Rocking motion helps calm agitation, creating a soothing effect




Posterior tilt aids in pressure redistribution



Front pivot seat tilt enhances visual orientation and encourages self-mobility with proper foot-on-floor positioning

Patient, Facility, & Caregiver Benefits

- Caregiver activated tilt allows for easy frequent repositioning
- Design provides safety and stability helping to reduce risk of injury
- Rocking provides stimulation; excellent for memory care



Comfort Tension Seating[®] improves postural support and sitting tolerance

*Enhanced mobility and
freedom for memory care*

SPECIFICATIONS

<i>Degrees of Tilt</i>	<i>up to 33°</i>
<i>Weight Capacity</i>	<i>250 lbs. (144 kg)</i>
<i>Seat Widths</i>	<i>15.5", 16", 17", 18.5", 20" (40cm, 45cm, 50cm, 55cm)</i>
<i>Seat Depths</i>	<i>17", 18.5", 20" (43cm, 46cm, 50cm)</i>
<i>Seat Surface Heights</i>	<i>Adjustable 14", 15", 16", 17", 18", 19" (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)</i>
<i>Wheels</i>	<i>5" (12cm) Total Lock Front Caster Wheels 20" or 22" (50cm or 55cm) Rear Mag Wheels 5" (12cm) Fixed Rear Caster Wheels</i>
<i>Lifetime Warranty Frame</i>	<i>10 years on 16 Gauge tubular steel frame</i>
<i>Pressure-Reducing Seat & Back Support</i>	<i>Integrated with Comfort Tension Seating[®]</i>
<i>Standard Padding</i>	<i>Head support cushion, lateral shoulder supports, lower side cushions, seat and back cushion with cover</i>