

**10D** Rotating Green  
Laser Lights

**Newest** Painless Fat  
Removal Choice

From **LuxMaster Slim**



LuxMaster Slim



## Fat loss in a healthier way

»» As the indispensable part of human's endocrine system, fat cells play a critical role in sensing and responding to changes in systemic energy balance and the metabolic cycle. **LuxMaster Slim** does not kill them like liposuction and cryolipolysis. It may sound reasonable to kill the fat cells in weight loss but actually it is harmful and has side-effects. When the fat cells are destroyed, the future fat will be stored and accumulated in other area and may cause asymmetric accumulation of fat.

»» **LuxMaster Slim** opens a new possibility in cellulite reduction. By using the low level laser, fat cells can be shrinked and then consumed as energy and flushed out naturally via the body's lymphatic system when combined with light exercise and sufficient water intake. It is non-invasive and no side-effects reported.

# Low Level Laser Therapy for Fat Loss

**Safest** ←

**Comfortable** ←

**No Downtime** ←



→ **No Freezing**

→ **No Bruising**

→ **No Swelling**

**Fast results (most patients see results in 4 weeks or less)**

# LuxMaster *Slim*



**10** Dancing Laser Light

**360°** Scanning on Larger Area

Shorten Treatment Time and  
Improve the Fat Reduction Efficiency

# LuxMaster *Slim*

## Upgraded Laser Head

-Multi diode laser foldable on both wings can focus multiple lasers and make the slim treatment more flexible

## High Quality Diode Laser

- Imported from Japan
- Double lens to increase the spot area of the laser and concentrate the energy
- High-temperature-resistant materials for continuous running

## Industrial-grade Elevator

-Make the entire laser head fully automatically up and down to adjust to the height of the bed

## 12.1 Inch high-sensitivity Touch Screen

- User-friendly
- More Convenient





# About the Treatment

Luxmaster Slim provides a quick and no downtime treatment. We can resume normal daily activities immediately after the treatment.

We just need to simply lie under the laser on the front for (20)minutes and then on the back for (20) minutes.

Easy!!



# LuxMaster *Slim*



## TECHNICAL PARAMETERS

Maximum reach of laser head	110cm
Angle adjustable of laser wings	100 degree
Weight of laser head	12kg
Maximum reach of Elevator	500mm
Size of screen	12.1 inches
Power of diode	200mw
Wavelength of diode	532nm
Voltage	90v-240v
Power	120w
Number of diode	10pcs



**10D**

**LOW LEVEL LASER THERAPY**

**For Fat Loss**