

METAL UPRIGHT AIR WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- During postoperative or rehab phase

FEATURES AND BENEFITS

- Malleable struts allow you to customized fit
- Rocker bottom sole promotes a more natural gait
- Pneumatic pump allows air compression to be adjusted as needed for a more customized fit.



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

| Size | Men's shoe size US | Women's shoe size US | Part# |
|---------|--------------------|----------------------|---------|
| X-Small | up to 4 | up to 5.5 | MAW0200 |
| Small | 4.5 – 7 | 6 – 8 | MAW0400 |
| Medium | 7.5 – 10.5 | 8.5 – 11.5 | MAW0600 |
| Large | 10.5 – 12.5 | 11.5 – 13.5 | MAW0800 |
| X-Large | 12.5 & up | 13.5 & up | MAW1000 |

FORMFIT PRO ANKLE

Össur Formfit® Pro Ankle is a compressive ankle sleeve with 2 integrated malleoli pads, designed to worn after mild ankle injuries. The dynamic 3-D MotionTech™ knit is supportive and meets certified medical grade compression requirements. Proven superior to leading brands for moisture wicking and breathability, your ankle stays cool during wear. CoolVent™ technology in the area in front of the ankle further enhances breathability and comfort.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

INDICATIONS FOR USE

For Ankle conditions that may benefit from compression, such as:

- Ankle injury rehabilitation
- Mild sprains and strains
- Mild instability


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your ankle cool during wear
- Certified medical grade compression (20-30 mmHg)
- CoolVent™ technology in the area in front of the ankle for added breathability and comfort
- 2 integrated malleoli SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application

REQUIRED MEASUREMENTS



Measure ankle circumference diagonally over the base of the heel and across the ankle bone. Ankle in 90 degrees.



BLUE

| Size | Ankle circumference | | Blue | | Black | |
|------|---------------------|----------|----------|----------|----------|----------|
| | cm | inch | Left | Right | Left | Right |
| XS | 18 - 20 | 7 - 7 | FP201201 | FP201101 | FP200201 | FP200101 |
| S | 20 - 23 | 7 - 9 | FP201202 | FP201102 | FP200202 | FP200102 |
| M | 23 - 26 | 9 - 10¼ | FP201203 | FP201103 | FP200203 | FP200103 |
| L | 26 - 29 | 10¼ - 11 | FP201204 | FP201104 | FP200204 | FP200104 |
| XL | 29 - 32 | 11 - 12 | FP201205 | FP201105 | FP200205 | FP200105 |